# **1,25(OH)<sub>2</sub>-VIT.D-RIA-CT**

***KIP1929***

---

**LOT** : 110217/1

Read entire protocol before use.

# 1,25(OH)<sub>2</sub>-VIT.D-RIA-CT

## I. INTENDED USE

Radioimmunoassay for the *in vitro* quantitative measurement of human 1,25(OH)<sub>2</sub>-Vitamin D (1,25(OH)<sub>2</sub>-Vit.D) in serum and plasma.

## II. GENERAL INFORMATION

- A. **Proprietary name :** DIAsource 1,25(OH)<sub>2</sub>-VIT.D-RIA-CT Kit
- B. **Catalog number :** KIP1929 : 48 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.  
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

**For technical assistance or ordering information contact :**  
Tel: +32 (0) 10 84.99.11 Fax: +32 (0) 10 84.99.90

## III. CLINICAL BACKGROUND

### A. Biological activities

Vitamin D<sub>3</sub> is mainly synthesized in the skin from 7- dehydrocholesterol and is partially from dietary origin. In the liver, Vitamin D<sub>3</sub> is hydroxylated on carbon 25 to produce the obligatory intermediate 25-OH-D<sub>3</sub>. 25-OH-D<sub>3</sub> must be metabolized further before it can carry out the functions of Vitamin D on intestine, kidney and bone. This subsequent reaction takes place exclusively in the kidney in the non-pregnant mammal. Thus 25-OH-D<sub>3</sub> is further hydroxylated in the 1 $\alpha$ -position to produce 1 $\alpha$ ,25 dihydroxyvitamin D<sub>3</sub> (1 $\alpha$ ,25-(OH)<sub>2</sub>D<sub>3</sub>).

In addition to renal tissue, placenta of pregnant women and macrophage cells in case of sarcoidosis can also produce some amount of 1 $\alpha$ ,25-(OH)<sub>2</sub>D<sub>3</sub>. 1 $\alpha$ ,25-(OH)<sub>2</sub>D<sub>3</sub> is the active form of Vitamin D with regard to the known functions whereas 25-OH-D<sub>3</sub> and Vitamin D<sub>3</sub> itself can be excluded as being physiologically functional. Furthermore since 1 $\alpha$ ,25-(OH)<sub>2</sub>D<sub>3</sub> is produced in the kidney and has some of its functions in the bone and intestine, it must be considered as a hormone. This hormone stimulates the intestinal absorption of both calcium and phosphorus. It also stimulates bone resorption and mineralization thereby preventing the development of rickets and osteomalacia.

1 $\alpha$ ,25-(OH)<sub>2</sub>D<sub>3</sub> might also be active in other tissues responsible for Calcium transport (placenta, kidney, mammary gland, ...) and endocrine glands such as parathyroid glands. 1 $\alpha$ ,25-(OH)<sub>2</sub>D<sub>3</sub> is rapidly metabolized and its lifetime is approximately 2-4 h in plasma. Its main metabolite is calcitroic acid, a C-23 carboxylic derivative essentially without any biological activity. In addition to this pathway, 1 $\alpha$ ,25-(OH)<sub>2</sub>D<sub>3</sub> undergoes 24-hydroxylation to produce 1,24,25-trihydroxy-Vitamin D<sub>3</sub>. This compound has less biological activity than its parent and this metabolism is considered as a minor pathway.

The levels of 1 $\alpha$ ,25-(OH)<sub>2</sub>D<sub>3</sub> in plasma or serum is 100 to 1000 less than that of 25-OH-D<sub>3</sub>. Due to its low concentrations and the presence of many similar metabolites, the measurement of 1 $\alpha$ ,25-(OH)<sub>2</sub>D<sub>3</sub> requires extraction and separation either by HPLC or by column chromatography.



### B. Clinical application

The measurement of circulating 1 $\alpha$ ,25-(OH)<sub>2</sub>D<sub>3</sub> is indicated in several disorders affecting calcium metabolism such as : sarcoidosis, renal failure, hyper and hypo-parathyroidism, rickets, tumor-associated hypercalcemia, Vitamin-resistant dysfunction and treatment with anti-convulsive medication.

#### IV. PRINCIPLES OF THE METHOD

Only samples and controls, not the calibrators, are extracted with a mix of solvents and applied on cartridges to separate 1,25(OH)<sub>2</sub> Vitamin-D from other Vitamin-D metabolites. After elution of samples and controls, the calibrators, samples and controls are incubated in coated tubes. A fixed amount of <sup>125</sup>I labelled 1,25(OH)<sub>2</sub> Vitamin D competes with the 1,25(OH)<sub>2</sub> Vitamin D to be measured present in the sample or in the calibrator for a fixed amount of antibody sites immobilized on the wall of a polystyrene tube. After an overnight incubation at room temperature, an aspiration step terminates the competition reaction. The tubes are then washed with washing solution and aspirated. A calibration curve is plotted and the 1,25(OH)<sub>2</sub> Vitamin D concentrations of the samples are determined by dose interpolation from the calibration curve.

#### V. REAGENTS PROVIDED

Reagents	48 Test Kit	Colour Code	Reconstitution
 Tubes coated with anti 1,25(OH) <sub>2</sub> -Vitamin D	1 x 48	green	Ready for use
 TRACER: <sup>125</sup> Iodine labelled 1,25(OH) <sub>2</sub> -Vitamin D (HPLC grade) in phosphate buffer with bovine casein and gentamycin.	1 vial lyophilised 75 kBq	red	Add 26 ml reconstitution solution
 Calibrators - N = 1 to 5 (see exact values on vial labels) in phosphate buffer with bovine casein and gentamycin	5 vials lyophilised	yellow	Add 2 ml elution solution
 Wash solution (TRIS-HCl)	1 vial 10 ml	brown	Dilute 70 x with distilled water (use a magnetic stirrer).
 Controls - N = 1 or 2 in human plasma with gentamycin	2 vials lyophilised	silver	Add 2 ml distilled water
 Reconstitution Solution: phosphate buffer with bovine casein and azide (<0.1%)	1 vial 30 ml	black	Ready for use
 Elution Solution: phosphate buffer with bovine casein, methanol and azide (<0.1%)	1 vial 30 ml	green	Ready for use
 Bond Elut Silica cartridges	20		Store at R.T.

**Note :** Use elution solution for calibrator 0 and for dilution of samples with values above the highest calibrator (dilute after separation step).

#### VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

- Distilled water
- Diisopropylether (p.a.)
- Cyclohexane (p.a.)
- Ethyl acetate (p.a.)
- Ethanol absolute (p.a.)
- Dichloromethane (p.a.)
- Pipettes for delivery of: 200 µl, 500 µl, 1 ml and 2 ml (the use of accurate pipettes with disposable plastic tips is recommended)
- Glass tubes (12 x 75 mm) for extraction and for elution. (closed with a cap for the extraction step)
- Glass tubes (16 x 100 mm) or (12 x 120 mm), or polypropylene tubes (falcon 2097), for the washing of the cartridges.
- Vortex mixer
- Magnetic stirrer
- Centrifuge operating at 800 g.

- Tube shaker (1200 rpm)
- 5 ml automatic syringe (Cornwall type) for washing
- Aspiration system (optional)
- Any gamma counter capable of measuring <sup>125</sup>I may be used (minimal yield 70%).

#### VII. REAGENT PREPARATION

- Calibrators:** Reconstitute the calibrators with 2 ml elution solution (**just before the incubation step**).
- Controls:** Reconstitute the controls with 2 ml distilled water.
- <sup>125</sup>I 1,25(OH)<sub>2</sub>-Vitamin.D :** Reconstitute with 26 ml of reconstitution solution.
- Working Wash solution:** Prepare an adequate volume of Working Wash solution by adding 69 volumes of distilled water to 1 volume of Wash Solution (70x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.
- Extraction solvent :** 2 ml for each control or sample to be tested, are needed. **Prepare a fresh solution** of diisopropylether, cyclohexane, ethyl acetate, (50, 40, 10 v/v).
- Washing solvent :** 1 ml for each control or sample to be tested, are needed. **Prepare a fresh solution** of diisopropylether, cyclohexane, ethyl acetate, ethanol absolute (50, 40, 10, 1 v/v).

#### VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C; except the cartridges which must be stored at room temperature.
- The calibrators and controls are very unstable, use them immediately after reconstitution, freeze immediately in aliquots and keep them at -20°C for 3 months. Avoid subsequent freeze-thaw cycles.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.
- Use freshly prepared extraction solvent and washing solvent, do not store them.**
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

#### IX. SPECIMEN COLLECTION AND PREPARATION

- Serum and plasma samples must be kept at 2-8°C.
- If the test is not run within 24 hrs, storage in aliquots, at -20°C is recommended.
- Avoid subsequent freeze-thaw cycles.
- After thawing, the samples should be vortexed and centrifuged.
- Serum or plasma (EDTA and heparin) provides similar results.

#### X. PROCEDURE

- Handling notes**  
Do not use the kit or components beyond expiry date.  
Do not mix materials from different kit lots.  
Bring all the reagents to room temperature prior to use.  
Thoroughly mix all reagents and samples by gentle agitation or swirling.  
Use a clean disposable pipette tip for addition of each different reagent and sample in order to avoid cross-contamination. High precision pipettes or automated pipetting equipment will improve the precision.  
Respect the incubation times.  
Prepare a calibration curve for each run, do not use data from previous runs.

#### B. Procedure

- Extraction step : ! Only for controls and samples.**
  - Label glass tubes (12x75 mm) for extraction: 2 controls and up to 16 samples.
  - Add 0.5 ml control or sample in the respective tubes.
  - Dispense 2 ml extraction solvent in each tube.
  - Tubes are closed with a cap and placed on a shaker for 1 hour at 1200 rpm.
  - Centrifuge each tube for 5 minutes at room temperature (at 800 g).
  - Supernatants are needed for the next step of separation.
- Separation step : ! Only for controls and samples**
  - Label glass tubes (16 x 100 mm) or (12 x 120 mm), or polypropylene tubes (falcon 2097), for washing cartridges: 2 controls and up to 16 samples.
  - Put one "Bond Elut" cartridge in each tube.
  - Apply 1.6 ml of supernatant (2 x 0.8 ml), obtained after extraction step, on cartridge.

- Then, wash cartridges with 1 ml washing solvent (cfr reagent preparation). ! Be careful never apply vacuum on cartridges, just let solvent draw by gravity.
- Add 300 µl dichloromethane on each cartridge, let draw by gravity.
- Add 300 µl of distilled water on each cartridge.
- Centrifuge each tube for 5 minutes at room temperature (at 800 g).
- Label glass tubes (12 x 75 mm) for elution of 1,25(OH)<sub>2</sub>-Vitamin D. After centrifugation, transfer cartridges in the corresponding glass tubes.
- Apply 400 µl elution solution on each cartridge to elute 1,25 (OH)<sub>2</sub>-Vitamin D and centrifuge 5 minutes at room temperature (at 800 g).
- Vortex** the eluted fraction.

**Note :** After this step, samples must be incubated in coated tubes as soon as possible to avoid degradation.

### III. Incubation step :

- Label coated tubes in duplicate for each calibrator, control and sample. For the determination of total counts, label 2 normal tubes
- Briefly vortex calibrators (use elution solution as zero calibrator), extracted controls and samples and dispense 150 µl of each into the respective tubes.
- Dispense 500 µl of <sup>125</sup>Iodine labelled 1,25(OH)<sub>2</sub>-Vitamin D into each tube, including the uncoated tubes for total counts.
- Shake the tube rack gently by hand to liberate any trapped air bubbles.
- Incubate overnight at room temperature
- Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
- Wash tubes with 2 ml Working Wash solution (except total counts) and aspirate (or decant). Avoid foaming during the addition of the Working Wash solution.
- Aspirate (or decant) the content of each tube (except total counts).
- Wash tubes again with 2 ml Wash solution (except total counts) and aspirate (or decant).
- After the last washing, let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
- Count tubes in a gamma counter for 60 seconds.

### XI. CALCULATION OF RESULTS

- Calculate the mean of duplicate determinations.
- Calculate the bound radioactivity as a percentage of the binding determined at the zero calibrator point (0) according to the following formula :

$$B/B_0(\%) = \frac{\text{Counts (Calibrator or sample)}}{\text{Counts (Zero Calibrator)}} \times 100$$

- Using a 3 cycle semi-logarithmic or logit-log graph paper, plot the (B/B<sub>0</sub>(%)) values for each calibrator point as a function of the 1,25(OH)<sub>2</sub>-Vitamin D concentration of each calibrator point. Reject obvious outliers.
- Computer assisted methods can also be used to construct the calibration curve. If automatic result processing is used, a 4-parameter logistic function curve fitting is recommended.
- By interpolation of the sample (B/B<sub>0</sub> (%)) values, determine the 1,25(OH)<sub>2</sub>-Vitamin D concentrations of the samples from the calibration curve.
- For each assay, the percentage of total tracer bound in the absence of unlabelled 1,25(OH)<sub>2</sub>-Vitamin D (B<sub>0</sub>/T) must be checked.

### XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

1,25(OH) <sub>2</sub> -Vitamin D	cpm	B/Bo (%)
Total count	51237	
Calibrator		
0.0 pg/ml	20425	100.0
4.7 pg/ml	17833	87.3
12.2 pg/ml	15921	77.9
60.4 pg/ml	10602	51.9
160.0 pg/ml	7150	35.0
411.0 pg/ml	4587	22.5

### XIII. PERFORMANCE AND LIMITATIONS

#### A. Detection limit

Twenty zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent concentration two standard deviations below the average counts at zero binding, was 1.6 pg/ml.

#### B. Specificity

The percentage of cross-reaction estimated by comparison of the concentration yielding a 50% inhibition are respectively:

Compound	Cross-Reactivity (%)
1,25(OH) <sub>2</sub> -Vitamin.D3	100
1,25(OH) <sub>2</sub> -Vitamin.D2	70
25OH-Vitamin-D3	0.0004
24,25(OH) <sub>2</sub> -Vitamin.D3	0.0027
25,26(OH) <sub>2</sub> -Vitamin.D3	0.007

**Note:** this table shows the cross-reactivity for the anti 1,25(OH)<sub>2</sub>-Vitamin D

#### C. Precision

##### INTRA-ASSAY PRECISION

##### INTER-ASSAY PRECISION

Serum	N	INTRA-ASSAY PRECISION		INTER-ASSAY PRECISION			
		<X> ± SD (pg/ml)	CV (%)	Serum	N	<X> ± SD (pg/ml)	CV (%)
A	20	24.5 ± 2.3	9.3	A	10	13.6 ± 1.7	12.7
B	20	77.3 ± 3.5	4.5	B	10	33.4 ± 3.6	11.3

SD: Standard Deviation; CV: Coefficient of variation

#### 7. D. Accuracy

##### DILUTION TEST

Sample	Dilution	Theoretical Concent. (pg/ml)	Measured Concent. (pg/ml)
	1/1	-	70.0
	1/2	35.0	35.7
	1/4	17.5	14.5
	1/8	8.8	7.8
	1/16	4.4	4.6

The sample was diluted with Elution solution.

##### RECOVERY TEST

Sample	added 1,25(OH) <sub>2</sub> -Vitamin D (pg/ml)	Recovered 1,25(OH) <sub>2</sub> -Vitamin D (pg/ml)	Recovered (%)
Serum n°1	25	23.8	95
	50	47.5	95
	100	100.2	100
Serum n° 2	25	29.6	118
	50	47.9	96
	100	90.4	90

#### Conversion factor :

From pg/ml to pmol/L : x 2.4  
From pmol/L to pg/ml : x 0.42

To the best of our knowledge, no international reference material exists for this parameter.

### XIV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises.

#### XV. REFERENCE INTERVALS

These values are given only for guidance; each laboratory should establish its own normal range of values.

The observed ranges are based on 2.5% to 97.5% percentiles.

Population	Range (pg/ml)	Mean	SD	n
Normal subjects	19.6 – 54.3	35.3	10.6	51

#### XVI. PRECAUTIONS AND WARNINGS

##### Safety

For *in vitro* diagnostic use only.

This kit contains <sup>125</sup>I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations. This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HbsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

#### XVII. BIBLIOGRAPHY

- Bouillon R.A., Auwerx J.D., Lissens W.D. and Pelemans W.K. (1987) **Vitamin D status in elderly; season substrate deficiency causes 1,25-dihydroxycholecalciferol deficiency.** Am. J. Clin. Nutr., 45:755-763
- Iqbal, S.J. (1994). **Vitamin D metabolism and the clinical aspects of measuring metabolites.** Ann. Clin. Biochem., 31:109-124
- Mawer E.B. (1980). **Clinical implications of measurements of circulating vitamin D metabolites.** Clinics in Endocr. Metabol., 9:63-79

- Jongen M.J.M., Van Ginkel F.C., Vander Vijgh W.J.F., Kuiper S., Netelenbos J.C. and Lips P; (1984). **An international comparison of Vitamin D metabolites measurements.** Clin. Chem., 30:399-403
- Deluca H.F. (1979). **The Vitamin D system in the regulation of calcium and phosphorus metabolism.** Nutritional Rev., 37:161-193
- Haussler, M.R., McCain, T.A. (1977). **Basic and clinical concepts related to Vitamin D metabolism and action.** N. Engl. J. Med., 297:974-983

#### XVIII. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS µl	CALIBRATORS µl	SAMPLE (S) CONTROLES µl
<b>EXTRACTION</b>			
Calibrators	-	-	-
Samples / Controls	-	-	500
Extraction solvent	-	-	2000
Shaking	1 hour at 1200 rpm		
Centrifugation	5 minutes at 800 g		
<b>SEPARATION</b>			
Supernatant from extraction step	-	-	1600
<b>CARTRIDGE</b>			
Supernatant	1600 µl		
Washing Solvent	1000 µl		
Dichloromethane	300 µl		
Distilled water	300 µl		
Centrifugation	5 minutes at 800 g		
Elution solution	400 µl		
Centrifugation	5 minutes at 800 g		
	<b>Vortex</b>		
<b>INCUBATION</b>			
Calibrators	-	150	-
Extracted samples	-	-	150
Tracer	500	500	500
Incubation	Overnight at R.T.		
Separation	-	Aspirate (or decant)	
Washing Solution	-	2 ml	
Separation	-	Aspirate (or decant)	
Washing Solution	-	2 ml	
Separation	-	Aspirate (or decant)	
Counting	Count tubes for 60 seconds in a gamma counter.		

DIASource Catalogue Nr : KIP1929	P.I. Number : 1700602/en	Revision nr : 110217/1
-------------------------------------	-----------------------------	---------------------------

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	<b>Used symbols</b>
	Consult instructions for use
	Storage temperature
	Use by
<b>LOT</b>	Batch code
<b>REF</b>	Catalogue number
<b>CONTROL</b>	Control
<b>I V D</b>	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
WASH SOLN CONC	Wash solution concentrated
CAL 0	Zero calibrator
CAL N	Calibrator #
CONTROL N	Control #
Ag 1251	Tracer
Ab 1251	Tracer
Ag 1251 CONC	Tracer concentrated
Ab 1251 CONC	Tracer concentrated
	Tubes
INC BUF	Incubation buffer
ACETONITRILE	Acetonitrile
SERUM	Serum
DIL SPE	Specimen diluent
DIL BUF	Dilution buffer
ANTISERUM	Antiserum
IMMUNOADSORBENT	Immunoabsorbent
DIL CAL	Calibrator diluent
REC SOLN	Reconstitution solution
PEG	Polyethylene glycol
EXTR SOLN	Extraction solution
ELU SOLN	Elution solution
GEL	Bond Elut Silica cartridges
PRE SOLN	Pre-treatment solution
NEUTR SOLN	Neutralization solution
TRACEUR BUF	Tracer buffer
	Microtiterplate
Ab HRP	HRP Conjugate
Ag HRP	HRP Conjugate
Ab HRP CONC	HRP Conjugate concentrate
Ag HRP CONC	HRP Conjugate concentrate
CONJ BUF	Conjugate buffer
CHROM TMB CONC	Chromogenic TMB concentrate
CHROM TMB	Chromogenic TMB solution
SUB BUF	Substrate buffer
STOP SOLN	Stop solution
INC SER	Incubation serum
BUF	Buffer
Ab AP	AP Conjugate
SUB PNPP	Substrate PNPP
BIOT CONJ CONC	Biotin conjugate concentrate
AVID HRP CONC	Avidine HRP concentrate
ASS BUF	Assay buffer
Ab BIOT	Biotin conjugate
Ab	Specific Antibody
SAV HRP CONC	Streptavidin HRP concentrate
NSB	Non-specific binding
2nd Ab	2nd Antibody
ACID BUF	Acidification Buffer
DIST	Distributor

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