Free TESTO-RIA-CT

KIPI19000

LOT : 110228/1

Free TESTOSTERONE-RIA-CT

Radioimmunoassay for the Quantitative Determination of Free Testosterone
in Human Serum
KIPI19000
IN VITRO DIAGNOSTIC USE

en

DIAsource ImmunoAssays SA - Rue du Bosquet 2, B-1348 Louvain-la-Neuve, Belgium - Tel: +32 10 84 99 11 - Fax : +32 10 84 99 90

1. INTENDED USE

For **IN VITRO** determination of Free Testosterone (FT) levels in hirsutism and hypogonadism.

Free testosterone diffuses through cell membranes and binds to specific receptor proteins (androgen receptors); the Testosterone-receptor complexes act as transcriptional modulators on cis-regulatory regions of many genes.

Excess of Androgens in women causes hirsutism and signs of virilization; Testosterone level in serum has to be determined before and after ovarian and adrenal stimulation and suppression to identify the source of excessive hormone production.

Primary and secondary hypogonadism in men result in clinical hypoandrogenization, correlated with the degree of gonadal failure in Testosterone production. The determination of serum Testosterone together with that of LH allows the correct assessment of those conditions.

The diagnosis of true anorchia also requires to discriminate this condition from cryptorchidism. Under prolonged hCG stimulation, Testosterone levels remain very low in true anorchia while cryptorchid testes can respond to stimulation.

Androgen resistance syndromes, due to X linked androgen receptor gene deficiencies, are made of various degrees of sexual ambiguity. Whatever the severity of the phenotypical abnormalities, serum Testosterone is systematically high in regards to elevated LH serum levels in these conditions.


Testosterone assays include total testosterone (direct, extraction, coated tubes) and free testosterone determinations.

Total Testosterone in plasma includes free Testosterone and Testosterone bound to SHBG, albumin, CBG. The mean percentage of each in normal men is 2.7, 32, 65 and <0.1 respectively.

Solvents break the protein binding in extraction assays whereas blocking agents release Testosterone from proteins in direct assays. The advantage of a free testosterone assay is that free testosterone concentrations are in equilibrium with testosterone bound to receptors in the organs.

2. PRINCIPLE OF THE METHOD

The Free Testosterone (FT) CT RIA obeys the law of mass action according to the following equation :

Since the concentrations of ${}^{125}\text{I} - \text{FT}$ and coated antibodies are constant, the advancing state of the equation depends on the concentration of FT. The amount of ${}^{125}\text{I} - \text{FT}$ bound to the coated tube is inversely proportional to the concentration of FT in the sample.

Following the incubation, the tube is aspirated to remove excess unbound labelled T.

Patient sample concentrations are read from a calibration curve.

3. MATERIAL PROVIDED AND STORAGE :

Stored at 2 - 8°C, the material can be used up to the expiration date printed on each label.

- 3.1.

--

 2 x 48 Polystyrene tubes (12 x 75 mm) coated with anti-Testosterone polyclonal antibodies. Systematically allow the coated tubes to reach room temperature before use.
- 3.2.

Ag	125I
----	------

 yellow, 42 ml
1 bottle of ${}^{125}\text{I}$ -labelled FREE TESTOSTERONE analogue in protein based buffer containing < 0.1 % NaN_3 as preservative.
Each bottle contains less than 185 Kbcq (5 μCi)
- 3.3.

CAL	N
-----	---

 0.5 ml in each vial - N=0 to 6
7 vials of FREE TESTOSTERONE in human serum containing preservative ($\text{NaN}_3 < 0.1 \%$).
The concentrations are printed on the labels.

- 3.4.

CONTROL	N
---------	---

 0.5 ml in each vial - N=1 or 2
2 vials of human serum containing preservative ($\text{NaN}_3 < 0.1 \%$). The control sera are to be assayed along with the patient samples. The ranges for the control sera are printed on the vial labels.

- 3.5.

WASH	SOLN	CONC
------	------	------

 70 x concentrated, 10 ml
1 bottle concentrated buffered solution containing sodium azide ($\text{NaN}_3 < 0.1 \%$). Pour the solution in 700 ml of distilled water.

4. MATERIAL REQUIRED BUT NOT PROVIDED :

- bench surfaces, protected by absorbent paper to reduce the effects of radioactive spillage.
- waste disposal containers, appropriately labelled and suitable for solid or liquid radioactive materials.
- manual or automated precision micropipettes for dispensing samples or reagents without cross-contamination.
- absorbent paper.
- vacuum pump, connected through a trap, for aspiration.
- water bath.
- a gamma scintillation counter
- appropriate graph paper for plotting the results.

5. METHODOLOGY

5.1. Collection and handling of blood samples :

The blood sample can be collected into a dry tube.

After separation from the red blood cells, serum samples can be assayed immediately, within 24 hours if stored at 2 - 8°C, or later, after a period of up to several months if stored at -20°C. Repeatedly freezing and thawing must be avoided.

5.2. Assay procedure :

Reagents stored at 2°- 8° C. must be brought at room temperature prior to use. Do not mix reagents of different lots. Label the tubes for T (« Total Counts » do not use coated tubes) calibrators, samples and controls. Calibrators and controls should be mixed before use by inverting or swirling rather than vortexing.

Perform the assay in duplicate. Calibrators, controls and samples must be assayed at the same time.

1. Calibrator curve :

Pipette 50 μl of each calibrator into the corresponding tubes.

2. Samples and control sera :

Pipette 50 μl of each sample or control serum into the corresponding tubes.

3. Add 400 μl of ${}^{125}\text{I} - \text{TESTOSTERONE}$ analogue tracer to each tube. Vortex and cover.

4. Incubate 2 hours at $37 \pm 2^\circ\text{C}$.

5. **Carefully** aspirate or decant (before to decant, add 2 ml of washing solution to each tube) the solution of all tubes. (Except total counts tubes).

6. Add 2 ml of washing solution to each tube. Aspirate or decant carefully.

7. Count the radioactivity fixed in each tube for at least 60 seconds

5.3. Data processing :

Determine the mean count rate for each set of duplicate tubes.

Calculate the ratio B/B0 as follows :

$$B/B0 \% = \frac{\text{Cal or Smp cpm}}{\text{B0 (Cal 0) cpm}} \times 100$$

Draw the calibrator curve on semilogarithmic paper by plotting the ratio B/B0 % (linear scale) obtained for each calibrator versus its respective concentration expressed in pg/ml (logarithmic scale). FREE TESTOSTERONE concentrations in samples can be read directly from the calibrator curve.

If a computer is used to calculate the results, the data can be fitted to the appropriate equation : smoothed spline.

5.4. Example of a typical assay :

	Contents (pg/ml)	cpm 1st duplicate	cpm 2nd duplicate	Mean count rate	B/Bo (%)	Free Testosterone (pg/ml)
Total counts	-	52039	51647	51843	-	-
Cal 0	0	25839	25961	25900	100	-
Cal 1	0.3	21086	21170	21128	81.6	-
Cal 2	1	16509	16203	16356	63.2	-
Cal 3	3	12437	12428	12433	48	-
Cal 4	10	8317	7916	8117	31.3	-
Cal 5	30	5017	4711	4864	18.8	-
Cal 6	90	2611	2528	2570	9.9	-
C 1 low	1.5 – 2.9	13461	13557	13508	52.2	2.3
C 2 high	15 – 29	5736	5002	5369	20.8	21
Sample 1		17478	16742	16975	65.5	0.86
Sample 2		7538	7423	7481	28.9	11.5
Sample 3		4681	4645	4663	24.3	33

Example of a typical assay, do not use for calculations

6. PERFORMANCE CHARACTERISTICS :

6.1. Specificity

Steroid	% Cross-reactivity
Testosterone	100
5 α DHT	0.006
androstenedione	0.02
β estradiol	0.0003
DHEA-S	0.000001
Androsterone, Corticosterone, 11 DOC, estriol, estrone, progesterone, DHEA	N.D.

6.2. Minimum detectable concentration of FREE TESTOSTERONE :

The minimum detectable concentration has been assayed at 0.13 pg/ml and corresponds to the concentration given by two standard deviations below the mean cpm of 20 replicate determinations of the zero calibrator.

6.3. Reproducibility :

	Mean value (pg/ml)	Within assay variation (% CV) 10 replicates	Between assay variation (% CV) 7 Separate assays in duplicate
Pool 1	0.73	11.4	18.07
Pool 2	10.89	5.7	6.72
Pool 3	33.94	9.3	8.46

7. LIMITATION OF THE PROCEDURE

- The results obtained from this or any other diagnostic kit should be used and interpreted only in the context of an overall clinical picture.
- Do not use lipemic, haemolyzed, icteric or turbid specimens.
- Do not use plasma samples

8. EXPECTED VALUES

It is recommended that each laboratory establishes its own reference values.

Revision date: 2011-02-28

Age group	MALES (pg/ml)		FEMALES (pg/ml)	
	Median	Range*	Median	Range*
<20	-	0,2 – 42,5	0,95	ND – 3,09
20 - 39	22	8,9 – 42,5	1	ND – 3,09
40 - 59	16,2	6,6 – 30,0	0,9	ND – 2,60
>60	11,9	4,9 – 21,6	0,7	ND – 1,80

* 95% percentiles

9. WARNING AND PRECAUTIONS

For IN VITRO DIAGNOSTIC use only

CAUTION : Radioactive material

This kit contains ¹²⁵I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

WARNING : Sodium azide

Some components contain sodium azide as preservative agent (NaN₃ < 0,1%). Dispose of the reagents by flushing with large amount of water through the plumbing system.

WARNING : Potentially infectious material

Handle all components (and all patient samples) as if capable of transmitting viral diseases such as hepatitis B and C and the acquired immunodeficiency syndrome (AIDS).

Source material derived from human body fluids or organs and used in the preparation of this kit were tested and found negative for HBsAg and anti-HCV by immunoassay. However, no known test can guarantee that such material does not contain the causative agent of viral hepatitis.

Likewise, all human materials used in the preparation of this kit were screened for the presence of antibodies against HIV-1 and -2 by enzyme-immunoassay and were found negative. However, absence of this antibody cannot guarantee the absence of the viral agent responsible for the acquired immunodeficiency syndrome.

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
Data wydania: 2011-02-28

	Used symbols
	Consult instructions for use
	Storage temperature
	Use by
LOT	Batch code
REF	Catalogue number
CONTROL	Control
I V D	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
WASH SOLN CONC	Wash solution concentrated
CAL 0	Zero calibrator
CAL N	Calibrator #
CONTROL N	Control #
Ag 1251	Tracer
Ab 1251	Tracer
Ag 1251 CONC	Tracer concentrated
Ab 1251 CONC	Tracer concentrated
	Tubes
INC BUF	Incubation buffer
ACETONITRILE	Acetonitrile
SERUM	Serum
DIL SPE	Specimen diluent
DIL BUF	Dilution buffer
ANTISERUM	Antiserum
IMMUNOADSORBENT	Immunoabsorbent
DIL CAL	Calibrator diluent
REC SOLN	Reconstitution solution
PEG	Polyethylene glycol
EXTR SOLN	Extraction solution
ELU SOLN	Elution solution
GEL	Bond Elut Silica cartridges
PRE SOLN	Pre-treatment solution
NEUTR SOLN	Neutralization solution
TRACEUR BUF	Tracer buffer
µP	Microtiterplate
Ab HRP	HRP Conjugate
Ag HRP	HRP Conjugate
Ab HRP CONC	HRP Conjugate concentrate
Ag HRP CONC	HRP Conjugate concentrate
CONJ BUF	Conjugate buffer
CHROM TMB CONC	Chromogenic TMB concentrate
CHROM TMB	Chromogenic TMB solution
SUB BUF	Substrate buffer
STOP SOLN	Stop solution
INC SER	Incubation serum
BUF	Buffer
Ab AP	AP Conjugate
SUB PNPP	Substrate PNPP
BIOT CONJ CONC	Biotin conjugate concentrate
AVID HRP CONC	Avidine HRP concentrate
ASS BUF	Assay buffer
Ab BIOT	Biotin conjugate
Ab	Specific Antibody
SAV HRP CONC	Streptavidin HRP concentrate
NSB	Non-specific binding
2nd Ab	2nd Antibody
ACID BUF	Acidification Buffer
DIST	Distributor

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	Used symbols
	Consult instructions for use
	Storage temperature
	Use by
LOT	Batch code
REF	Catalogue number
CONTROL	Control
I V D	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
WASH SOLN CONC	Wash solution concentrated
CAL 0	Zero calibrator
CAL N	Calibrator #
CONTROL N	Control #
Ag 1251	Tracer
Ab 1251	Tracer
Ag 1251 CONC	Tracer concentrated
Ab 1251 CONC	Tracer concentrated
	Tubes
INC BUF	Incubation buffer
ACETONITRILE	Acetonitrile
SERUM	Serum
DIL SPE	Specimen diluent
DIL BUF	Dilution buffer
ANTISERUM	Antiserum
IMMUNOADSORBENT	Immunoabsorbent
DIL CAL	Calibrator diluent
REC SOLN	Reconstitution solution
PEG	Polyethylene glycol
EXTR SOLN	Extraction solution
ELU SOLN	Elution solution
GEL	Bond Elut Silica cartridges
PRE SOLN	Pre-treatment solution
NEUTR SOLN	Neutralization solution
TRACEUR BUF	Tracer buffer
	Microtiterplate
Ab HRP	HRP Conjugate
Ag HRP	HRP Conjugate
Ab HRP CONC	HRP Conjugate concentrate
Ag HRP CONC	HRP Conjugate concentrate
CONJ BUF	Conjugate buffer
CHROM TMB CONC	Chromogenic TMB concentrate
CHROM TMB	Chromogenic TMB solution
SUB BUF	Substrate buffer
STOP SOLN	Stop solution
INC SER	Incubation serum
BUF	Buffer
Ab AP	AP Conjugate
SUB PNPP	Substrate PNPP
BIOT CONJ CONC	Biotin conjugate concentrate
AVID HRP CONC	Avidine HRP concentrate
ASS BUF	Assay buffer
Ab BIOT	Biotin conjugate
Ab	Specific Antibody
SAV HRP CONC	Streptavidin HRP concentrate
NSB	Non-specific binding
2nd Ab	2nd Antibody
ACID BUF	Acidification Buffer
DIST	Distributor

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