IGF-1-ELISA

KAP1581

History

Summary of change:

Previous Version: 200224-1	Current Version: 220516																																																																																																
Old DIAsource logo 	New DIAsource logo 																																																																																																
V. REAGENTS PROVIDED	V. REAGENTS PROVIDED																																																																																																
<table border="1"> <tr> <td>CONJ</td> <td>BUF</td> <td>1 vial 11.5 ml</td> <td>Ready for use</td> </tr> <tr> <td colspan="4">Conjugate buffer: Phosphate buffer with caseine and thymol</td> </tr> <tr> <td>CAL</td> <td>0</td> <td>1 vial lyophilized</td> <td>Add 1 ml reconstitution solution</td> </tr> <tr> <td colspan="4">Zero calibrator in phosphate buffer with bovine casein</td> </tr> <tr> <td>CAL</td> <td>N</td> <td>5 vials lyophilized</td> <td>Add 1 ml reconstitution solution</td> </tr> <tr> <td colspan="4">Calibrators - N = 1 to 5 (see exact values on vial labels) in phosphate buffer with bovine casein</td> </tr> <tr> <td>STOP</td> <td>SOLN</td> <td>1 vial 12 ml</td> <td>Ready for use</td> </tr> <tr> <td colspan="4">Stop solution: HCl 1M</td> </tr> <tr> <td>CHROM</td> <td>TMB</td> <td>1 vial 25 ml</td> <td>Ready for use</td> </tr> <tr> <td colspan="4">Chromogenic TMB Solution (Tetramethylbenzidine)</td> </tr> <tr> <td>NEUTR</td> <td>SOLN</td> <td>1 vial 30 ml</td> <td>Ready for use</td> </tr> <tr> <td colspan="4">Neutralization Solution containing phosphate buffer with bovine casein and azide (<0.1%)</td> </tr> </table>	CONJ	BUF	1 vial 11.5 ml	Ready for use	Conjugate buffer: Phosphate buffer with caseine and thymol				CAL	0	1 vial lyophilized	Add 1 ml reconstitution solution	Zero calibrator in phosphate buffer with bovine casein				CAL	N	5 vials lyophilized	Add 1 ml reconstitution solution	Calibrators - N = 1 to 5 (see exact values on vial labels) in phosphate buffer with bovine casein				STOP	SOLN	1 vial 12 ml	Ready for use	Stop solution: HCl 1M				CHROM	TMB	1 vial 25 ml	Ready for use	Chromogenic TMB Solution (Tetramethylbenzidine)				NEUTR	SOLN	1 vial 30 ml	Ready for use	Neutralization Solution containing phosphate buffer with bovine casein and azide (<0.1%)				<p>Removal of the column relative to color code</p> <table border="1"> <tr> <td>CONJ</td> <td>BUF</td> <td>1 vial 13 ml</td> <td>Ready for use</td> </tr> <tr> <td colspan="4">Conjugate buffer: Phosphate buffer with caseine and thymol</td> </tr> <tr> <td>CAL</td> <td>0</td> <td>1 vial 1 ml</td> <td>Ready for use</td> </tr> <tr> <td colspan="4">Zero calibrator in phosphate buffer with ovalbumin and gentamycin</td> </tr> <tr> <td>CAL</td> <td>N</td> <td>5 vials lyophilized</td> <td>Add 1 ml reconstitution solution</td> </tr> <tr> <td colspan="4">Calibrators - N = 1 to 5 in phosphate buffer with ovalbumin and gentamycin. (see exact values on vial labels)</td> </tr> <tr> <td>STOP</td> <td>SOLN</td> <td>1 vial 13 ml</td> <td>Ready for use</td> </tr> <tr> <td colspan="4">Stop solution: 0.2M H₂SO₄</td> </tr> <tr> <td>CHROM</td> <td>TMB</td> <td>2 vials 13 ml</td> <td>Ready for use</td> </tr> <tr> <td colspan="4">Chromogenic TMB Solution (Tetramethylbenzidine)</td> </tr> <tr> <td>NEUTR</td> <td>SOLN</td> <td>1 vial 30 ml</td> <td>Ready for use</td> </tr> <tr> <td colspan="4">Neutralization Solution containing phosphate buffer with bovine casein</td> </tr> </table>	CONJ	BUF	1 vial 13 ml	Ready for use	Conjugate buffer: Phosphate buffer with caseine and thymol				CAL	0	1 vial 1 ml	Ready for use	Zero calibrator in phosphate buffer with ovalbumin and gentamycin				CAL	N	5 vials lyophilized	Add 1 ml reconstitution solution	Calibrators - N = 1 to 5 in phosphate buffer with ovalbumin and gentamycin. (see exact values on vial labels)				STOP	SOLN	1 vial 13 ml	Ready for use	Stop solution: 0.2M H ₂ SO ₄				CHROM	TMB	2 vials 13 ml	Ready for use	Chromogenic TMB Solution (Tetramethylbenzidine)				NEUTR	SOLN	1 vial 30 ml	Ready for use	Neutralization Solution containing phosphate buffer with bovine casein			
CONJ	BUF	1 vial 11.5 ml	Ready for use																																																																																														
Conjugate buffer: Phosphate buffer with caseine and thymol																																																																																																	
CAL	0	1 vial lyophilized	Add 1 ml reconstitution solution																																																																																														
Zero calibrator in phosphate buffer with bovine casein																																																																																																	
CAL	N	5 vials lyophilized	Add 1 ml reconstitution solution																																																																																														
Calibrators - N = 1 to 5 (see exact values on vial labels) in phosphate buffer with bovine casein																																																																																																	
STOP	SOLN	1 vial 12 ml	Ready for use																																																																																														
Stop solution: HCl 1M																																																																																																	
CHROM	TMB	1 vial 25 ml	Ready for use																																																																																														
Chromogenic TMB Solution (Tetramethylbenzidine)																																																																																																	
NEUTR	SOLN	1 vial 30 ml	Ready for use																																																																																														
Neutralization Solution containing phosphate buffer with bovine casein and azide (<0.1%)																																																																																																	
CONJ	BUF	1 vial 13 ml	Ready for use																																																																																														
Conjugate buffer: Phosphate buffer with caseine and thymol																																																																																																	
CAL	0	1 vial 1 ml	Ready for use																																																																																														
Zero calibrator in phosphate buffer with ovalbumin and gentamycin																																																																																																	
CAL	N	5 vials lyophilized	Add 1 ml reconstitution solution																																																																																														
Calibrators - N = 1 to 5 in phosphate buffer with ovalbumin and gentamycin. (see exact values on vial labels)																																																																																																	
STOP	SOLN	1 vial 13 ml	Ready for use																																																																																														
Stop solution: 0.2M H ₂ SO ₄																																																																																																	
CHROM	TMB	2 vials 13 ml	Ready for use																																																																																														
Chromogenic TMB Solution (Tetramethylbenzidine)																																																																																																	
NEUTR	SOLN	1 vial 30 ml	Ready for use																																																																																														
Neutralization Solution containing phosphate buffer with bovine casein																																																																																																	
XII. TYPICAL DATA	XII. TYPICAL DATA																																																																																																
The following data are for illustration only and should never be used instead of the real time calibration curve.	The following data are for illustration only and should never be used instead of the real time calibration curve.																																																																																																
<table border="1"> <thead> <tr> <th>Calibrators</th> <th>ng/ml</th> <th>OD units</th> </tr> </thead> <tbody> <tr> <td>S0</td> <td>0.0</td> <td>2.704</td> </tr> <tr> <td>S1</td> <td>15.0</td> <td>2.278</td> </tr> <tr> <td>S2</td> <td>23.8</td> <td>1.944</td> </tr> <tr> <td>S3</td> <td>59.0</td> <td>1.714</td> </tr> <tr> <td>S4</td> <td>183.0</td> <td>1.104</td> </tr> <tr> <td>S5</td> <td>774.0</td> <td>0.504</td> </tr> </tbody> </table>	Calibrators	ng/ml	OD units	S0	0.0	2.704	S1	15.0	2.278	S2	23.8	1.944	S3	59.0	1.714	S4	183.0	1.104	S5	774.0	0.504	<table border="1"> <thead> <tr> <th>Calibrators</th> <th>ng/ml</th> <th>OD units</th> </tr> </thead> <tbody> <tr> <td>S0</td> <td>0.0</td> <td>2.841</td> </tr> <tr> <td>S1</td> <td>18.0</td> <td>2.560</td> </tr> <tr> <td>S2</td> <td>28.6</td> <td>2.269</td> </tr> <tr> <td>S3</td> <td>56.0</td> <td>1.820</td> </tr> <tr> <td>S4</td> <td>154.0</td> <td>1.091</td> </tr> <tr> <td>S5</td> <td>643.0</td> <td>0.494</td> </tr> </tbody> </table>	Calibrators	ng/ml	OD units	S0	0.0	2.841	S1	18.0	2.560	S2	28.6	2.269	S3	56.0	1.820	S4	154.0	1.091	S5	643.0	0.494																																																						
Calibrators	ng/ml	OD units																																																																																															
S0	0.0	2.704																																																																																															
S1	15.0	2.278																																																																																															
S2	23.8	1.944																																																																																															
S3	59.0	1.714																																																																																															
S4	183.0	1.104																																																																																															
S5	774.0	0.504																																																																																															
Calibrators	ng/ml	OD units																																																																																															
S0	0.0	2.841																																																																																															
S1	18.0	2.560																																																																																															
S2	28.6	2.269																																																																																															
S3	56.0	1.820																																																																																															
S4	154.0	1.091																																																																																															
S5	643.0	0.494																																																																																															
XIII. PERFORMANCE AND LIMITATIONS	XIII. PERFORMANCE AND LIMITATIONS																																																																																																
A. Detection Limit	A. Detection Limit																																																																																																
The LoB (Limit of Blank) was calculated by measuring the blank 20 times and was calculated as the mean - 1.65 standard deviations of the distribution of these values. The LoB was calculated to be 4.5 ng/ml.	The LoB (Limit of Blank) was calculated by measuring the blank 20 times and was calculated as the mean - 1.65 standard deviations of the distribution of these values. The LoB was calculated to be 4.0 ng/ml.																																																																																																
B. Specificity	B. Specificity																																																																																																
B2. Interferences	B2. Interferences																																																																																																
The effect of potential interfering substances on samples using the DIAsource IGF-1-ELISA kit was evaluated. Different levels of Haemoglobin, Bilirubin, Triglyceride were tested on samples with different IGF-1 concentrations. Our acceptance criteria was to have interference of less than 10%. Hemolyzed and icteric samples have to be avoided. However, there is no restriction for the measurement of lipemic samples.	The effect of potential interfering substances on samples using the DIAsource IGF-1-ELISA kit was evaluated. Different levels of Haemoglobin, Bilirubin, Triglyceride were tested on samples with different IGF-1 concentrations. Our acceptance criteria was to have interference of less than 10%. Hemolyzed and icteric samples have to be avoided. However, there is no restriction for the measurement of lipemic samples.																																																																																																

Compound	IGF-1 concentration before addition of interferent (ng/ml)	Interferent added concentration (mg/dl)	IGF-1 concentration after addition of interferent	Mean variation (%)
Haemoglobin	69.8	500	58.5	-9.9
		250	60.1	
	99.9	500	100.8	
		250	89.4	
Bilirubin	69.8	100	56.1	-13.0
		50	64.1	
	99.9	100	83.3	
		50	92.1	
Triglyceride	45.6	2500	47.9	4.4
		1250	45.8	
	64.6	2500	68.2	
		1250	68.8	

C. Precision

Intra assay				Inter assay			
Sample	n	$\bar{x} \pm \text{sd}$ ng/ml	cv %	Sample	n	$\bar{x} \pm \text{sd}$ ng/ml	cv %
A	10	51,1 ± 3,1	6,1	C	15	58,4 ± 7,2	12,3
B	10	479,2 ± 42,5	8,9	D	15	127,6 ± 13,8	10,8

D. Accuracy

Sample (ng)	Added IGF-1 (ng)	Expected (ng/ml)	Measured (ng/ml)	Recovery (%)
16	15	31	32	103,2
16	30	46	48	104,3
16	100	116	120	103,4
16	250	266	256	96,2
16	400	416	431	103,6
16	600	616	681	110,6

Sample (ng)	Added IGF-1 (ng)	Expected (ng/ml)	Measured (ng/ml)	Recovery (%)
17	15	32	32	100
17	30	47	45	95,7
17	100	117	128	109,4
17	250	267	288	107,9
17	400	417	508	121,8
17	600	617	785	127,2

Sample A			Sample B		
Dil	Expected ng/ml	Measured ng/ml	Dil	Expected ng/ml	Measured ng/ml
1/1	540,0	540,0	1/1	875,0	875,0
1/2	270,0	266,0	1/2	438,0	440,0
1/4	135,0	115,0	1/4	219,0	183,0
1/8	68,0	62,0	1/8	109,0	79,0
1/16	34,0	37,0	1/16	55,0	46,1
1/32	17,0	25,3	1/32	27,0	31,1

Compound	IGF-1 concentration before addition of interferent (ng/ml)	Interferent added concentration (mg/dl)	IGF-1 concentration after addition of interferent	Mean variation (%)
Haemoglobin	320.2	500	316.0	8.4
		250	274.2	
	100.7	500	85.8	
		250	103.9	
Bilirubin	320.2	100	347.6	7.2
		50	316.3	
	100.7	100	105.2	
		50	115.3	
Triglyceride	320.2	2500	392.0	7.0
		1250	320.0	
	100.7	2500	103.9	
		1250	98.6	

C. Precision

Intra assay				Inter assay			
Sample	n	$\bar{x} \pm \text{sd}$ ng/ml	cv %	Sample	n	$\bar{x} \pm \text{sd}$ ng/ml	cv %
A	20	159,2 ± 12,3	7,7	C	10	124,4 ± 7,5	6,1
B	20	140,6 ± 9,7	6,9	D	10	187,4 ± 14,0	7,5

D. Accuracy

Sample (ng)	Added IGF-1 (ng)	Expected (ng/ml)	Measured (ng/ml)	Recovery (%)
1	15	45	51	112.5
1	30	60	60	99.5
1	100	151	130	86.4
1	250	280	290	103.5

Sample (ng)	Added IGF-1 (ng)	Expected (ng/ml)	Measured (ng/ml)	Recovery (%)
2	15	47	49	106
2	30	62	58	94.0
2	100	150	123	82.6
2	250	282	281	99.8

Sample A			Sample B		
Dil	Expected ng/ml	Measured ng/ml	Dil	Expected ng/ml	Measured ng/ml
1/1	478.8	478.8	1/1	490.4	490.4
1/2	239.4	214.0	1/2	245.2	184.0
1/4	119.7	114.5	1/4	122.6	108.3
1/8	59.8	59.2	1/8	61.3	55.4
1/16	29.9	34.7	1/16	30.7	29.7
1/32	15.0	18.6	1/32	15.3	18.3

Read entire protocol before use.

IGF1-ELISA

I. INTENDED USE

Immunoenzymetric assay for the *in vitro* quantitative measurement of human Insulin-like Growth Factor-I (IGF-I) in serum.

II. GENERAL INFORMATION

- A. **Proprietary name:** DIAsource IGF1-ELISA Kit
- B. **Catalogue number:** KAP1581: 96 tests
- C. **Manufactured by:** DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact:
Tel: +32 (0)10 84.99.11 Fax: +32 (0)10 84.99.91

III. CLINICAL BACKGROUND

A. Biological activities

Insulin-like growth factor I (IGF-I) or Somatomedin-C (SM-C) is a basic 70 amino acid single chain polypeptide (MW: 7649 Da) similar to proinsulin (50% sequence homology), and to the other well-characterized member of the somatomedin family: IGF II (67AA, 70 % sequence homology with IGF-I). IGF-I is the most important factor, which mediates the growth promoting actions of growth hormone, a pituitary hormone with highly fluctuating blood levels due to pulsatile release. The blood concentration of IGF-I is more stable due to the binding to carrier proteins. The concentration of the predominant binding protein (MW 53000) as well as the production of IGF-I, are regulated by growth hormone. IGF-I is produced by the liver, and other tissues, and it has endocrine, paracrine and autocrine activities. It stimulates growth and regulates differentiation of various tissues, displays insulin-like activities and promotes cartilage growth. Although GH is the most important factor controlling IGF-I secretion and concentration, other factors are also determinant: the age (with a peak at adolescence), the sex, the nutritional status, and other hormones (oestrogen, thyroxin, prolactin, ...). Specific trophic stimuli mainly control IGF-I secretion in the local microenvironment of a particular organ (paracrine activities), while blood IGF-I concentration is the most important variable for balanced systemic growth (endocrine activities).

B. Clinical applications

Growth retardation: Growth retardation may be due to several causes, among which deficient GH production (hypopituitarism), which is associated with low IGF-I blood levels. Because of the difficulties to get interpretable results from GH measurements (by dynamic multiple or stimulation tests), the determination of the stable IGF-I concentration in plasma is often considered as a simple screening test to evaluation "GH impregnation" of the patient before deciding more extensive investigations. In several clinical situations with impaired growth, low IGF-I levels may be observed despite normal or high GH production (i.e. malnutrition, chronic diseases states, some genetic dwarfs like Pygmies, ...). Interestingly, children with discrete GH neuro-secretory dysfunction may display low IGF I values despite normal GH levels by conventional testing. The results of IGF-I assay must be interpreted cautiously by considering the normal variations of IGF-I during childhood and adolescence (see Rosenfeld et al).

Acromegaly: IGF-I levels are elevated in acromegaly (excess production of GH) and may serve as an indicator of disease severity. Results are more readily interpreted because the normal values are more easily defined in adults. IGF-I measurements are also useful to monitor treatment.

Research: The IGF-I ELISA kit is an invaluable tool to study the modifications of this growth factor during physiologic (i.e. pregnancy) or pathologic (i.e. diabetes) situations, and the local regulation of IGF-I production in relation to its paracrine and autocrine activities (wound healing, organ regeneration, neoplastic growth, foetal development, gonadal regulation, etc).

IV. PRINCIPLES OF THE METHOD

The DIIAsource IGF1-ELISA is a solid phase Enzyme Amplified Sensitivity Immunoassay performed on polystyrene microtiterplates.

In the present kit, DIIAsource has introduced a pre-treatment step in order to improve the clinical performance of the assay. It is well established that the binding proteins interfere with the ELISA assay for IGF1. The pre-treatment step used by DIIAsource is the acid-ethanol procedure of Daughaday et al.(8)

A fixed amount of IGF1-labelled with horseradish peroxidase (HRP), compete with unlabelled IGF1 present in the calibrators, controls and samples for a limited number of binding sites on a specific antibody. After 1 hour incubation at room temperature, the microtiterplate is washed to stop the competition reaction. The Chromogenic solution (TMB) is added and incubated for 15 min. The reaction is stopped with the addition of Stop Solution and the microtiterplate is then read at the appropriate wavelength. The amount of substrate turnover is determined colourimetrically by measuring the absorbance, which is inversely proportional to the IGF1 concentration.

A calibration curve is plotted and IGF1 concentration in samples is determined by interpolation from the calibration curve.

V. REAGENTS PROVIDED

Reagents	96 tests Kit	Reconstitution			
 Microtiterplate with 96 anti IGF1 coated wells	96 wells	Ready for use			
<table border="1" data-bbox="60 801 255 840"> <tr> <td>Ag</td> <td>HRP</td> <td>CONC</td> </tr> </table> Conjugate: HRP labelled IGF1 in phosphate buffer with caseine and thymol	Ag	HRP	CONC	1 vial 0.2 ml	Dilute 100 x with conjugate buffer
Ag	HRP	CONC			
<table border="1" data-bbox="60 929 191 967"> <tr> <td>CONJ</td> <td>BUF</td> </tr> </table> Conjugate buffer: Phosphate buffer with caseine and thymol	CONJ	BUF	1 vial 13 ml	Ready for use	
CONJ	BUF				
<table border="1" data-bbox="60 1034 191 1072"> <tr> <td>PRE</td> <td>SOLN</td> </tr> </table> Pre-treatment Solution containing ethanol	PRE	SOLN	1 vial 20 ml	Ready for use	
PRE	SOLN				
<table border="1" data-bbox="60 1140 191 1178"> <tr> <td>REC</td> <td>SOLN</td> </tr> </table> Reconstitution Solution containing ethanol	REC	SOLN	1 vial 10 ml	Ready for use	
REC	SOLN				
<table border="1" data-bbox="60 1245 231 1283"> <tr> <td>NEUTR</td> <td>SOLN</td> </tr> </table> Neutralization Solution containing phosphate buffer with bovine casein	NEUTR	SOLN	1 vial 30 ml	Ready for use	
NEUTR	SOLN				
<table border="1" data-bbox="60 1350 231 1388"> <tr> <td>CAL</td> <td>0</td> </tr> </table> Zero calibrator in phosphate buffer with ovalbumin and gentamycin	CAL	0	1 vial 1 ml	Ready for use	
CAL	0				
<table border="1" data-bbox="60 1456 231 1494"> <tr> <td>CAL</td> <td>N</td> </tr> </table> Calibrators - N = 1 to 5 in phosphate buffer with ovalbumin and gentamycin. (see exact values on vial labels)	CAL	N	5 vials lyophilized	Add 1 ml reconstitution solution	
CAL	N				
<table border="1" data-bbox="60 1606 215 1644"> <tr> <td>CONTROL</td> <td>N</td> </tr> </table> Controls - N = 1 or 2 in human plasma with thymol	CONTROL	N	2 vials lyophilized	Add 1 ml distilled water	
CONTROL	N				
<table border="1" data-bbox="60 1711 279 1749"> <tr> <td>WASH</td> <td>SOLN</td> <td>CONC</td> </tr> </table> Wash Solution (Tris-HCl)	WASH	SOLN	CONC	1 vial 10 ml	Dilute 200 x with distilled water (use a magnetic stirrer).
WASH	SOLN	CONC			
<table border="1" data-bbox="60 1783 207 1821"> <tr> <td>CHROM</td> <td>TMB</td> </tr> </table> Chromogenic TMB Solution (Tetramethylbenzidine)	CHROM	TMB	2 vials 13 ml	Ready for use	
CHROM	TMB				
<table border="1" data-bbox="60 1888 199 1926"> <tr> <td>STOP</td> <td>SOLN</td> </tr> </table> Stop solution: 0.2M H ₂ SO ₄	STOP	SOLN	1 vial 13 ml	Ready for use	
STOP	SOLN				

- Note:** 1. Use the zero calibrator for sample dilutions.
2. 1 ng of the calibrator preparation is equivalent to 1 ng of the NIBSC 1stWHO IS 02/254.

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

- High quality distilled water
- Pipettes for delivery of: 50 µl, 100 µl, 400 µl, 500 µl, 600 µl and 1 ml (the use of accurate pipettes with disposable plastic tips is recommended)
- Optional: microcentrifuge tubes in polypropylene for pre-treatment of samples (for example: Eppendorf 3810)
- Polypropylene tubes for neutralization step (for example: Falcon 352063)
- Vortex mixer
- Magnetic stirrer
- Microcentrifuge or Centrifuge
- Washer for microtiterplates
- Horizontal microtiterplate shaker capable of 500 rpm
- Microtiterplate reader capable of reading at 450 nm and 650 nm (monochromatic reading)

VII. REAGENT PREPARATION

A. Calibrators: Reconstitute the calibrators with 1 ml reconstitution solution.

B. Controls: Reconstitute the controls with 1 ml distilled water.

C. Working IGF1-HRP conjugate: Prepare an adequate volume of conjugate solution by adding for example: 20 µl of the 100 x concentrated IGF1-HRP conjugate to 2 ml of conjugate buffer. Use a vortex to homogenize. Extemporaneous preparation is recommended.

D. Working Wash solution: Prepare an adequate volume of Working Wash solution by adding 199 volumes of distilled water to 1 volume of Wash Solution (200x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kit components are stable until the expiry date, indicated on the vial label, if kept at 2 to 8°C.
- Unused strips must be stored, at 2-8°C, in a sealed bag containing a desiccant until expiration date.
- After reconstitution, calibrators and controls are stable for one week at 2 to 8°C. For longer storage periods, aliquots should be made and kept at -20°C for maximum 3 months. Avoid successive freezing and thawing.
- The concentrated Wash Solution is stable at room temperature until expiration date.
- Freshly prepared Working Wash solution should be used on the same day.
- The Working IGF1-HRP conjugate is stable for 4 hours at room temperature, avoid direct sunlight.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- Serum must be kept at 2 - 8°C.
- If the test is not run within 24 hours, storage in aliquots at -20°C is recommended. Avoid subsequent freeze thaw cycles.
- Prior to use, all samples should be at room temperature. It is recommended to vortex the samples before use.
- Do not use haemolysed samples.

X. PROCEDURE

A. Handling notes

Do not use the kit or components beyond expiry date.

Do not mix materials from different kit lots.

Bring all the reagents to room temperature prior to use.

Thoroughly mix all reagents and samples by gentle agitation or swirling. Perform calibrators, controls and samples in duplicate. Vertical alignment is recommended.

Use a clean plastic container to prepare the Wash Solution.

In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.

For the dispensing of the Chromogenic Solution and the Stop Solution avoid pipettes with metal parts.

High precision pipettes or automated pipetting equipment will improve the precision.

Respect the incubation times.

To avoid drift, the time between pipetting of the first calibrator and the last sample must be limited to the time mentioned in section XIII paragraph E (Time delay).

Prepare a calibration curve for each run, do not use data from previous runs.

Dispense the Chromogenic Solution within 15 minutes following the washing of the microtiterplate.

During incubation with Chromogenic Solution, avoid direct sunlight on the microtiterplate.

B.1. Pre-treatment step for microfuge tubes

! Do not extract calibrators

1. Label one microcentrifuge tube (for extraction) and one polypropylene tube (for neutralization) for each sample and control. Dispense 100 µl of each sample and control into microfuge tube.
2. Add 400 µl of pre-treatment solution into this tube.
3. Close the tubes, vortex and incubate 30 minutes at room temperature.
4. Centrifuge for 2 minutes at 10000 rpm.
5. Take 100 µl of the supernatant and transfer it into the polypropylene labelled tube.
6. Add 600 µl of the neutralisation solution to this tube.
7. Vortex each tube.

B.2. Pre-treatment step for other polypropylene tubes

! Do not extract calibrators

1. Label two polypropylene tubes (one for extraction and the other for neutralization) for each sample and control.
2. Dispense 100 µl of each sample and control into extraction tube.
3. Add 400 µl of pre-treatment solution into this tube.
4. Close the tubes, vortex and incubate 30 minutes at room temperature.
5. Centrifuge for 15 minutes at > or = 3000 rpm.
6. Take 100 µl of the supernatant and transfer it into the polypropylene labelled tube.
7. Add 600 µl of the neutralisation solution to this tube.
8. Vortex each tube.

C. Procedure

1. Select the required number of wells for the run. The unused wells should be resealed in the bag with a desiccant and stored at 2-8°C.
2. Secure the wells into the holding frame.
3. Pipette 50 µl of each Calibrator, extracted control and extracted sample into the appropriate wells.
4. Pipette 100 µl of IGF1-HRP conjugate solution into all the wells.
5. Incubate for 1 hour at room temperature on a shaker at 500 rpm.
6. Aspirate the liquid from each well.
7. Wash the plate 3 times by:
 - dispensing 0.4 ml of Wash Solution into each well
 - aspirating the content of each well
8. Pipette 200 µl of the Chromogenic solution into each well within 15 minutes following the washing step.
9. Incubate the microtiterplate for 15 minutes at room temperature on the shaker (500 rpm), avoid direct sunlight
10. Pipette 100 µl of Stop Solution into each well.
11. Read the absorbances at 450 nm (reference filter 630 nm or 650 nm) within 1 hour and calculate the results as described in section XI

XI. CALCULATION OF RESULTS

1. Read the plate at 450 nm against a reference filter set at 650 nm (or 630 nm).
2. Calculate the mean of duplicate determinations.
3. Calculate for each calibrator, control and sample:

$$B/B0 (\%) = \frac{OD (\text{Calibrator, Control or Sample})}{OD (\text{Zero Calibrator})} \times 100$$

4. Using either linear-linear or semi-logarithmic graph paper, plot the (B/B0(%)) values for each calibrator point as a function of the IGF1 concentration of each calibrator point. Reject obvious outliers.
5. Computer assisted methods can also be used to construct the calibration curve. If automatic result processing is used, a 4-parameter logistic function curve fitting is recommended.

By interpolation of the sample (B/B0 (%)) values, determine the IGF1 concentrations of the samples from the calibration curve.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

Calibrators	ng/ml	OD units
S0	0.0	2.841
S1	18.0	2.560
S2	28.6	2.269
S3	56.0	1.820
S4	154.0	1.091
S5	643.0	0.494

XIII. PERFORMANCE AND LIMITATIONS

A. Detection Limit

The LoB (Limit of Blank) was calculated by measuring the blank 20 times and was calculated as the mean - 1.65 standard deviations of the distribution of these values. The LoB was calculated to be 4.0 ng/ml.

B. Specificity

B1. Cross-reaction

The percentages of cross-reaction estimated by comparison of the concentration yielding a 50% inhibition are respectively:

Compound	Cross reactivity %
IGF-II	0.09
Insuline	ND
GH	ND

B2. Interferences

The effect of potential interfering substances on samples using the DIAsource IGF-1-ELISA kit was evaluated. Different levels of Haemoglobin, Bilirubin, Triglyceride were tested on samples with different IGF-1 concentrations. Our acceptance criteria was to have interference of less than 10%.

Hemolyzed and icteric samples have to be avoided. However, there is no restriction for the measurement of lipemic samples.

Compound	IGF-1 concentration before addition of interferent (ng/ml)	Interferant added concentration (mg/dl)	IGF-1 concentration after addition of interferent	Mean variation (%)
Haemoglobin	320.2	500	316.0	8.4
		250	274.2	
		500	85.8	
Bilirubin	100.7	250	103.9	7.2
		100	347.6	
		50	316.3	
Triglyceride	320.2	100	105.2	7.0
		50	115.3	
		2500	392.0	
	100.7	1250	320.0	7.0
		2500	103.9	
		1250	98.6	

C. Precision

Intra assay				Inter assay			
Sample	n	x ± sd ng/ml	cv %	Sample	n	x ± sd ng/ml	cv %
A	20	159.2 ± 12.3	7.7	C	10	124.4 ± 7.5	6.1
B	20	140.6 ± 9.7	6.9	D	10	187.4 ± 14.0	7.5

D. Accuracy

RECOVERY TEST

Sample (ng)	Added IGF-1 (ng)	Expected (ng/ml)	Measured (ng/ml)	Recovery (%)
1	15	45	51	112.5
1	30	60	60	99.5
1	100	151	130	86.4
1	250	280	290	103.5

Sample (ng)	Added IGF-1 (ng)	Expected (ng/ml)	Measured (ng/ml)	Recovery (%)
2	15	47	49	106
2	30	62	58	94.0
2	100	150	123	82.6
2	250	282	281	99.8

DILUTION TEST

Sample A			Sample B		
Dil	Expected ng/ml	Measured ng/ml	Dil	Expected ng/ml	Measured ng/ml
1/1	478.8	478.8	1/1	490.4	490.4
1/2	239.4	214.0	1/2	245.2	184.0
1/4	119.7	114.5	1/4	122.6	108.3
1/8	59.8	59.2	1/8	61.3	55.4
1/16	29.9	34.7	1/16	30.7	29.7
1/32	15.0	18.6	1/32	15.3	18.3

Samples were diluted after extraction with zero calibrator

Conversion factor:

From ng/ml to nmol/L: x 0.13
From nmol/L to ng/ml: x 7.65

E. Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 30 minutes after the calibrators have been added to the coated wells.

Time delay	
0 minute	30 minutes
994,0 ng/ml	1057,0 ng/ml
458,0 ng/ml	476,0 ng/ml

XIV. INTERNAL QUALITY CONTROL


- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Controls which contain azide will interfere with the enzymatic reaction and cannot be used.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises.
- It is recommended that Controls be routinely assayed as unknown samples to measure assay variability. The performance of the assay should be monitored with quality control charts of the controls.
- It is good practise to check visually the curve fit selected by the computer.

XV. REFERENCE INTERVALS

These values are given only for guidance; each laboratory should establish its own normal range of values.

Age group (years)	Males (ng/ml)			Females (ng/ml)		
	n	Median	Range	n	Median	Range
0 - 5	30	37	14 - 154	50	57	21 - 262
6 - 8	20	93	45 - 213	18	207	89 - 485
9 - 11	20	147	53 - 453	19	354	99 - 708
12 - 15	24	313	103 - 753	31	274	100 - 744
16 - 20	30	203	99 - 655	29	248	73 - 522
21 - 25	9	167	115 - 304	9	171	83 - 511
26 - 39	24	162	87 - 415	20	165	101 - 267
40 - 54	35	156	69 - 343	17	121	60 - 271
> 55	25	117	33 - 232	14	109	69 - 189

These ranges are expressed as 2.5% to 97.5% percentiles.

XVI. PRECAUTIONS AND WARNINGS

Safety

For *in vitro* diagnostic use only.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections.

Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with all reagents, Stop Solution contains H₂SO₄. In case of contact, wash thoroughly with water.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVII. BIBLIOGRAPHY

1. DAUGHADAY W.H. and ROTWEIN P. (1989)
Insulin-like growth factors I and II. peptide, messenger ribonucleic acid and gene structures, serum and tissue concentrations.
Endocrine Rev., 10(1): 68-91.
2. FROESCH E.R. and ZAPF J. (1985)
Insulin-like growth factors and insulin: comparative aspects.
Diabetologia, 28: 485-493.
3. FURLANETTO R.W., UNDERWOOD L., VAN WYK J.J., D'ERCOLE A.J. (1977)
Estimation of somatomedin-C levels in normals and patients with pituitary disease by radioimmunoassay.
J. Clin. Invest. 60: 648.
4. ROSENFELD R.G., WILSON D.M., LEE P.D.K. and HINTZ R.L. (1986)
Insulin-like growth factors I and II in evaluation of growth retardation.
J. Pediatr., 109: 428.
5. RUDMAN D., KUTNER M.H. and CHAWLA R.K. (1985)
The short child with subnormal plasma somatomedin-C. Pediatric Res., 19(10): 975-980.
6. ALBERTSSON-WIKLAND K. and HALL K. (1987)
Growth hormone treatment in short children: relationship between growth and serum insulin-like growth factor I and II levels.
J. Clin. Endocrinol Metab., 65: 671.
7. WASS J.A.H., CLEMMONS D.R., UNDERWOOD L.E., BARROW L., BESSER G.M. and VAN WYK J.J. (1982)
Changes in circulating somatomedin-C levels in bromocriptine treated acromegaly.
Clin. Endocrinol., 17: 369-377.
8. DAUGHADAY W.H., MARIZ I.K. and BLETHEN S.L. (1980)
Inhibition of access of bound somatomedin to membrane receptor and immunobinding sites - a comparison of radioreceptor and radioimmunoassay of somatomedin in native and acid-ethanol extracted serum.
J. Clin. Endocrinol. Metab., 51: 781.
9. BREIERB.H., GALLAHER B.W. and GLUCKMAN P.D. (1991)
Radioimmunoassay for insulin-like growth factor-I: solutions to some potential problems and pitfalls.
Journal of Endocrinology, 128: 347-357.
10. SCHOUTEN J.S.A.G. and al. (1993)
IGF1: a prognostic factor of knee osteoarthritis.
British J. Rheumatol., 32: 274-280.
11. GRONBAEK H., SKJAERBAEK C., NIELSEN B., FRYSTYK J., FOEGH M.L., FLYVBJERG A., ORSKOV H. (1995)
Growth hormone and insulin-like growth factor-I: a suggested role in renal transplantation and graft vessel disease.
Transplant. Proceed., 27/3: 2133-2136.
12. TSITOURAS P.D., ZHONG Y.G., SPUNGEN A.M., BAUMAN W.A. (1995)
Serum testosterone and growth hormone insulin-like growth factor-I in adults with spinal cord injury.
Hormone and metabolic Research, 27/6 : 287-292.
13. KOCH A., DORR H.G., GERLING S., BEHRENS R., BOHLES H.J. (1995)
Effect of growth hormone on IDF-I levels in patients with growth hormone deficiency and Wilson disease.
Hormone Research, 44/1: 40-44.

XVIII. SUMMARY OF THE PROTOCOL

	CALIBRATORS µl	SAMPLE(S) CONTROLS µl
PRE-TREATMENT		
Samples, controls	-	100
Pre-treatment solution	-	400
Incubation	30 minutes RT	
Centrifugation	2 minutes at 10000 rpm (microfuge tubes) or 15 minutes at > or = 3000 rpm (polypropylene tubes)	
Supernatant	-	100
Neutralization Solution	-	600
Shaking	Vortex	
INCUBATION		
Calibrators (0 to 5)	50	-
Pre-treated Samples, controls	-	50
Diluted Conjugate	100	100
Incubate for 1 hour at room temperature on a shaker at 500 rpm. Aspirate the contents of each well. Wash 3 times with 400 µl of Wash Solution and aspirate.		
Chromogenic Solution	200	200
Incubate for 15 min at room temperature on a shaker at 500 rpm		
Stop Solution	100	100
Read on a microtiterplate reader and record the absorbance of each well at 450 nm (versus 630 or 650 nm)		

Other translations of this Instruction for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

DIAsource Catalogue Nr: KAP1581	Revision nr: 220516
------------------------------------	------------------------

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Lapso de tiempo	
minuto 0	30 minutos
994,0 ng/ml	1057,0 ng/ml
458,0 ng/ml	476,0 ng/ml

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
- DAUGHADAY W.H. and ROTWEIN P. (1989) **Insulin-like growth factors I and II. peptide, messenger ribonucleic acid and gene structures, serum and tissue concentrations.** Endocrine Rev., 10(1): 68-91.
- FROESCH E.R. and ZAPF J. (1985) **Insulin-like growth factors and insulin: comparative aspects.** Diabetologia, 28: 485-493.

3. FURLANETTO R.W., UNDERWOOD L., VAN WYK J.J., D'ERCOLE A.J. (1977)
Estimation of somatomedin-C levels in normals and patients with pituitary disease by radioimmunoassay.
J. Clin. Invest. 60: 648.
4. ROSENFELD R.G., WILSON D.M., LEE P.D.K. and HINTZ R.L. (1986)
Insulin-like growth factors I and II in evaluation of growth retardation.
J. Pediatr., 109: 428.
5. RUDMAN D., KUTNER M.H. and CHAWLA R.K. (1985)
The short child with subnormal plasma somatomedin-C. Pediatric Res., 19(10): 975-980.
6. ALBERTSSON-WIKLAND K. and HALL K. (1987)
Growth hormone treatment in short children: relationship between growth and serum insulin-like growth factor I and II levels.
J. Clin. Endocrinol Metab., 65: 671.
7. WASS J.A.H., CLEMMONS D.R., UNDERWOOD L.E., BARROW L., BESSER G.M. and VAN WYK J.J. (1982)
Changes in circulating somatomedin-C levels in bromocriptine treated acromegaly.
Clin. Endocrinol., 17: 369-377.
8. DAUGHADAY W.H., MARIZ I.K. and BLETHEN S.L. (1980)
Inhibition of access of bound somatomedin to membrane receptor and immunobinding sites - a comparison of radioreceptor and radioimmunoassay of somatomedin in native and acid-ethanol extracted serum.
J. Clin. Endocrinol. Metab., 51: 781.
9. BREIERB.H., GALLAHER B.W. and GLUCKMAN P.D. (1991)
Radioimmunoassay for insulin-like growth factor-I: solutions to some potential problems and pitfalls.
Journal of Endocrinology, 128: 347-357.
10. SCHOUTEN J.S.A.G. and al. (1993)
IGF1: a prognostic factor of knee osteoarthritis.
British J. Rheumatol., 32: 274-280.
11. GRONBAEK H., SKJAERBAEK C., NIELSEN B., FRYSTYK J., FOEGH M.L., FLYVBJERG A., ORSKOV H. (1995)
Growth hormone and insulin-like growth factor-I: a suggested role in renal transplantation and graft vessel disease.
Transplant. Proceed., 27/3: 2133-2136.
12. TSITOURAS P.D., ZHONG Y.G., SPUNGEN A.M., BAUMAN W.A. (1995)
Serum testosterone and growth hormone insulin-like growth factor-I in adults with spinal cord injury.
Hormone and metabolic Research, 27/6: 287-292.
13. KOCH A., DORR H.G., GERLING S., BEHRENS R., BOHLES H.J. (1995)
Effect of growth hormone on IGF-I levels in patients with growth hormone deficiency and Wilson disease.
Hormone Research, 44/1: 40-44.

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Opóźnienie	
0 minut	30 minut
994,0 ng/ml	1057,0 ng/ml
458,0 ng/ml	476,0 ng/ml

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
XVII. BIBLIOGRAFIA

1. DAUGHADAY W.H. and ROTWEIN P. (1989)
Insulin-like growth factors I and II. peptide, messenger ribonucleic acid and gene structures, serum and tissue concentrations.
Endocrine Rev., 10(1): 68-91.
2. FROESCH E.R. and ZAPF J. (1985)
Insulin-like growth factors and insulin: comparative aspects.
Diabetologia, 28: 485-493.
3. FURLANETTO R.W., UNDERWOOD L., VAN WYK J.J., D'ERCOLE A.J. (1977)
Estimation of somatomedin-C levels in normals and patients with pituitary disease by radioimmunoassay.
J. Clin. Invest. 60: 648.
4. ROSENFELD R.G., WILSON D.M., LEE P.D.K. and HINTZ R.L. (1986)
Insulin-like growth factors I and II in evaluation of growth retardation.
J. Pediatr., 109: 428.
5. RUDMAN D., KUTNER M.H. and CHAWLA R.K. (1985)
The short child with subnormal plasma somatomedin-C. Pediatric Res., 19(10): 975-980.
6. ALBERTSSON-WIKLAND K. and HALL K. (1987)
Growth hormone treatment in short children: relationship between growth and serum insulin-like growth factor I and II levels.
J. Clin. Endocrinol Metab., 65: 671.
7. WASS J.A.H., CLEMMONS D.R., UNDERWOOD L.E., BARROW L., BESSER G.M. and VAN WYK J.J. (1982)
Changes in circulating somatomedin-C levels in bromocriptine treated acromegaly.
Clin. Endocrinol., 17: 369-377.
8. DAUGHADAY W.H., MARIZ I.K. and BLETHEN S.L. (1980)
Inhibition of access of bound somatomedin to membrane receptor and immunobinding sites - a comparison of radioreceptor and radioimmunoassay of somatomedin in native and acid-ethanol extracted serum.
J. Clin. Endocrinol. Metab., 51: 781.
9. BREIERB.H., GALLAHER B.W. and GLUCKMAN P.D. (1991)
Radioimmunoassay for insulin-like growth factor-I: solutions to some potential problems and pitfalls.
Journal of Endocrinology, 128: 347-357.
10. SCHOUTEN J.S.A.G. and al. (1993)
IGF1: a prognostic factor of knee osteoarthritis.
British J. Rheumatol., 32: 274-280.
11. GRONBAEK H., SKJAERBAEK C., NIELSEN B., FRYSTYK J., FOEGH M.L., FLYVBJERG A., ORSKOV H. (1995)
Growth hormone and insulin-like growth factor-I: a suggested role in renal transplantation and graft vessel disease.

Transplant. Proceed., 27/3: 2133-2136.

12. TSITOURAS P.D., ZHONG Y.G., SPUNGEN A.M., BAUMAN W.A. (1995)
Serum testosterone and growth hormone insulin-like growth factor-I in adults with spinal cord injury.
Hormone and metabolic Research, 27/6: 287-292.
13. KOCH A., DORR H.G., GERLING S., BEHRENS R., BOHLES H.J. (1995)
Effect of growth hormone on IGF-I levels in patients with growth hormone deficiency and Wilson disease.
Hormone Research, 44/1: 40-44.

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