17 OH-PROGESTERONE ELISA

KAPD1292

History

Summary of change :

| Previous Version : | Current Version : | | | | | | | | | | | | | | | | | | | | | | | | |
|--|--|------------------------|------------------------|------|---------------------------|------|--------------------------|-----|--------------------------|------|--------------------------|------|---|------------|------------------------|------------------------|------|---------------------------|------|--------------------------|------|--------------------------|------|--------------------------|------|
| 200831 | 210208 | | | | | | | | | | | | | | | | | | | | | | | | |
| <p>2 PRINCIPLE OF THE TEST The DIAsource 17-α-OH Progesterone ELISA Kit is a solid phase enzyme-linked immunosorbent assay (ELISA), based on the principle of competitive binding. The microtiter wells are coated with a polyclonal antibody directed towards an antigenic site on the 17 α OHP molecule. Endog-enous 17-α-OHP of a patient sample competes with a 17-α-OHP-horseradish peroxidase conjugate for binding to the coated antibody. After incubation the unbound conjugate is washed off. The amount of bound peroxidase conjugate is inversely proportional to the concentration of 17-α-OHP in the sample. After addition of the substrate solution, the intensity of colour developed is inversely proportional to the concentration of 17-α-OHP in the patient sample.</p> | <p>2 PRINCIPLE OF THE TEST The DIAsource 17-OH Progesterone ELISA is a solid phase enzyme-linked immunosorbent assay (ELISA), based on the principle of competitive binding. The microtiter wells are coated with a polyclonal [rabbit] antibody directed towards antigenic sites of the 17 α OHP molecule. Samples are pre-incubated in the coated wells. During the second incubation, 17-α-OHP in the added sample competes with the added enzyme conjugate, which is 17 α OHP conjugated to horseradish peroxidase, for binding to the coated antibody. After a washing step to remove all unbound substances, the solid phase is incubated with the substrate solution. The colorimetric reaction is stopped by addition of stop solution, and optical density (OD) of the resulting yellow product is measured. The intensity of colour is inversely proportional to the concentration of the analyte in the sample. A standard curve is constructed by plotting OD values against concentrations of standards, and concentrations of unknown samples are determined using this standard curve.</p> | | | | | | | | | | | | | | | | | | | | | | | | |
| <p>4.1 Contents of the Kit No indication of calibration for the calibrators</p> | <p>4.1 Contents of the Kit The calibrators are calibrated against the following reference material: Certified Reference Material Cerilliant H-085</p> | | | | | | | | | | | | | | | | | | | | | | | | |
| <p>4.4 Preparation of Reagents The diluted Wash Solution is stable for 2 weeks at room temperature.</p> | <p>4.4 Preparation of Reagents The diluted Wash Solution is stable for 1 weeks at room temperature.</p> | | | | | | | | | | | | | | | | | | | | | | | | |
| <p>6.3 Calculation of Results</p> <table border="1"> <thead> <tr> <th>Calibrator</th> <th>Optical Units (450 nm)</th> </tr> </thead> <tbody> <tr> <td>Calibrator 0 (0 ng/mL)</td> <td>1,89</td> </tr> <tr> <td>Calibrator 1 (0,15 ng/mL)</td> <td>1,51</td> </tr> <tr> <td>Calibrator 2 (0,5 ng/mL)</td> <td>1,1</td> </tr> <tr> <td>Calibrator 3 (1,5 ng/mL)</td> <td>0,69</td> </tr> <tr> <td>Calibrator 4 (3,0 ng/mL)</td> <td>0,46</td> </tr> </tbody> </table> | Calibrator | Optical Units (450 nm) | Calibrator 0 (0 ng/mL) | 1,89 | Calibrator 1 (0,15 ng/mL) | 1,51 | Calibrator 2 (0,5 ng/mL) | 1,1 | Calibrator 3 (1,5 ng/mL) | 0,69 | Calibrator 4 (3,0 ng/mL) | 0,46 | <p>6.3 Calculation of Results</p> <table border="1"> <thead> <tr> <th>Calibrator</th> <th>Optical Units (450 nm)</th> </tr> </thead> <tbody> <tr> <td>Calibrator 0 (0 ng/mL)</td> <td>2.15</td> </tr> <tr> <td>Calibrator 1 (0,15 ng/mL)</td> <td>1.77</td> </tr> <tr> <td>Calibrator 2 (0,5 ng/mL)</td> <td>1.28</td> </tr> <tr> <td>Calibrator 3 (1,5 ng/mL)</td> <td>0.77</td> </tr> <tr> <td>Calibrator 4 (3,0 ng/mL)</td> <td>0.49</td> </tr> </tbody> </table> | Calibrator | Optical Units (450 nm) | Calibrator 0 (0 ng/mL) | 2.15 | Calibrator 1 (0,15 ng/mL) | 1.77 | Calibrator 2 (0,5 ng/mL) | 1.28 | Calibrator 3 (1,5 ng/mL) | 0.77 | Calibrator 4 (3,0 ng/mL) | 0.49 |
| Calibrator | Optical Units (450 nm) | | | | | | | | | | | | | | | | | | | | | | | | |
| Calibrator 0 (0 ng/mL) | 1,89 | | | | | | | | | | | | | | | | | | | | | | | | |
| Calibrator 1 (0,15 ng/mL) | 1,51 | | | | | | | | | | | | | | | | | | | | | | | | |
| Calibrator 2 (0,5 ng/mL) | 1,1 | | | | | | | | | | | | | | | | | | | | | | | | |
| Calibrator 3 (1,5 ng/mL) | 0,69 | | | | | | | | | | | | | | | | | | | | | | | | |
| Calibrator 4 (3,0 ng/mL) | 0,46 | | | | | | | | | | | | | | | | | | | | | | | | |
| Calibrator | Optical Units (450 nm) | | | | | | | | | | | | | | | | | | | | | | | | |
| Calibrator 0 (0 ng/mL) | 2.15 | | | | | | | | | | | | | | | | | | | | | | | | |
| Calibrator 1 (0,15 ng/mL) | 1.77 | | | | | | | | | | | | | | | | | | | | | | | | |
| Calibrator 2 (0,5 ng/mL) | 1.28 | | | | | | | | | | | | | | | | | | | | | | | | |
| Calibrator 3 (1,5 ng/mL) | 0.77 | | | | | | | | | | | | | | | | | | | | | | | | |
| Calibrator 4 (3,0 ng/mL) | 0.49 | | | | | | | | | | | | | | | | | | | | | | | | |

| | | | | | |
|-----------|------------------------------|------|-----------|---|------|
| | Calibrator 5 (7,5 ng/mL) | 0,28 | | Calibrator 5 (7,5 ng/mL) | 0.25 |
| | Calibrator 6 (20 ng/mL) | 0,18 | | Calibrator 6 (20 ng/mL) | 0.12 |
| 9 | ASSAY CHARACTERISTICS | | 9 | ASSAY CHARACTERISTICS Complete new data | |
| 12 | REFERENCES | | 12 | REFERENCES References updated | |

17- α -OH Progesterone-ELISA

en

KAPD1292

IN VITRO DIAGNOSTIC USE

DIAsource ImmunoAssays S.A. - Rue du Bosquet 2, B-1348 Louvain-la-Neuve, Belgium - Tel: +32 10 84 99 11 - Fax : +32 10 84 99 91

1 INTRODUCTION

The DIAsource 17-OH Progesterone Enzyme Immunoassay Kit provides materials for the quantitative determination of 17- α -OH Progesterone in serum or plasma (EDTA-, heparin- or citrate plasma). This assay is intended for in vitro diagnostic use only.

The steroid 17- α -Hydroxyprogesterone (17- α -OHP) is produced by both the adrenal cortex and gonads. Even though 17- α -OHP has relatively low progestational activity, it is of intense clinical interest because it is the immediate precursor to 11-desoxycortisol (Cpd-S). Because Cpd-S is produced by 21-hydroxylation of 17- α -OHP, measurement of 17- α -OHP is a useful indirect indicator of 21-hydroxylase activity. In congenital 21-hydroxylase deficiency, the most common variety of congenital adrenal hyperplasia (CAH), 17- α -OHP is secreted in abundant excess. It is moderately elevated in the 11- β -hydroxylase deficiency as well. Measurement of 17- α -OHP is therefore valuable in the initial diagnosis of CAH.

1.1 CLINICAL PHYSIOLOGY

Adult non-pregnant women:

In adult non-pregnant women in the childbearing age group, 17- α -OHP concentrations vary over the menstrual cycle with luteal phase concentrations being higher than follicular phase concentrations. This is because 17- α -OHP is secreted in parallel with progesterone from maturing follicles or from the corpus luteum. There is also a diurnal variation of 17- α -OHP concentrations.

This rhythm is parallel with adrenal cortisol secretion such that maximum 17- α -OHP concentrations are measured in samples obtained between midnight and 8:00 am.

Adult males:

There is little information available on the systematic variability of 17- α -OHP concentration in adult males.

Pregnant women and newborn children:

The steroid 17- α -OHP is produced in large amounts by the fetus and the adrenals. It is secreted in abundance into both the fetal and maternal circulation. The maternal concentrations of 17- α -OHP increase very sharply after 32 weeks gestational age to about 4-fold above basal concentrations at term.

1.2 CLINICAL APPLICATIONS

Congenital adrenal hyperplasia:

The principal application of the 17- α -OHP Elisa is in the diagnosis of CAH in newborns with ambiguous genitalia and in virilized adolescent girls. Since 17- α -OHP is the immediate precursor to 11-desoxycortisol, basal 17- α -OHP concentrations are sharply elevated in patients with 21-hydroxylase deficiency and to a lesser degree in patients with 11-hydroxylase deficiency.

Because 17- α -OHP concentrations are so markedly elevated in newborns and adolescent girls afflicted with CAH, a single basal measurement is all that is normally required to make the diagnosis.

Late onset adrenal hyperplasia:

More recently, 17- α -OHP concentrations have been utilized in the evaluation of androgenized women where late onset 21-hydroxylase deficiency is suspected. This condition is clinically very subtle and since the presentation is the same as classical polycystic ovarian disease, basal plasma 17- α -OHP concentrations, unlike classical congenital adrenal hyperplasia, are normal. The diagnosis is made by administration of an ACTH stimulation test.

Other applications:

Measurement of 17- α -OHP concentrations is also utilized in evaluation of both men and women with acne vulgaris, male pattern baldness and in some subtle forms of infertility. Experiences with these applications are very limited.

2 PRINCIPLE OF THE TEST

The DIAsource 17-OH Progesterone ELISA is a solid phase enzyme-linked immunosorbent assay (ELISA), based on the principle of competitive binding.

The microtiter wells are coated with a polyclonal [rabbit] antibody directed towards antigenic sites of the 17 α OHP molecule.

Samples are pre-incubated in the coated wells.

During the second incubation, 17- α -OHP in the added sample competes with the added enzyme conjugate, which is 17 α OHP conjugated to horseradish peroxidase, for binding to the coated antibody.

After a washing step to remove all unbound substances, the solid phase is incubated with the substrate solution. The colorimetric reaction is stopped by addition of stop solution, and optical density (OD) of the resulting yellow product is measured. The intensity of colour is inversely proportional to the concentration of the analyte in the sample.

A standard curve is constructed by plotting OD values against concentrations of standards, and concentrations of unknown samples are determined using this standard curve.

3 PRECAUTIONS

- This kit is for in vitro diagnostic use only. For professional use only.
- All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.

- Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.
- The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided.
- Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
- Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.
- Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
- Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
- Allow the reagents to reach room temperature (21°C to 26°C) before starting the test. Temperature will affect the absorbance readings of the assay. However, values for the patient samples will not be affected.
- Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
- Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
- Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
- Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
- Do not use reagents beyond expiry date as shown on the kit labels.
- All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiter plate readers.
- Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
- Avoid contact with *Stop Solution* containing 0.5 M H₂SO₄. It may cause skin irritation and burns.
- Some reagents contain Proclin, BND and MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
- TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
- Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
- For information on hazardous substances included in the kit please refer to Material Safety Data Sheets. Material Safety Data Sheets for this product are available upon request directly from DIAsource ImmunoAssays S.A.

4 KIT COMPONENTS

4.1 Contents of the Kit

- | |
|-----|
| 111 |
|-----|

 12x8 (break apart) strips, 96 wells
 Wells coated with Anti-17- α -OH Progesterone antibody (polyclonal)
- | | |
|-----|---|
| CAL | N |
|-----|---|

 N=1 to 6, Reference Calibrator Set, 6 vials, 1 ml each
 The calibrators are calibrated against the following reference material: Certified Reference Material Cerilliant H-085
 Contains non-mercury preservative
 See exact values on vial labels
 Conversion : ng/mL x 3.03 = nmol/L
- | | |
|-----|---|
| CAL | 0 |
|-----|---|

 Zero Calibrator, 1 vial, 1 ml
 Contains non-mercury preservative
- | | |
|---------|---|
| CONTROL | N |
|---------|---|

 2 vials, 1 mL each, ready to use;
 For control values and ranges please refer to vial label or QC-Datasheet
 Contains non-mercury preservative
- | | |
|----|-----|
| Ag | HRP |
|----|-----|

 1 vial, 25 ml, ready to use
 17- α -OH Progesterone conjugated to horseradish peroxidase
 Contains non-mercury preservative
- | | |
|-------|-----|
| CHROM | TMB |
|-------|-----|

 (Substrate Solution) 1 vial, 25ml, ready to use
 TMB
- | | |
|------|------|
| STOP | SOLN |
|------|------|

 1 vial, 14 ml
 contains 0.5M H₂SO₄
 Avoid contact with the stop solution. It may cause skin irritations and burns.
- | | | |
|------|------|------|
| WASH | SOLN | CONC |
|------|------|------|

 (40x concentrated), 1 vial, 30 ml
 See "Preparation of Reagents"

Note: Additional *Calibrator 0* for sample dilution is available upon request.

4.2 Equipment and material required but not provided

- A microtiter plate calibrated reader (450±10 nm)
- Calibrated variable precision micropipettes.
- Absorbent paper
- Aqua dest.
- Semi-logarithmic graph paper

4.3 Storage and stability of the Kit

When stored at 2-8°C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date. Opened reagents must be stored at 2-8°C. Microtiter wells must be stored at 2-8°C. Once the foil bag has been opened, care should be taken to close it tightly again. Opened kits retain activity for two months if stored as described above.

4.4 Preparation of Reagents

Allow all reagents and required number of strips to reach room temperature prior to use.

Wash Solution

Add deionized water to the 40 x concentrated Wash Solution.

Dilute 30 mL of concentrated Wash Solution with 1170 mL deionized water to a final volume of 1200 mL.

The diluted Wash Solution is stable for 1 weeks at room temperature.

4.5 Disposal of the Kit

The disposal of the kit must be made according to the national regulations. Special information for this product is given in the Material Safety Data Sheets.

4.6 Damaged Test Kits

In case of any severe damage to the test kit or components, DIAsource ImmunoAssays S.A. has to be informed in writing, at the latest, one week after receiving the kit. Severely damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed according to the official regulations.

5 SPECIMEN

Serum or plasma (EDTA-, heparin- or citrate plasma) can be used in this assay.

Do not use haemolytic, icteric or lipaemic specimens.

Please note: Samples containing sodium azide should not be used in the assay.

5.1 Specimen Collection

Serum:

Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

Plasma:

Whole blood should be collected into centrifuge tubes containing anti-coagulant (e.g. Sarstedt Monovette with the appropriate plasma preparation) and centrifuged immediately after collection.

5.2 Specimen Storage

Specimens should be capped and may be stored for up to 7 days at 2-8°C prior to assaying.

Specimens held for a longer time should be frozen only once at -20°C prior to assay. Thawed samples should be inverted several times prior to testing.

5.3 Specimen Dilution

If in an initial assay, a specimen is found to contain more than the highest calibrator, the specimens can be diluted with *Calibrator 0* and reassayed as described in Assay Procedure.

For the calculation of the concentrations this dilution factor has to be taken into account.

Example:

- a) Dilution 1:10: 10 µL Serum + 90 µL *Calibrator 0* (mix thoroughly)
- b) Dilution 1:100: 10 µL dilution a) 1:10 + 90 µL *Calibrator 0* (mix thoroughly).

6 TEST PROCEDURE

6.1 General Remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each calibrator, control or sample in order to avoid cross contamination.
- Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.

6.2 Assay Procedure

Each run must include a calibrator curve.

1. Secure the desired number of Microtiter wells in the holder.
2. Dispense **25 µL** of each *Calibrator*, *Control* and samples with new disposable tips into appropriate wells.
3. Incubate for **5 minutes** at room temperature
4. Dispense **200 µL Enzyme Conjugate** into each well.
5. Thoroughly mix for 10 seconds. It is important to have a complete mixing in this step.
6. Incubate for **60 minutes** at room temperature.
7. Briskly shake out the contents of the wells.

8. Rinse the wells 3 times with diluted Wash Solution (400 μL per well). Strike the wells sharply on absorbent paper to remove residual droplets.
Important note:
The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!
9. Add **200 μL** of *Substrate Solution* to each well.
10. Incubate for **30 minutes** at room temperature.
11. Stop the enzymatic reaction by adding **100 μL** of *Stop Solution* to each well.
12. Read the OD at **450 \pm 10 nm** with a microtiter plate reader **within 10 minutes** after adding the *Stop Solution*.

6.3 Calculation of Results

1. Calculate the average absorbance values for each set of calibrators, controls and patient samples.
2. Using semi-logarithmic graph paper, construct a calibrator curve by plotting the mean absorbance obtained from each calibrator against its concentration with absorbance value on the vertical(Y) axis and concentration on the horizontal (X) axis.
3. Using the mean absorbance value for each sample determine the corresponding concentration from the calibrator curve.
4. Automated method: The results in the Instructions For Use have been calculated automatically using a 4 Parameter curve fit. (4 Parameter Rodbard or 4 Parameter Marquardt are the preferred methods.) Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this calibrator curve. Samples with concentrations higher than that of the highest calibrator have to be further diluted. For the calculation of the concentrations this dilution factor has to be taken into account.

Below is listed a typical example of a calibrator curve with the 17- α -OH Progesterone ELISA.

| Calibrator | Optical Units (450 nm) |
|---------------------------|------------------------|
| Calibrator 0 (0 ng/mL) | 2.15 |
| Calibrator 1 (0,15 ng/mL) | 1.77 |
| Calibrator 2 (0,5 ng/mL) | 1.28 |
| Calibrator 3 (1,5 ng/mL) | 0.77 |
| Calibrator 4 (3,0 ng/mL) | 0.49 |
| Calibrator 5 (7,5 ng/mL) | 0.25 |
| Calibrator 6 (20 ng/mL) | 0.12 |

7 EXPECTED VALUES

It is strongly recommended that each laboratory should determine its own normal and abnormal values.

In a study conducted with newborns and children, using the DIAsource 17-OH Progesterone ELISA, the following values are observed:

| Newborns (boys and girls) | n | | Range (min - max) (ng/mL) | Mean (ng/mL) | Median (ng/mL) | 2.5 th - 97.5 th Percentile (ng/mL) |
|---------------------------|-----------------------------------|-----------------------------------|------------------------------|-----------------|-------------------|--|
| | 26 | 1 st month after birth | | 0 - 17.3 | 7.2 | 6.7 |
| 43 | 2 nd month after birth | | 0.32 - 13.7 | 4.9 | 4.6 | 1.6 - 9.8 |
| 21 | 3 rd month after birth | | 0.06 - 4.2 | 2.3 | 2.3 | 0.5 - 4.1 |
| 12 | 4 th month after birth | | 0.2 - 4.6 | 2.1 | 2.3 | 0.2 - 4.3 |

| | n | Age (years) | Range (min - max) (ng/mL) | Mean (ng/mL) | Median (ng/mL) | 2.5 th - 97.5 th Percentile (ng/mL) |
|-------------------|----|-------------|------------------------------|-----------------|-------------------|--|
| Children | 75 | 1 - 10 | 0.03 - 2.85 | 1.04 | 0.88 | 0.08 - 2.58 |
| Adolescent | 3 | 11 - 14 | 0.06 - 1.38 | 0.65 | 0.50 | 0.07 - 1.34 |
| | 10 | 15 - 18 | 0.41 - 2.35 | 1.24 | 1.26 | 0.42 - 2.26 |

In a study conducted with apparently normal healthy donors, using the DIAsource 17-OH Progesterone ELISA, the following values are observed:

| | | |
|----------------------|------------------|-------------------|
| Adult females | Follicular phase | 0.1 - 0.8 ng/mL |
| | Luteal phase | 0.6 - 2.3 ng/mL |
| | Ovulation | 0.3 - 1.4 ng/mL |
| | Post ACTH | < 3.2 ng/mL |
| | Third trimester | 2.0 - 12 ng/mL |
| | Postmenopausal | 0.13 - 0.51 ng/mL |
| Adult males | | 0.5 - 2.1 ng/mL |

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated to other clinical observations and diagnostic tests.

8 QUALITY CONTROL

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results. Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid.

In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

9 ASSAY CHARACTERISTICS

9.1 Assay Dynamic Range

The range of the assay is between 0.042 – 20 ng/mL.

9.2 Specificity of Antibodies (Cross Reactivity)

The following substances were tested for cross reactivity of the assay:

| Substance | Conc. Range (ng/mL) | Mean Bias (ng/mL) | Mean Bias % |
|------------------------|---------------------|-------------------|-------------|
| 17-Benzoate Estradiol | 2 - 2000 | 0.05 | < 0.01 |
| 17-Cypionate Estradiol | 2 - 2000 | < 0.01 | < 0.01 |
| 17-Valerate Estradiol | 2 - 2000 | < 0.01 | < 0.01 |
| Aldosterone | 2 - 2000 | 0.06 | < 0.01 |
| Androstenedione | 2 - 2000 | < 0.01 | < 0.01 |
| Corticosterone | 2 - 200 | < 0.01 | < 0.01 |
| Cortisol | 2 - 200 | < 0.01 | < 0.01 |
| Cortisone | 2 - 200 | < 0.01 | < 0.01 |
| DHEA | 2 - 2000 | < 0.01 | < 0.01 |
| DHEA-S(a) | 2 - 2000 | < 0.01 | < 0.01 |
| Estradiol | 2 - 2000 | < 0.01 | < 0.01 |
| Estriol | 2 - 2000 | < 0.01 | < 0.01 |
| Estrone | 2 - 2000 | < 0.01 | < 0.01 |
| Progesterone | 2 - 20 | 0.05 | 1.31 |
| Testosterone | 2 - 2000 | < 0.01 | < 0.01 |

9.3 Analytical Sensitivity

The analytical sensitivity of the Diasource ELISA was calculated by subtracting 2 standard deviations from the mean OD of 20 replicate analyses of Standard 0 and was found to be 0.013 ng/mL.

The Limit of Blank (LoB) is 0.006 ng/mL.

The Limit of Detection (LoD) is 0.042 ng/mL.

The Limit of Quantification (LoQ) is 0.156 ng/mL...

9.4 Precision

9.4.1 Intra Assay Variation

The within assay variability is shown below:

| Sample | n | Mean (ng/mL) | CV (%) |
|-----------------------|----|--------------|--------|
| 1 (Li-Heparin plasma) | 10 | 0.23 | 4.6 |
| 2 (Serum) | 10 | 2.10 | 4.0 |
| 3 (Serum) | 10 | 7.74 | 3.0 |
| 4 (Citrate plasma) | 10 | 12.11 | 4.2 |

9.4.2 Inter Assay Variation

The between assay variability is shown below:

| Sample | n | Mean (ng/mL) | CV (%) |
|-----------------------|----|--------------|--------|
| 1 (Li-Heparin plasma) | 30 | 0.21 | 8.7 |
| 2 (Serum) | 30 | 1.99 | 6.3 |

| | | | |
|--------------------|----|-------|-----|
| 3 (Serum) | 30 | 7.41 | 6.3 |
| 4 (Citrate plasma) | 30 | 11.53 | 6.2 |

9.5 Recovery

Samples have been spiked by adding 17- α -OHP solutions with known concentrations. The recovery (%) was calculated by multiplying the ratio of measured and expected values with 100.

| | Li-Heparin plasma | Citrate plasma | EDTA plasma | Serum |
|-----------------------|-------------------|----------------|-------------|-------|
| Concentration (ng/mL) | 0.17 | 0.46 | 0.72 | 1.92 |
| Average Recovery (%) | 97.9 | 96.0 | 95.7 | 93.8 |
| Range of Recovery (%) | from | 92.6 | 89.2 | 93.3 |
| | to | 103.8 | 103.2 | 98.6 |

9.6 Linearity

Samples were measured undiluted and in serial dilutions with standard 0. The recovery (%) was calculated by multiplying the ratio of expected and measured values with 100.

| | Citrate plasma | Serum | EDTA plasma | Li-Heparin plasma |
|-----------------------|----------------|-------|-------------|-------------------|
| Concentration (ng/mL) | 1.31 | 6.99 | 10.90 | 15.20 |
| Average Recovery (%) | 104.1 | 104.0 | 108.0 | 105.7 |
| Range of Recovery (%) | from | 93.9 | 92.2 | 102.4 |
| | to | 111.4 | 108.8 | 113.8 |

10 LIMITATIONS OF USE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice.

Any improper handling of samples or modification of this test might influence the results.

10.1 Interfering Substances

Haemoglobin (up to 4 mg/mL), Bilirubin (up to 0.5 mg/mL) and Triglyceride (up to 7.5 mg/mL) have no influence on the assay results.

10.2 Drug Interferences

| Substance | Concentration | Mean Bias |
|----------------------|---------------|-----------|
| | ng/mL | ng/mL |
| Coumestrol | 5 – 500 | -0.34 |
| Daidzein | 5 – 500 | 0.59 |
| Ethisterone | 23.4 – 2340 | -0.48 |
| Fulvestrant | 5 – 500 | -0.11 |
| Genistein | 5 – 500 | 0.28 |
| Levonorgestrel | 0.2 – 20 | 0.27 |
| Mifepristone | 23.4 – 2340 | -0.53 |
| Prednisolone | 0.2 – 20 | -0.12 |
| Prednisone | 0.2 – 20 | -0.59 |
| Secoisolariciresinol | 5 – 500 | 0.23 |

Daidzein will increase the measured 17-OH Progesterone concentration in a sample on average by more than 0.5 ng/mL. Mifepristone and Prednisolone will decrease the measured 17-OH Progesterone concentration in a sample on average by more than 0.5 ng/mL. All other tested substances will not change the 17-OH Progesterone concentration by more than ± 0.5 ng/mL.

10.3 High-Dose-Hook Effect

"High Dose Hook Effect" is not detected in the range between 0 - 640 ng/mL.

11 LEGAL ASPECTS

11.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DIASource.

11.2 Therapeutical Consequences

Therapeutical consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutical consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutical consequences.

11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2. are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

12 REFERENCES

1. Chrousos GP et al. Late onset 21-hydroxylase deficiency mimicking idiopathic hirsutism or polycystic ovarian disease. *Annals Intern Med.* 1982; 96,143-148.
2. Choi J-H, and Yoo H-W. Management issues of congenital adrenal hyperplasia during transition from pediatric to adult care. *Korean J. Pediatr.* 2017; 60(2), 31-37.
3. Eisenhoefer G et al. Reference intervals for plasma concentrations of adrenal steroids measured by LC-MS/MS: Impact of gender, age, oral contraceptives, body mass index and blood pressure status. *Clin. Chim Acta.* 2017; 470, 115-124.
4. Honour JW. 17-Hydroxyprogesterone in children, adolescents and adults. *Annals of Clinical Biochemistry,* 2014; 51(4), 424-440.
5. Hughes IA, Riad-Fahmy D, and Griffith K. Plasma 17OH-progesterone concentrations in newborn infants. *Archives of Disease in Childhood.* 1979; 54, 347-349.
6. Ishii T et al. Guidelines for diagnosis and treatment of 21-hydroxylase deficiency (2014 revision). *Clin. Pediatr. Endocrinol.* 2015; 24(3), 77-105.
7. Kang MJ et al. Relationships of basal level of serum 17-hydroxyprogesterone with that of serum androstenedione and their stimulated responses to a low dose of ACTH in young adult patients with congenital adrenal hyperplasia due to 21-hydroxylase deficiency. *J. Korean Med. Sci.* 2011; 26(11), 1454-60.
8. Maas KH et al. Relationship between 17-Hydroxyprogesterone responses to human chorionic gonadotropin and markers of ovarian follicle morphology in women with polycystic ovary syndrome. *J. Clin. Endocrinol. Metab.* 2015; 100(1), 293-300.

Other translations of this Instructions for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

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