CE
0344

HBsAg Elisa

KAPG4SGE3/KAPG4SGE11

History

Summary of change :

Previous Version : 190318/1	Current Version : 200224/1
Multilanguage IFU	Addition of the following sentence at the end of the English IFU: "Other translations of this Instructions for Use can be downloaded from our website: https://www.diasource-diagnostics.com/ "

HBsAg Elisa

en

For in-vitro qualitative detection of hepatitis B surface antigen (HBsAg) in human serum or plasma.

KAPG4SGE3 : 96 tests / KAPG4SGE11 : 480 tests

IN VITRO DIAGNOSTIC USE **IVD**

 DIAsource ImmunoAssays SA - Rue du Bosquet 2, B-1348 Louvain-la-Neuve, Belgium
Tel: +32 10 84 99 11 - Fax : +32 10 84 99 90

1) INTENDED USE

HBsAg Elisa is an enzyme immunoassay diagnostic kit for in-vitro qualitative detection of hepatitis B surface antigen (HBsAg) in human serum or plasma (heparin, citrate or EDTA).

2) SUMMARY AND TEST EXPLANATION

The hepatitis B surface antigen (HBsAg) is the first marker that appears in the blood following infection with hepatitis B virus (HBV) some days or weeks before clinical symptoms manifest. It is a lipoprotein polypeptide which constitutes the external envelope of the HB virus. The detection of HBsAg in human serum or plasma indicates an ongoing HBV infection, either acute or chronic. Testing of additional HBV markers is needed to define the specific disease state. HBsAg assays are used not only to diagnose HBV infections but also to monitor the course of the disease and the efficacy of antiviral therapy.

HBsAg Elisa is a fast test for the qualitative detection of the presence of HBsAg in serum or plasma (heparin, citrate or EDTA) specimen. The test utilizes monoclonal and polyclonal (anti-guinea pig) antibodies to selectively detect elevated levels of HBsAg in serum or plasma.

Specimens which are non-reactive by **HBsAg Elisa** are considered negative for HBsAg. Specimens with positive reaction should be retested in duplicate.

In case of a reactive repeat reaction, the specimen should be confirmed for HBsAg reactivity with validated confirmatory reagents .

Only confirmed positive specimens are considered to contain HBsAg.

3) TEST DESCRIPTION

HBsAg Elisa is a solid-phase enzyme immunoassay (ELISA= enzyme-linked immuno-sorbent assay) based on the "sandwich principle".

The solid phase of the microtiter plate is made of polystyrene wells coated with mouse monoclonal antibodies specific for HBsAg; whereas guinea pig polyclonal antibody purified by affinity chromatography is used to prepare the anti-HBs•peroxidase (horseradish) conjugate in the liquid-phase.

When a serum or plasma specimen containing HBsAg is added to the anti-HBs antibody-coated wells together with the peroxidase conjugated anti-HBs antibody and incubated, an antibody-HBsAg-antibody-peroxidase complex will form on the wells.

After washing the microtiter plate to remove unbound material, a solution of TMB substrate is added to the wells and incubated. A color develops in proportion to the amount of HBsAg bound to Anti-HBs. The peroxidase-TMB reaction is stopped by addition of sulfuric acid. The optical density of developed color is read with a suitable photometer at 450nm with a selected reference wavelength within 620 to 690nm*¹.

The test principle is shown also in the following figure.

A Specimen containing HBsAg:




1. Plate well (Anti-HBs) + specimen (HBsAg)+ Anti-HBs•peroxidase → Anti-HBs•HBsAg•(Anti-HBs•peroxidase) sandwich complex
2. Sandwich complex + TMB substrate solution → Light blue to blue color
3. Add sulfuric acid to stop the color development → Read OD at 450nm (reference wavelength 620-690 nm*¹)

B Specimen without HBsAg:

1. Plate well (Anti-HBs) + specimen (no HBsAg) + Anti-HBs•peroxidase → Anti-HBs (on the well)
2. Anti-HBs (on the well) + TMB substrate solution → Colorless to light blue color
3. Add sulfuric acid to stop the color development → Read OD at 450nm (reference wavelength 620-690 nm^{*1})

4) DESCRIPTION OF MATERIALS PROVIDED

- Item 1 - 6 in the following reagent table should be refrigerated at 2-8°C . Washing Solution D (20X) and stop solution can be stored at 2-30°C.

ITE MS	Components	Description	Qt. per 96 tests	Qt. per 480 tests
(1)	Anti-HBs Plate 	Microtiter plate coated with mouse monoclonal anti-HBs.	1 plate	5 plates
(2)	Anti-HBs Peroxidase Solution 	Polyclonal Anti-HBs HRPO conjugate, diluted in buffer with protein stabilizers. Preservatives: 0.003% Gentamycin and 0.01% Thimerosal. Dye: phenol red.	1 bottle, 8 ml	1 bottle, 25 ml
(3)	HBsAg Positive Control 	Inactivated human serum positive for HBsAg but non-reactive for anti-HCV and anti-HIV1/2, diluted in buffer with protein stabilizers. Preservatives: 0.003% Gentamycin and 0.01% Thimerosal.	1 bottle, 1.5 ml	1 bottle, 3 ml
(4)	HB Negative Control 	Serum non-reactive for HBV markers, anti-HCV and anti-HIV1/2, diluted in buffer with protein stabilizers. Preservatives: 0.003% Gentamycin and 0.01% Thimerosal.	1 bottle, 2 ml	1 bottle, 3 ml
(5)	Chromogenic TMB concentrate 	3,3',5,5'-tetramethylbenzidine (TMB) in an organic base.	1 bottle, 12 ml	1 bottle, 35 ml
(6)	Substrate buffer 	Acetic Acid Buffer with Urea Hydrogen Peroxydase.	1 bottle, 12 ml	1 bottle, 35 ml
(7)	Conc. Washing Solution (20X) 	Concentrated Phosphate buffer with Tween-20	1 bottle, 58 ml	1 bottle, 250ml
(8)	Stop Solution 	2N sulfuric acid	1 bottle, 12 ml	1 bottle, 50 ml

- OTHER REQUIRED MATERIALS, BUT NOT PROVIDED

ITEMS	Components
(1)	50µl, 100µl micropipettes and tips are needed
(2)	Water bath or incubator with temperature control at +37°C.
(3)	Plate washing equipment.
(4)	ELISA Microwell Reader: Dual wavelength 450nm with 620-690nm as reference wavelength, bandwidth 10nm ^{*1} .
(5)	Fully automatic EIA micro-plate analyzer is optional. User should validate the automatic EIA micro-plate analyzer in combination with the kit.

4.1) Storage condition and Stability of the kit and components

Kit/components	Storage temp.	State	Stability
HBsAg Elisa KIT	+2 to +8°C	Original	18 months
		Once open	1 month
HBsAg Positive Control	+2 to +8°C	Original	18 months
		Once open	1 month
HB Negative Control	+2 to +8°C	Original	18 months
		Once open	1 month
Anti-HBs Plate	+2 to +8°C	Original	24 months
		Once open	1 month
Anti-HBs • HRPO Conjugate Solution	+2 to +8°C	Original	18 months
		Once open	1 month
Concentrated Washing Solution (20X)	Room temp.	Original	24 months
		Once open	1 month
20X Diluted Washing Solution	Room temp.	Diluted	2 days
	+2 to +8°C	Diluted	1 week
Chromogenic TMB concentrate	+2 to +8°C	Original	24 months
		Once open	1 month
Substrate buffer	+2 to +8°C	Original	24 months
		Once open	1 month
TMB substrate mixture	Room temp.	Mixture	6 hours
2N Sulfuric Acid	Room temp.	Original	24 months
		Once open	1 month

5) Instructions for Use

5.1) Warnings:

- 5.1.1) This reagent kit is for professional use only.
- 5.1.2) **IVD** This reagent kit is for *in vitro* diagnosis only.
- 5.1.3) Bring all kit reagents and samples to room temperature (+20 to +30°C) and mix carefully before use.
- 5.1.4) Do not use reagent beyond its expiration date.
- 5.1.5) Do not interchange reagents between different lots.
- 5.1.6) Do not put pipette in mouth.
- 5.1.7) Do not smoke or eat in areas where specimens or reagents are handled.
- 5.1.8) All kit components and specimens should be regarded as potential health hazards. It should be used and discarded according to your laboratory's safety procedures. Such safety procedures probably include the wearing of protective gloves and avoiding the use of aerosols.
- 5.1.9) Potential infectious specimens and non-acid containing spills or leakages should be wiped up thoroughly with 5% sodium hypochlorite or treated in accordance with your practice for potential bio-hazard control.
- 5.1.10) Prior to disposing used specimens and kit reagents as general waste; it should be treated in accordance with the local practice of potential bio-hazardous waste or treated as follows:
Both liquid and solid waste should be autoclaved at +121°C for at least 30 minutes.
Solid waste can also be incinerated.
Non-acidic liquid waste can be treated with sodium hypochlorite diluted to a final concentration of 1%.
Acidic liquid wastes must be neutralized before treatment with sodium hypochlorite as mentioned above and should stand for 30 minutes to obtain effective disinfection.
- 5.1.11)  Stop Solution is an irritant to skin, eyes, respiratory tract and mucous membranes. Avoid contact of the stop solution with skin and mucous membranes. In case of contact, flush immediately with abundant amounts of water. In case of inhalation, find fresh air and seek medical attention in case of pain.
- 5.1.12)   Chromogenic TMB concentrate contains organic solvent, which is flammable : danger of serious irreversible effects through inhalation, in contact with skin and if swallowed. Chromogenic TMB concentrate contains dimethyl sulfoxide, an irritant to skin and mucous membranes.

- 5.1.13) Although all human sourced material are tested non reactive for Anti-HCV and Anti-HIV, and inactivated at +56°C for one hour, the reagent shall be handled as potential infectious material.*²

5.2) Residual risks

- 5.2.1) Incorrect operation could lead to wrong results.
- 5.2.2) Using the kit after expiration date might lead to wrong results.
- 5.2.3) Inefficient Washer and ELISA Reader could affect the final result.
- 5.2.4) Storage at higher temperature could affect the shelf life of the product.
- 5.2.5) The product is for single use only. Reuse could lead to wrong results.
- 5.2.6) Any change to the test procedure may lead to incorrect results which could cause diagnostic error.

5.3) Specimen Collection and Preparation for Analysis

- 5.3.1) No special preparation of the patient is required prior to blood collection. Blood should be collected by approved medical techniques.
- 5.3.2) Either serum or plasma specimens can be used with this test kit. Whole blood specimen should be separated as soon as possible in order to avoid hemolysis. Any particulates (e.g. fibrin clots, erythrocytes) contained in the specimen should be removed prior to use.
- 5.3.3) Specimens must be stored at +2 to +8°C and avoid heat-inactivation to minimize deterioration. For long-term storage, they should be frozen below -20°C. Storage in self-defrosting freezer is not recommended.
- 5.3.4) Frozen specimens must be thoroughly thawed and mixed homogenously before test.
- 5.3.5) Avoid multiple freeze-thaw procedures.

WARNING

- 1. The specimen must not contain any compounds of AZIDE, which inhibits the peroxidase activity.
- 2. Incompletely coagulated sera and microbial-contaminated specimens should not be used.

5.4) Storage conditions and Stability of Reagents

- 5.4.1) The kit must be stored at +2 to +8°C. Do not freeze.
- 5.4.2) Strips of the plate should be used within one month after opening the original aluminum foil bag. The unused strips should be kept in the aluminum foil bag and taped tightly.
- 5.4.3) Return reagents to +2 to +8°C immediately after use.
- 5.4.4) Washing Solution (20X) Concentrate can be stored at room temperature to avoid crystallization, because the kits are stored and shipped at +2 to +8°C. If crystals have been precipitated before use, warm up the solution at 37°C in a water bath until the crystals dissolve.

5.5) Plate Washing Procedure

- 5.5.1) Preparation of washing solution:
Dilute Washing Solution (20X) Concentrate with distilled or de-ionized water to obtain a 1:20 dilution. Do not use tap water.
- 5.5.2) Plate washing:
 - (a) For plate washer with overflow aspirating function: 6 cycles with at least 0.5 ml washing buffer per well per cycle.
 - or
 - (b) For plate washer without overflow aspirating function: 8 cycles with at least 0.35 ml washing buffer per well per cycle.
- 5.5.3) Blot dry by inverting the plate and tapping firmly onto absorbent paper. Too much residual wash buffer will cause false results.

WARNING : Improper washing will cause false results.

5.6) Test Procedure

5.6.1) Bring all reagents and specimens to room temperature (+20 to +30°C) before assay. Adjust water bath or incubator to +37±1°C.

5.6.2) Reserve one well for blank. Add 50 µl of each control or specimen to appropriate wells of the microtiter plate (3 Negative Controls and 2 Positive Controls).

NOTE:

- a. Use a clean pipette tip for each sampling to avoid cross-contamination.
- b. Each plate needs respective negative controls, positive controls and blank well.
- c. Do not use any cut-off value established for other plates of HBsAg Elisa.

5.6.3) Add 50 µl of Anti-HBs Peroxidase Solution to each well except the blank.

NOTE: Do not touch the wall of the plate wells to prevent contamination.

5.6.4) Gently tap the plate.

5.6.5) Remove the protective backing from the Adhesive Slip and press it onto the reaction plate, so that it is tightly sealed.

5.6.6) Incubate the reaction plate at +37±1°C in a water bath or incubator for **80 minutes**.

5.6.7) At the end of the incubation period, remove and discard the Adhesive Slip and wash the plate in accordance with 5.4) Plate Washing Procedure.

5.6.8) Select one of the following two methods for color development:

- A. Mix equal volumes of **Chromogenic TMB concentrate and Substrate Buffer** in a clean container immediately prior to use. Add **100 µl** of the mixture solution to each well including the blank well.
- B. Add **50µl** of **Chromogenic TMB concentrate** first, and then add **50 µl** of **Substrate buffer** into each well including the blank. Mix well carefully.

NOTE: Chromogenic TMB concentrate should be colorless to light blue; otherwise, it should be discarded. The mixture of Chromogenic TMB concentrate and Substrate buffer should be used within 6 hours after mix. The mixture should be kept away from intense light.

5.6.9) Cover the plate with a black cover and incubate at **room temperature** for **30 minutes**.

5.6.10) Stop the reaction by adding **100 µl** of stop solution to each well including the blank.

5.6.11) Determine the absorbance of Controls and test specimens within 30 minutes with a precision spectrophotometer at 450/620-690nm (450nm reading wavelength with 620 – 690 nm reference wavelength) *1.

Use the blank well to blank spectrophotometer.

NOTE: The color of the blank should be colorless to light yellow; otherwise, the test results are invalid.
Substrate blank : absorbance value must be less than 0.100.

5.7) Calculations of Results

5.7.1) Calculation of the NC (Mean Absorbance of Negative Control).

Example: Sample No.	Absorbance
1	0.022
2	0.025
3	0.023

$$NC = (0.022 + 0.025 + 0.023) / 3 = 0.023$$

NCx should be ≤0.1, otherwise, the test is invalid.

5.7.2) Calculation of the PC (Mean Absorbance of Positive Control)

Example: Sample No.	Absorbance
1	1.432
2	1.508

$$PC = (1.432 + 1.508) / 2 = 1.470$$

PC should be ≥0.6, otherwise, the test is invalid.

5.7.3) Calculation of the **P - N Value**

$$P - N = PC - NC$$

Example: NC = 0.024

PC = 1.470

$$P - N = 1.470 - 0.024 = 1.446$$

P - N Value must be ≥0.5, otherwise, the test is invalid.

5.7.4) Calculation of the Cutoff Value

$$\text{Cutoff Value} = NC + 0.025$$

Example:

$$\text{Cutoff Value} = 0.023 + 0.025 = 0.048$$

5.8) Validity of Test Runs

5.8.1) NC should be ≤ 0.1 ; otherwise, the test is invalid.

5.8.2) PC should be ≥ 0.6 , otherwise, the test is invalid.

5.8.3) P - N Value must be ≥ 0.5 , otherwise, the test is invalid.

NOTE: Negative Control: absorbance value must be less than or equal to 0.100 after subtracting the blank.

5.9) Interpretation of Results

5.9.1) Specimens with absorbance values **LESS** than the **Cutoff Value** are **NON-REACTIVE** and are considered **NEGATIVE for HBsAg**.

5.9.2) Specimens with absorbance value **GREATER** than or **EQUAL** to the **Cutoff Value** are considered **INITIALLY REACTIVE**. The original specimens must be retested in duplicate.

5.9.3) If both absorbance values in the retest are **less** than the cutoff value, the specimens are considered **NEGATIVE for HBsAg**.

If in the retest at least one of the two absorbance values is **GREATER** than or **EQUAL** to the **Cutoff Value** then the specimens are considered as **repeated HBsAg positive**. The repeated positive specimen shall be confirmed with Antisurgen B.

5.10) Troubleshooting

If the result cannot be reproduced, perform a preliminary troubleshooting by checking the possibilities listed below:

5.10.1) Improper washing procedure.

5.10.2) Contamination with positive specimen.

5.10.3) Wrong volume of sample, conjugate or substrates.

5.10.4) Contamination of the well rim with conjugate.

5.10.5) Improper specimen, such as hemolyzed serum or plasma, specimen containing sediments and specimen not thoroughly mixed before use.

5.10.6) Wrong incubation time or temperature.

5.10.7) Obstructed or partial obstructed washer aspirate/dispense head and needles.

5.10.8) Insufficient aspiration.

5.11) Limitations and Interferences

5.11.1) This reagent kit is to be used for un-pooled human serum or plasma only.

5.11.2) The reagent kit has not been validated for use with cadaveric samples.

5.11.3) A negative HBsAg result without other evidence should not be used to exclude an HBV infection.

5.11.4) Interfering Substances:

The following results were obtained in respective experiments:

1. No interferences with different anticoagulants such as lithium heparin, EDTA, citrate have been observed.

2. Heat-treated specimens (+60°C, 10 hours) exhibited diminished HBsAg titer.

3. No cross reactivity was detected using specimens deriving from patients with a) other infections by HAV, EBV, CMV, HSV, VZV, Lyme Borreliosis, HCV, HIV, b) other disease states such as chronic renal failure, hemodialysis, autoimmune hepatitis, liver cirrhosis, and c) presence of certain antibodies like HAMA, GAD, IA2, APS).

4. Samples containing potential interfering substances [e.g. triglycerides (lipemia), hemoglobin (hemolysis), bilirubin (icterus), monoclonal immunoglobulin components, elevated levels of autoimmune antibodies (rheumatoid factor-RF, antinuclear antibodies-ANA, or antimitochondrial antibodies-AMA)] and samples from pregnant women did not interfere with the HBsAg Elisa assay.

5.12) Performance Characteristics

5.12.1) Diagnostic Specificity

Results from the European Performance Evaluation for DIAsource HBsAg Elisa - Reactivity of HBV Negative "Donor" and "Clinical" Specimens.

Total No. of Specimens	DIAsource HBsAg Elisa					
	N	Neg	*IR	**RR	Confirmed	False Positive
HBV negative (clinical specimen)	213	211	2	2	0	2
HBV negative (donor specimen)	5501	5479	22	22	0	22
Total	5714	5690	24	24	0	24

*IR: initial reactive **RR: repeat reactive

Diagnostic specificity = $5690/5714 = 99.58\%$

5.12.2) Diagnostic Sensitivity

1. The diagnostic sensitivity determined in the European performance evaluations yielded the following results:

Sample	No. of sample	Reactive	Sensitivity
HBsAg positive sera	400	400	100%

Diagnostic sensitivity = $400/400 = 100\%$

2. The evaluation of the HBsAg Sensitivity Panel PHA807 (BBI, USA) yielded the following results:

Panel	No. of panel members	Results
BBI HBsAg Sensitivity Panel PHA807	PHA807-01 ~ 21	Equivalent to or better than the competitors.

3. The evaluation of the HBsAg Mixed Titer Performance Panel PHA205 (BBI, USA) yielded the following results:

Panel	No. of panel members	Results
BBI HBsAg Mixed Titer Performance Panel PHA205	PHA205-01 ~ 25	Equivalent to or better than the competitors.

4. The HBsAg Mixed Titer Performance Panel VHA601 (BBI, USA) evaluation yielded the following results:

Panel	No. of panel members	Results
BBI HBsAg Mixed Titer Performance Panel VHA601	VHA601-01 ~ 06	Met BBI's expected results.

5. The evaluation of the HBsAg Mixed Titer Performance Panel QHA711 (BBI, USA) yielded the following results:

Panel	No. of panel members	Results
BBI HBsAg Mixed Titer Performance Panel QHA711	QHA711-01 ~ 06	Met BBI's expected results.

6. INTS HBsAg Subtype Panel Test Result:

The results of the evaluation demonstrate that all known HBsAg serotypes (HBsAg subtypes) available in the INTS HBsAg subtype panel can be detected by both HBsAg Elisa and the reference assay up to $10^4x \sim 10^6x$ dilutions.

7. Summary table of the evaluation of all tested seroconversion panels:

PPanel ID	HBsAg Positive Result From Date of First Blood Collection			
	Reference HBsAg Assay (A in days)	HBsAg Elisa (B in days)	Difference in Bleed Days (A-B)	Difference in Bleed #s (panel member No.)
PHM907(ay)	50	50	0	0
PHM916(ay)	62	65	-3	-1
PHM920(ad)	26	26	0	0
PHM921(ad)	0	0	0	0
PHM930(ad)	3	3	0	0
PHM933(ad)	7	9	-2	-1
PHM934(ad)	0	0	0	0

PHM935A (ad)	9	21	-12	-3
PHM935B (ad)	231	217	14	-1
NabiSB0411	6	6	0	0
PHM903(ad)	6	10	-4	-1
PHM904(ad)	7	18	-11	-1
PHM906(ad)	137	0	137	1
PHM909(ad)	9	9	0	0
PHM910(ad)	35	18	17	1
PHM912(ad)	42	0	42	7
PHM914(ad)	0	146	-146	-1
PHM915(ind)	12	28	-16	-5
PHM917(ind)	36	>43	->7	->2
PHM918 (ad)	7	7	0	0
PHM923 (ay)	15	21	-6	-1
PHM924 (ad)	29	35	6	1
PHM925(ind)	8	14	-6	-1
PHM926(ind)	13	15	-2	-1
PHM927(ind)	4	7	-3	-1
PHM928 (ad)	9	9	0	0
PHM929 (ad)	14	18	-2	-1
PHM931(ind)	19	21	-2	-1
PHM932 (ad)	61	61	0	0
PHM935A (ad)	21	28	-7	-2
Over all	-----	-----	Sum = -13	-----
		Average	-13/30 = 0.4 days	-----

Summary of the evaluation of all tested seroconversion panels:

The results have shown that the HBsAg Elisa is nearly equivalent to the reference assay. On average the HBsAg Elisa test is only 0.4 days later in the 30 tested seroconversion panels in comparison to the reference assay.

5.12.3) Mutants

1. Mutant Panel

A Hepatitis B Pre Core Mutant panel of Teragenix, USA was tested. HBsAg Elisa can detect the results while panel's dilution ratio is 1:2 or 1:5.

This panel consists of 3 samples Core Promoter (A1762/G1764) wild type; PreCore codon 28 mutant/wild type infection, of 4 samples Core Promotor (A1762/G1764) wild type; PreCore codon 28 mutant, of 2 samples Core Promotor (T1762/A1754) MUTANT/(A1762/G1754) wild type mixed infection.

2. Mutant samples

The detection rate of HBsAg Elisa is 17/21 which is better than reference kit's (14/21).

Mutation	HbsAg Elisa	Reference kit	Mutation	HbsAg Elisa	Reference kit
	Int	Int		Int	Int
113A114	positive	positive	C137W	positive	negative
P120L	positive	positive	C139Y	positive	positive
R122I	negative	negative	K141E	positive	positive
R122T	positive	positive	P142L/G145R	negative	negative
E122I	positive	positive	P142L	positive	positive
T123N	positive	positive	D144G	positive	negative
122RA123	negative	negative	G145K	positive	negative
C124R	negative	negative	C147T/C148T	positive	positive
T131I	positive	positive	E164D	positive	positive
T131P	positive	positive	S174N	positive	positive
Y134S	positive	positive			

Mutation	HbsAg Elisa	Reference kit	Mutation	HbsAg Elisa	Reference kit
Detection Rate	8/11	8/11	Detection Rate	9/10	6/10

5.12.4) 25 same day fresh positive samples

HBsAg Elisa detected 25 same day fresh positive samples. The criteria of the chapter 3.1.9 of the 2009/886/EC was fulfilled.

5.12.5) Analytical Sensitivity

Standard Material	Calculated Sensitivity
WHO 2 nd international Standard ^{**}	0.03 IU/mL
PEI standard Subtype ad	0.08 PEI U/mL
PEI standard Subtype ay	0.09 PEI U/mL

^{**} WHO Second International Standards for HBsAg, subtype adw2, genotype A, NIBSC code: 00/588

5.12.6) Precision

1. Intra-assay reproducibility

Intra-assay precision was determined using one positive control sample and two patient serum samples of different HBsAg concentration (slightly above the cutoff level and at medium level) which were analyzed in replicates of 20 in a single run over 3 days.

Run	Sample ID	Mean	Standard Deviation	Coefficient of Variation
1	PC 1 :32	0.7541	0.0857	11.36
1	PS 1 :32	2.2766	0.1571	6.90
1	PS 1 :64	1.3159	0.1168	8.88
2	PC 1 :32	0.9671	0.0358	3.70
2	PS 1 :32	2.7325	0.1025	3.75
2	PS 1 :64	1.5822	0.0888	5.61
3	PC 1 :32	0.9669	0.0909	9.40
3	PS 1 :32	2.6088	0.2367	9.07
3	PS 1 :64	1.6502	0.1499	9.08

The calculated CV's ranged between 3.7 % and 11.36 %. (= acceptable value for an immunoassay in microtiter plate format).

2. Inter-assay reproducibility

The precision evaluation experiments were performed over 10 operating days is using five serum samples (with borderline positive and clearly above cutoff value HBsAg levels).

Sample ID	N	Mean	Standard Deviation	Coefficient of Variation
F1	10	0.0609	0.0185	30.41
F2	10	0.0770	0.0189	24.56
F3	10	0.1000	0.0245	24.51
F4	10	0.1786	0.0462	25.87
F5	10	0.6107	0.1412	23.12

The calculated CV's ranges between 30.4 for an HBsAg negative sample and 23.1 % for an HBsAg low positive sample (= acceptable values for the inter-assay imprecision of an immunoassay in microtiter plate format).


5.12.7) Antigen Excess/High-dose hook effect

This was performed testing a serum sample with a very high HBsAg concentration of 125 mg/L in serial dilution with the DIAsource HBsAg Elisa assay. No high-dose hook effect was observed.

5.12.8) Traceability

The DIAsource HBsAg Master Calibrator has been calibrated against the British Working Standard for HBsAg (NIBSC-Code: 01/476-006) using the **HBsAg Elisa** assay. The relative potency (ratio) of the British Working Standard for HBsAg versus the DIAsource HBsAg Master Calibrator is 4.081 (3.853-4.325 95% CI). The concentration of the Positive Control of **HBsAg Elisa** assay has been determined against the DIAsource HBsAg Master Calibrator and was established with 42 IU/ml \pm 20%.

5.13) Flow chart of the test procedure

6) NOTES

- *1 The reference wavelength of the photometer to be used can be 620 nm to 690 nm. However, the user should validate the photometer in combination with **HBsAg Elisa** before use.
- *2 Incomplete inactivation of hepatitis B virus after heat treatment at +60°C for 10 hours, J. Infect. Dis. 138:242-244.

7) BIBLIOGRAPHY

1. Blumberg BS, Alter HJ, Visnich S. A „new“ antigen in leukemia sera. JAMA, 1965;191:101- 106.
2. Dane DS, Cameron CH, Briggs M. Virus-like particles in serum of patients with Australia-antigen-associated hepatitis. Lancet 1970; 1: 695 - 698.
3. Aach RD, G ham JW, Paker CW. Detection of Australia antigen by radioimmunoassay. Proc Natl Acad Sci. USA 1971;68:1056-1060.
4. Kim CY, Tikes JG. Purification and biophysical characterization of hepatitis antigen. J Clin Invest. 1973; 52:1176-1186.
5. Wolters G. Kuijpers LP, Kacaki J, Schuurs AH, Enzyme linked immunosorbent assay for hepatitis B surface antigen. J Infect. Dis 1977;136:Suppl 311-377.
6. Shih JW, Cote PJ Jr, Dapolito GM, Gerin JL. Production of monoclonal antibodies against hepatitis B surface antigen (HBsAg) by somatic cell hybrids. J Virol Methods. 1980;1:257-273.
7. Hoofnagle JH, Di Bisceglie AM. Serologic diagnosis of acute and chronic viral hepatitis. Semin Liver Dis. 1991;11:73-83

Revision date : 2020-02-24

Other translations of this Instructions for Use can be downloaded from our website:
<https://www.diasource-diagnostics.com/>

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