# **IGFBP-1 Elisa**

***KAPME01***

# History

---

**Summary of change :**

<b>Previous Version :</b> 200224/1	<b>Current Version :</b> 210628
No Spanish version	Spanish version added

# IGFBP-1-ELISA

Enzyme Immunoassay for the Quantitative Determination of  
Human Insulin-like Growth Factor Binding Protein-1  
KAPME01  
IN VITRO DIAGNOSTIC USE

DIAsource ImmunoAssays S.A. - Rue du Bosquet 2, B-1348 Louvain-la-Neuve, Belgium - Tel: +32 10 84 99 11 - Fax : +32 10 84 99 91

## INTENDED USE

This enzyme immunoassay kit is suited for measuring IGFBP-1 in human serum and plasma.

## INTRODUCTION

The Insulin-like Growth Factors I and – II are free in body fluids and tissues but are bound to specific binding proteins. Until today seven different binding proteins (IGFBP-1 to -7) can be differentiated additionally several IGFBP-related proteins have also been detected. Bioavailability of IGF is regulated by these IGFBPs or their proteolytic cleavage which reduces affinity to IGF. The IGFBPs as well as their proteolytic fragments can also exert IGF-independent effects, like influencing cell migration or proliferation.

IGFBP-1 (Placental Protein 12) consists of 234 amino acids and has a molecular weight of approximately 25kDa. The coding DNA region is located on chromosome 7 [1, 2]. IGFBP-1 is mainly synthesized by foetal and adult liver tissue and decidual endometrium. Intensity of Expression varies enduring menstruation with a maximal expression in the late secretory phase [3, 4]. Further IGFBP-1 expression seems to be regulated by Insulin concentration, with Insulin inhibiting the expression. Insulin regulation results in diurnal fluctuations of up to factor 10 [5]. IGFBP-1 is posttranslational modified by phosphorylation of serine residues 101, 119 and 169. Phosphorylation has physiological relevance as it increases affinity of IGFBP-1 to IGF. In adult humans phosphorylated IGFBP-1 of the liver is the predominant form in circulation. IGFBP-1 produced by endometrial tissue is significantly less phosphorylated than the liver originated form [6].

In pregnancy IGFBP-1 maternal serum concentration increases significantly with maximal values in the second trimester or 22-23 week of gestation (75.8 ng/ml) [5] and decreases slowly until term. IGFBP-1 concentration is not only increased in maternal but also in foetal serum, extremely high concentrations are found in amniotic fluid. Here concentration can reach more than the 1000-fold of serum values [7]. Long-term changes of serum IGFBP-1 concentration can also be found in amniotic fluid. After birth the IGFBP-1 level of the child decreases until it reaches the low steady-state level of puberty and adulthood [8, 9].

Short term IGFBP-1 serum concentration is strongly influenced by nutrition level and therewith by insulin. Decreasing IGFBP-1 levels can be found enduring fasting or in diabetes; IGFBP-1 levels increase in case of intensive exercises [10-12].

Relevance of serum and amnion IGFBP-1 in diagnostics has been investigated in several areas.

### - Energy metabolism

Based on the influence of Insulin on IGFBP-1 serum concentrations IGFBP-1 is said to be a possible marker for insulin resistance.

In a small study Maddux et al were able to demonstrate with 23 non-diabetic patients, that IGFBP-1 serum concentration correlated very well with Glucose-uptake rate, even better than the HOMA index does [13].

### - Pregnancy

In pregnancy IGFBP-1 concentrations could be used for different applications. Specially, a potential diagnostic value was shown for membrane rupture and intra uterine growth disturbances. There are also data that IGFBP-1 in maternal serum might be a marker for trisomy 18, but here the IGFBP-1/IGFBP-2 ratio seems to be of value [14].

A significant difference in IGFBP-1 serum concentration was detected between healthy pregnant and diabetic and pre-eclamptic women (102,8 vs. 203,71 or 281,09 ng/ml respectively) [15].

Also the evaluation of IGFBP-1 as marker for membrane rupture showed a high specificity (97%) and sensitivity (75%) of IGFBP-1 in vaginal/cervical secrets. In case of intact membrane IGFBP-1 concentration was < 90ng/ml in the secretion. Enduring 8 hours after spontaneous or induced membrane rupture IGFBP-1 values increased significantly with a median concentration of 1900 ng/ml. In this study IGFBP-1 concentrations of >100ng/ml were set as threshold for detection of amnion fluid and therewith diagnosis of membrane rupture [16]. A positive predictive value of 97% indicates that IGFBP-1 is a suitable marker for premature membrane rupture [17].

## ASSAY PRINCIPLE

The DIAsource **IGFBP-1 Elisa KAPME01** is a so-called Sandwich-Assay. It utilizes two specific antibodies of high affinity. First the IGFBP-1 in the sample binds to the immobilized antibody on the microtiter plate. In the following step, the anti-IGFBP-1-Antibody binds in turn to the immobilised IGFBP-1. This is biotinylated and allows the binding of a streptavidin-peroxidase enzyme conjugate.

Subsequently, the peroxidase catalyzes an enzymatic reaction resulting in a blue coloration. The intensity of the blue color depends on the IGFBP-1 content of the sample. The reaction is stopped by the addition of stop solution and color intensity is quantified by measuring the absorption.

## WARNINGS AND PRECAUTIONS

### For In Vitro Diagnostic Use only. For Professional use only.

The DIIAsource kit is suitable only for in vitro diagnostics and not for internal use in humans and animals. Follow strictly the test protocol. Use the valid version of the package insert provided with the kit. Be sure that everything has been understood. DIIAsource will not be held responsible for any loss or damage (except as required by statute) howsoever caused, arising out of noncompliance with the instructions provided.

Do not use obvious damaged or microbial contaminated or spilled material.

**Caution: This kit contains material of human and/or animal origin. Therefore all components and patient's specimens should be treated as potentially infectious.**

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. The disposal of the kit components must be made according to the local regulations.

### Human Serum

Following components contain human serum: **Control Sera 1 / 2, Calibrators 0-6**

Source human serum for the control sera provided in this kit was tested and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV). No known methods can offer total security of the absence of infectious agents; therefore all components and patient's specimens should be treated as potentially infectious.

H- and P statements refer to undiluted substances and are for information purpose.

### Reagents

Contain as preservative **5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-4-isothiazolin-3-one** (<0.015%)

H317	May cause an allergic skin reaction.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P272	Contaminated work clothing should not be allowed out of the workplace.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P501	Dispose of contents/ container in accordance with local/ regional/ national/ international regulations.

### Substrate Solution

The TMB-Substrate contains 3,3',5,5' Tetramethylbencidine (<0.05%)

H315	Causes skin irritation.
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P305+P351+	IF IN EYES: Rinse cautiously with water for several minutes.
P338	Remove contact lenses, if present and easy to do. Continue rinsing.

### Stopping Solution

The Stopping solution contains 0.2 M acid sulphur acid (H<sub>2</sub>SO<sub>4</sub>)

H290	May be corrosive to metals.
H314	Causes severe skin burns and eye damage.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P301+P330+	IF SWALLOWED: rinse mouth.
P331	Do NOT induce vomiting.
P305+P351+	IF IN EYES: Rinse cautiously with water for several minutes.
P338	Remove contact lenses, if present and easy to do. Continue rinsing.
P309+P310	IF exposed or if you feel unwell: Immediately call a POISON CENTER or doctor/physician.

### General first aid procedures:

Skin contact: Wash affected area rinse immediately with plenty of water at least 15 minutes. Remove contaminated cloths and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing spread the eyelids.

Ingestion: After swallowing the product, if the affected person is conscious, rinse out the mouth with plenty of water: seek medical advice immediately.

## SAMPLES

### Sample type

Beside serum also EDTA- and Heparin Plasma can be used

### Specimen collection

Use standard venipuncture for the blood sampling. Haemolytic reactions have to be avoided. For blood sampling diurnal variations specially influenced by nutrition should be considered [5].

### Required sample volume: 20 µL

### Sample stability

- In firmly closable sample vials
- Storage at 20-25°C: 3 days
- Storage at -20°C: min. 2 years
- Freeze-thaw cycles max. 3

The storage of samples over a period of 2 years at -20°C, showed no influence on the reading. Freezing and thawing of samples should be minimized. 3 Freezing-Thawing showed no effect on samples.

### Interference

Triglyceride, bilirubin and hemoglobin in the sample do not interfere to a concentration of 100 mg/mL, 100 µg/mL or 5 mg/mL, respectively. However, the use of haemolytic, lipemic or icteric samples should be validated by the user.

### Sample dilution


- Dilution:  $\geq 1:16$  with Dilution Buffer
- Pipette **300 µL** Dilution Buffer (red colored) in PE-/PP-Tubes (application of a multi-stepper is recommended in larger series), add **20 µL Serum-** or **Plasma** (dilution factor 16). After mixing use **50 µL** diluted solution **within 1 hour per determination** in the assay.

## REAGENTS PROVIDED


The reagents listed below are sufficient for 96 wells including the calibration curve.

**Microtiter plate**, ready for use: **Microtiter plate** with 96 wells (8x12), divided up in 12 strips with 8 wells separately breakable, coated with anti-IGFBP-1 antibody and packed in a laminate bag.


**Calibrators 0-6**, lyophilized: contain native human IGFBP-1. Calibrators values are between 0 – 8 ng/ml (0, 0.1, 0.5, 1, 2, 4 and 8 ng/ml) IGFBP-1. Calibrators should be reconstituted with 500 µl Dilution Buffer each. Use 50 µl per well in the assay.


**Dilution Buffer**, 125 ml, ready for use, please use for reconstitution of samples, control and calibrators.


**Control 1 and 2**, lyophilized, contain human Serum and should be reconstituted in 250 µl Dilution Buffer each. The IGFBP-1 target values and the respective ranges are given on the vial label. The dilutions should be according to the dilution of the respected samples. Use 50 µl per well in the assay.


**Biotin Conjugate**, 6 ml, ready for use, contains biotinylated anti-human IGFBP-1 Antibody. Use 50 µl per well in the assay.


**Streptavidin HRP**, 12 ml, contains HRP (Horseradish-Peroxidase)-labelled Streptavidin. Ready for use. Use 100 µl per well in the assay.

**Washing Buffer**, 50 ml, 20 X concentrated solution. Dilute 1:20 with Aqua dest. The 1:20 diluted Washing Buffer is only limited stable. Please dilute only according to daily requirements.

**TMB-Substrate Solution** 12 ml, ready for use, horseradish-peroxidase-(HRP)-substrate, stabilised H<sub>2</sub>O<sub>2</sub>-Tetramethylbenzidine. Use 100 µl per well in the assay.

**Stopping Solution**, 12 ml, ready for use, 0.2 M sulphuric acid. Use 100 µl per well in the assay.

**Sealing tape** for covering of the microtiter plate, 2 x, adhesive.

## MATERIALS REQUIRED BUT NOT PROVIDED

Distilled (Aqua destillata) or deionized water for dilution of the Washing Buffer, 950 mL.  
Precision pipettes and multichannel pipettes with disposable plastic tips  
Polyethylene PE/Polypropylene PP tubes for dilution of samples  
Vortex-mixer  
Microtiter plate shaker (350 rpm)  
Microtiter plate washer (recommended)  
Micro plate reader ("ELISA-Reader") with filter for 450 and  $\geq 590$  nm

## TECHNICAL NOTES

### Storage Conditions

Store the kit at 2-8°C after receipt until its expiry date. The lyophilized reagents should be stored at -20 °C after reconstitution. Avoid repeated thawing and freezing.

### Storage Life

The shelf life of the components **after initial opening** is warranted for **4 weeks**, store the unused strips and microtiter wells **airtight** together with the desiccant at 2-8°C in the clip-lock bag, use in the frame provided. The **reconstituted components** calibrators 0-6 and Control Sera **1** and **2** must be stored at -20°C. For further use, thaw quickly but gently (avoid temperature increase above room temperature and avoid excessive vortexing). Up to 3 of the freeze-thaw cycles did not influence the assay. The 1:20 diluted Washing Buffer is stable for **4 weeks** at **2-8°C**

### Preparation of reagents

Bring all reagents to room temperature (20°C - 25°C) before use. Possible precipitations in the buffers have to be resolved before usage by mixing and / or warming. Reagents with different lot numbers cannot be mixed.

### Reconstitution

The **Calibrators 0 - 6** and **Controls 1 and 2** are reconstituted with the Dilution Buffer. It is recommended to keep reconstituted reagents at room temperature for 15 minutes and then to mix them thoroughly but gently (no foam should result) with a Vortex mixer.

### Dilution

After reconstitution dilute the Control Sera 1 and 2 with the Dilution Buffer in the same ratio as the sample.  
The required volume of Washing Buffer is prepared by 1:20 dilution of the provided 20 fold concentrate with Aqua dest.

### Assay Procedure

When performing the assay, Calibrators 0-6, Control Serum 1 and 2 and the samples should be pipette as fast as possible (e.g. <15 minutes). To avoid distortions due to differences in incubation times, Biotin Conjugate, Streptavidin HRP as well as the succeeding Substrate Solution should be added to the plate in the same order and in the same time interval as the samples. Stopping Solution should be added to the plate in the same order as Substrate Solution.

All determinations (Calibrators 0-6, Control Sera 1 and 2 and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

### Incubation

**Incubation at room temperature means: Incubation at 20 - 25°C.** The Substrate Solution, stabilized H<sub>2</sub>O<sub>2</sub>-Tetramethylbendidine, is photosensitive—store and incubation in the dark.

### Shaking

The incubation steps should be performed at mean rotation frequency of a particularly suitable microtiter plate shaker. We are recommending 350 rpm. Due to certain technical differences deviations may occur, in case the rotation frequency must be adjusted. Insufficient shaking may lead to inadequate mixing of the solutions and thereby to low optical densities, high variations and/ or false values, excessive shaking may result in high optical densities and/ or false values.

### Washing

Proper washing is of basic **importance** for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results calculations of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided Washing Buffer diluted to usage concentration. Washing volume per washing cycle and well must be 300  $\mu$ L at least.

The danger of handling with potentially infectious material must be taken into account.

When using an **automatic microtiter** plate washer, the respective instructions for use must be carefully followed. Device adjustments, e.g. for plate geometry and the provided washing parameters, must be performed. Dispensing and aspirating manifold must not scratch the inside well surface. Provisions must be made that the remaining fluid volume of every aspiration step is minimized. Following the last aspiration step of each washing cycle, this could be controlled, and possible remaining fluid could then be removed, by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

**Manual washing** is an adequate alternative option. Washing Buffer may be dispensed via a multistepper device, a multichannel pipette, or a squirt bottle. The fluid may be removed by dynamical swinging out the microtiter plate over a basin. If aspirating devices are used, care has to be taken that the inside well surface is not scratched. Subsequent to every single washing step, the remaining fluid should be removed by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

#### ASSAY PROCEDURE

Preparation of reagents		Reconstitution:	Dilution
	<b>Calibrators 0-6</b>	in <b>500 µL</b> Dilution Buffer	-
	<b>Control Serum 1</b>	in <b>250 µL</b> Dilution Buffer	<b>≥ 1:16</b> with Dilution buffer
	<b>Control Serum 2</b>	in <b>250 µL</b> Dilution Buffer	<b>≥ 1:16</b> with Dilution buffer
	<b>Washing Buffer</b>	-	<b>1:20</b> with distilled water
<b>Sample dilution:</b> with Dilution Buffer <b>≥ 1:16</b>			
Before assay procedure bring all reagents to room temperature <b>20-25°C</b> .			
<b>Assay Procedure in Double Determination:</b>			
Pipette	Reagents	Position	
50 µL	Biotin Conjugate	Pipette in <b>all</b> required number of wells	
50 µL	Calibrator <b>0 (0 ng/mL)</b>	A1/A2	
50 µL	Calibrator <b>1 (0.1 ng/mL)</b>	B1/B2	
50 µL	Calibrator <b>2 (0.5 ng/mL)</b>	C1/C2	
50 µL	Calibrator <b>3 (1 ng/mL)</b>	D1/D2	
50 µL	Calibrator <b>4 (2 ng/mL)</b>	E1/E2	
50 µL	Calibrator <b>5 (4 ng/mL)</b>	F1/F2	
50 µL	Calibrator <b>6 (8 ng/mL)</b>	G1/G2	
50 µL	Control Serum <b>1 (≥ 1:16 diluted)</b>	H1/G2	
50 µL	Control Serum <b>2 (≥ 1:16 diluted)</b>	A3/A4	
50 µL	Sample ( <b>≥ 1:16 diluted</b> )	in the rest of the wells according the requirements	
Cover the wells with the sealing tape.			
<b>Sample Incubation <u>with Shaking</u>: 1 h at 20°C - 25°C, 350 rpm</b>			
5 x 300 µL	Aspirate the contents of the wells and <b>wash</b> 5 x with 300 µL each Washing Buffer/ <b>well</b>		In each well
100 µL	Streptavidin HRP		In each well
Cover the wells with the sealing tape.			
<b>Incubation <u>with Shaking</u>: 30 Minutes at 20°C - 25°C, 350 rpm</b>			
5 x 300 µL	Aspirate the contents of the wells and <b>wash</b> 5 x with 300 µL each Washing Buffer/ <b>well</b>		In each well

100 µL	TMB-Substrate Solution	In each well
<b>Incubation: 15 Minutes in the Dark at 20°C - 25°C</b>		
100 µL	Stopping Solution	In each well
Measure the absorbance within <b>30 min at 450 nm</b> with $\geq 590$ nm as reference wavelength		

#### QUALITY CONTROL

Good laboratory practice requires that controls are included in each assay. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance. The test results are only valid if the test has been performed following the instructions. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state or local standards/laws. All standards and kit controls must be found within the acceptable ranges as stated on the QC Certificate. If the criteria are not met, the run is not valid and should be repeated. Each laboratory should use known samples as further controls.

#### Quality criteria

For the evaluation of the assay it is required that the absorbance values of the Calibrator 0 should be below 0.25, and the absorbance of Calibrator 6 should be above 1.00. Samples, which yield higher absorbance values than Calibrator 6, should be re-tested with a higher dilution.

#### EVALUATION OF RESULTS

##### Establishing of the calibration curve

The Calibrator provided contain the following concentrations of hIGFBP-1

Standard	0	1	2	3	4	5	6
ng/ml	0	0.1	0.5	1	2	4	8
OD (450-620 nm)	0.012	0.043	0.189	0.381	0.771	1.596	2.748

- 1) Plot the calibrator concentrations on the x-axis versus the mean value of the absorbance of the calibrator on the y-axis.
- 2) Recommendation: Calculation of the calibration curve should be done by using a computer program, because the curve is in general (without respective transformation) not ideally described by linear regression. **A higher-grade polynomial, or four parametric logistic (4-PL) curve fit or non-linear regression** are usually suitable for the evaluation (as might be spline or point-to-point alignment in individual cases).
- 3) The IGFBP-1 concentration in ng/mL of the samples can be calculated by **multiplication** with the respective **dilution factor**.

#### Example of a typical calibration curve

The following data is for demonstration only and cannot be used in place of data generation at the time of assay. The exemplary shown calibration curve in Figure 1 **cannot** be used for calculation of your test results. You have to establish a calibration curve for each test you conduct!


Figure 1. Exemplary calibration curve

### Exemplary calculation of IGFBP-1 concentrations

Sample dilution: 1:16

Measured extinction of your sample: 0.199

Your **measurement program** will calculate the IGFBP-1 concentration of the diluted sample automatically. You only have to determine the most suitable curve fit (here: polynomial 3<sup>rd</sup> degree). In this exemplary case the following equation is solved by the program to calculate the IGFBP-1 concentration in the sample:

$$y = -0.00269x^3 + 0.0205x^2 + 0,351x + 0.0097$$

$$0.199 = -0.00269x^3 + 0.0205x^2 + 0.351x + 0.0097$$

$$x = 0.5235$$

If the dilution factor (**1:16**) is taken into account the IGFBP-1 concentration of the undiluted sample is  $0.5235 \text{ ng/mL} \times 16 = 8.376 \text{ ng/mL}$

### Interpretation of results

The test results should not be the only base for therapeutic decisions. The results should be interpreted in regard to anamnesis, further clinical observations and results of other diagnostic investigations. Further, it is recommended to establish reference and cut-off values corresponding to the relevant group of patients for each laboratory.

### LIMITATION OF PROCEDURE

The DAsource IGFBP-1 ELISA kit, is based on murine antibodies. Generally, this technique is sensible to heterophilic antibodies as well as to human anti mouse antibodies in the sample. The influence of these antibodies is reduced by assay design, but cannot be excluded completely. Further, for interpretation of IGFBP-1 concentrations diagnostic sensitivity and specificity must be taken into account.

Also test duration has to be considered, this test system is not suitable for point-of care.

### REFERENCE VALUES

Concentrations of IGFBP-1 in human sera of 69 healthy adult donors were determined with the DAsource ELISA kit. Slight gender dependent differences were found, the concentrations of all samples varied from minimal 0.23 ng/ml to maximal 17.94 ng/ml (see table 1).

**Table 1.** IGFBP-1 Expectation values in sera of healthy adults (measured values in ng/ml)

Gender	No. of Samples	Average value	Median	Min. – Max.
female	33	4.79	4.24	0.23 – 16.07
male	36	5.22	2.71	0.42 – 17.94
total	69	5.01	2,77	0.23 – 17.94

### PERFORMANCE CHARACTERISTICS

#### Sensitivity

The analytical sensitivity was determined by measuring the null standard (calibrator 0) and calculating the corresponding concentration of the signal intensity of the calibrator 0 + 2SD. In three different assays analytical sensitivities from 0.03 to 0.08 ng/mL were measured with an average value of 0.055 ng/mL.

#### Specificity

The IGF/IGFBP-System consists of several related and homologous proteins. Thus the cross-reactivity of the DAsource IGFBP-1 ELISA kit to other IGFBPs was tested exemplarily with the two IGFBPs most abundant in circulation. Table 2 summarizes the results.

**Table 2.** Cross-reactivity to related IGFBPs. 500 ng of the indicated IGFBP were added to 1mL dilution buffer and measured as a sample. Shown is the concentration measured and the relative cross reactivity.

	Concentration [ng/ml]	Relative cross-reactivity [%]
IGFBP-2	0.0071	0.00142
IGFBP-3	0.0033	0.00066

## Reproducibility and Precision

### Intra-Assay-Variation

Serum samples were measured up to 20 times in the same assay and the coefficient was calculated based on the recalculated IGFBP-1 concentration. Exemplary data are shown in Table 3. The mean measured coefficient of variation (CV) is 6.52% (n =6).

**Table 3** Intra-Assay Variation

	Sample 1	Sample 2	Sample 3
Mean [ng/mL]	4.25	12.21	55.17
SD [ng/mL]	0.24	0.75	2.38
CV [%]	5.59	6.11	4.32
Number [n]	20	20	20

### Inter-Assay Variance


Serum samples were measured in independent assays. On average the coefficient of variation was 6.05% (SD 0.46). Exemplary results are shown in table 4.

**Table 4:** Inter-Assay Variation

	Sample 1	Sample 2	Sample 3	Sample 4	Sample 5	Sample 6
Mean [ng/mL]	12.50	8.52	6.42	18.62	4.43	9.01
SD [ng/mL]	0.17	0.24	0.25	0.52	0.18	0.18
CV [%]	1.33	2.80	3.90	2.80	4.11	2.05
Number [n]	5	5	5	5	5	5

## Linearity

The linearity of sample dilution is shown in Figure 1. Exemplarily three serum samples were diluted 1:5 up to 1:512 and the IGFBP-1 concentration was measured within each dilution. The results were analysed by linear regression analysis and revealed coefficients of determination of >0.98 for each of the samples.

**Figure 2.** Linearity of the sample dilution

## Recovery

Native IGFBP-1 (4.5 ng/mL) was added to human serum and quantified in comparison to the same amount of IGFBP-1 in buffer. The mean recovery detected was 82%. Also, the IGFBP-1 content of samples enriched with recombinant IGFBP-1 was measured and recovery calculated in comparison to enriched buffer (45 ng/mL) with a mean recovery of 94%.

## Trueness / Assay calibration

There is no IGFBP-1 reference preparation available. The DIAsource IGFBP-1 ELISA kit was calibrated against purified, native human IGFBP-1. Traceability / lot-to-lot consistency is accomplished by a panel of different serum samples.

## Interference

Interference of haemoglobin, bilirubin and triglycerides was tested by adding the indicated amount of these substances to human serum containing IGFBP-1. For comparison the same amount of buffer without any substance was also added to the serum and the relative recovery rate was calculated. Table 5 shows that on average neither haemoglobin, bilirubin nor

triglycerides exert significant influence on the measurement of IGFBP-1 in human serum. But generally measurement of IGFBP-1 in haemolytic, icteric or lipaemic samples should be avoided.

**Table 5.** Interference of haemoglobin, bilirubin and triglycerides on IGFBP-1 measurement in human serum. Shown is the relative recovery of IGFBP-1 measured in human serum samples after adding the indicated amount of possibly interfering substances in comparison to serum samples only buffer was added to.

Recovery [%]	Triglyceride 100 mg/mL	Bilirubin 100 µg/mL	Haemoglobin 5 mg/mL
Sample 1	142	72	122
Sample 2	125	46	119
Sample 3	82	113	97
Mean	116	77	113

## ASSAY COMPARISON

The DAsource IGFBP-1 ELISA kit was compared with two competitor assays. The correlation of these test systems and the DAsource assay was evaluated by Passing-Bablok regression and linear regression. The results are shown in Table 6.

**Table 6.** Assay comparison  
RSD.= residual standard deviation

	Sample [n]	Passing Bablok Regression	Linear Regression coefficient of regression
Competitor A	36	$y = -0.99 + 0.3x$ (RSD 2.97)	0.89
Competitor B	34	$y = 0.1 + 1.17x$ (RSD 0.84)	0.84

The analysis reveals that the DAsource test correlates well with both of the competitive tests. But regarding the absolute concentrations measured, there is a significant difference especially in comparison to Competitor A. This indicates differences in calibration or binding properties of the antibodies used.

## Instruction for scientific application

### SCIENTIFIC APPLICATION

This enzyme immunoassay kit is suited for measuring IGFBP-1 in human serum or Heparin and EDTA plasma or in other body fluids, for example amnion fluid, mother's milk, urine or saliva, as for scientific purposes. It is also suited to quantitate IGFBP-1 in cell culture media.

### Samples suitable for scientific application

Where required, depending on the expected IGFBP-1-values, the dilution with Dilution Buffer can be higher or lower (at least however 1:2.5). The IGFBP-1 concentration may be completely variable in different body fluids. Examples as well as dilution recommendations are given in table 7.

**Table 7** Results of matrix tests. Purified native IGFBP-1 was added to the respectively diluted samples. Enriched samples were measured without further dilution. Shown is the relative recovery of added IGFBP-1 of the value measured in enriched Dilution Buffer.

Samples	Expected IGFBP-1 [ng/ml]	Recovery IGFBP-1 [%]	Recommended Dilution
Amniotic Fluid	8.140 - 16.450	n.d.	≥1:5000 up to 1:25000
Mother's milk	5.12	96	1:10
Urine	0.07	90	1:2.5
Saliva	< 0.02	63	≥1:2.5
Bronchial Lavage	< 0.02	100	1:2.5
Sputum	< 0.02	100	1:20
Cell Culture Media	---	94	≥1:5

n.d.= not determined

### Species cross reactivity

Serum of the cited species was diluted and used as sample in this assay system.

**No signal** was detected in serum of the following species:

Rat, mouse, guinea pig, rabbit, pig, sheep, goat, cattle, horse, dog, cat, chicken or donkey.

## REFERENCES

1. Rechler, M.M., *Insulin-like growth factor binding proteins*. Vitam Horm, 1993. **47**: p. 1-114.
2. Allander, S.V., et al., *Structure and chromosomal localization of human insulin-like growth factor-binding protein genes*. Growth Regul, 1993. **3**(1): p. 3-5.
3. Rutanen, E.M., et al., *Synthesis of placental protein 12 by human decidua*. Endocrinology, 1985. **116**(4): p. 1304-9.
4. Julkunen, M., et al., *Identification by hybridization histochemistry of human endometrial cells expressing mRNAs encoding a uterine beta-lactoglobulin homologue and insulin-like growth factor-binding protein-1*. Mol Endocrinol, 1990. **4**(5): p. 700-7.
5. Khosravi, J., et al., *Immunoassay of serine-phosphorylated isoform of insulin-like growth factor (IGF) binding protein (IGFBP)-1*. Clin Biochem, 2007. **40**(1-2): p. 86-93.
6. Westwood, M., *Role of insulin-like growth factor binding protein 1 in human pregnancy*. Rev Reprod, 1999. **4**(3): p. 160-7.
7. Rutanen, E.M., H. Bohn, and M. SeVPala, *Radioimmunoassay of placental protein 12: levels in amniotic fluid, cord blood, and serum of healthy adults, pregnant women, and patients with trophoblastic disease*. Am J Obstet Gynecol, 1982. **144**(4): p. 460-3.
8. Drop, S.L., et al., *Immunoassay of a somatomedin-binding protein from human amniotic fluid: levels in fetal, neonatal, and adult sera*. J Clin Endocrinol Metab, 1984. **59**(5): p. 908-15.
9. Hall, K., G. Lundin, and G. Pova, *Serum levels of the low molecular weight form of insulin-like growth factor binding protein in healthy subjects and patients with growth hormone deficiency, acromegaly and anorexia nervosa*. Acta Endocrinol (Copenh), 1988. **118**(3): p. 321-6.
10. Busby, W.H., D.K. Snyder, and D.R. Clemmons, *Radioimmunoassay of a 26,000-dalton plasma insulin-like growth factor-binding protein: control by nutritional variables*. J Clin Endocrinol Metab, 1988. **67**(6): p. 1225-30.
11. Brismar, K., et al., *Insulin regulates the 35 kDa IGF binding protein in patients with diabetes mellitus*. J Endocrinol Invest, 1988. **11**(8): p. 599-602.
12. Suikkari, A.M., et al., *Prolonged exercise increases serum insulin-like growth factor-binding protein concentrations*. J Clin Endocrinol Metab, 1989. **68**(1): p. 141-4.
13. Maddux, B.A., et al., *IGF-binding protein-1 levels are related to insulin-mediated glucose disposal and are a potential serum marker of insulin resistance*. Diabetes Care, 2006. **29**(7): p. 1535-7.
14. Miell, J.P., et al., *The maternal insulin-like growth factor (IGF) and IGF-binding protein response to trisomic pregnancy during the first trimester: a possible diagnostic tool for trisomy 18 pregnancies*. J Clin Endocrinol Metab, 1997. **82**(1): p. 287-92.
15. Than, G.N., et al., *Levels of placenta-specific tissue protein 12 (VP12) in serum during normal pregnancy and in patients with trophoblastic tumour*. Arch Gynecol, 1983. **234**(1): p. 39-46.
16. Rutanen, E.M., F. Pekonen, and T. Karkkainen, *Measurement of insulin-like growth factor binding protein-1 in cervical/vaginal secretions: comparison with the ROM-check Membrane Immunoassay in the diagnosis of ruptured fetal membranes*. Clin Chim Acta, 1993. **214**(1): p. 73-81.
17. Martinez de Tejada, B., et al., *Can we improve the diagnosis of rupture of membranes? The value of insulin-like growth factor binding protein-1*. Bjog, 2006. **113**(9): p. 1096-9.

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### REFERENCIAS

1. Rechler, M.M., *Insulin-like growth factor binding proteins*. Vitam Horm, 1993. **47**: p. 1-114.
2. Allander, S.V., et al., *Structure and chromosomal localization of human insulin-like growth factor-binding protein genes*. Growth Regul, 1993. **3**(1): p. 3-5.
3. Rutanen, E.M., et al., *Synthesis of placental protein 12 by human decidua*. Endocrinology, 1985. **116**(4): p. 1304-9.
4. Julkunen, M., et al., *Identification by hybridization histochemistry of human endometrial cells expressing mRNAs encoding a uterine beta-lactoglobulin homologue and insulin-like growth factor-binding protein-1*. Mol Endocrinol, 1990. **4**(5): p. 700-7.
5. Khosravi, J., et al., *Immunoassay of serine-phosphorylated isoform of insulin-like growth factor (IGF) binding protein (IGFBP)-1*. Clin Biochem, 2007. **40**(1-2): p. 86-93.
6. Westwood, M., *Role of insulin-like growth factor binding protein 1 in human pregnancy*. Rev Reprod, 1999. **4**(3): p. 160-7.
7. Rutanen, E.M., H. Bohn, and M. SeVPala, *Radioimmunoassay of placental protein 12: levels in amniotic fluid, cord blood, and serum of healthy adults, pregnant women, and patients with trophoblastic disease*. Am J Obstet Gynecol, 1982. **144**(4): p. 460-3.
8. Drop, S.L., et al., *Immunoassay of a somatomedin-binding protein from human amniotic fluid: levels in fetal, neonatal, and adult sera*. J Clin Endocrinol Metab, 1984. **59**(5): p. 908-15.
9. Hall, K., G. Lundin, and G. Pova, *Serum levels of the low molecular weight form of insulin-like growth factor binding protein in healthy subjects and patients with growth hormone deficiency, acromegaly and anorexia nervosa*. Acta Endocrinol (Copenh), 1988. **118**(3): p. 321-6.
10. Busby, W.H., D.K. Snyder, and D.R. Clemmons, *Radioimmunoassay of a 26,000-dalton plasma insulin-like growth factor-binding protein: control by nutritional variables*. J Clin Endocrinol Metab, 1988. **67**(6): p. 1225-30.
11. Brismar, K., et al., *Insulin regulates the 35 kDa IGF binding protein in patients with diabetes mellitus*. J Endocrinol Invest, 1988. **11**(8): p. 599-602.
12. Suikkari, A.M., et al., *Prolonged exercise increases serum insulin-like growth factor-binding protein concentrations*. J Clin Endocrinol Metab, 1989. **68**(1): p. 141-4.
13. Maddux, B.A., et al., *IGF-binding protein-1 levels are related to insulin-mediated glucose disposal and are a potential serum marker of insulin resistance*. Diabetes Care, 2006. **29**(7): p. 1535-7.
14. Miell, J.P., et al., *The maternal insulin-like growth factor (IGF) and IGF-binding protein response to trisomic pregnancy during the first trimester: a possible diagnostic tool for trisomy 18 pregnancies*. J Clin Endocrinol Metab, 1997. **82**(1): p. 287-92.

15. Than, G.N., et al., *Levels of placenta-specific tissue protein 12 (VP12) in serum during normal pregnancy and in patients with trophoblastic tumour*. Arch Gynecol, 1983. **234**(1): p. 39-46.
16. Rutanen, E.M., F. Pekonen, and T. Karkkainen, *Measurement of insulin-like growth factor binding protein-1 in cervical/vaginal secretions: comparison with the ROM-check Membrane Immunoassay in the diagnosis of ruptured fetal membranes*. Clin Chim Acta, 1993. **214**(1): p. 73-81.
17. Martinez de Tejada, B., et al., *Can we improve the diagnosis of rupture of membranes? The value of insulin-like growth factor binding protein-1*. Bjog, 2006. **113**(9): p. 1096-9.

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