Angiotensin I RIA

KIPB3518

History

Summary of change :

Previous Version :	Current Version :
150226/5	200224/1
	Addition of the following sentence at the end of the English IFU: "Other translations of this Instructions for Use can be downloaded from our website: https://www.diasource-diagnostics.com/ "
No IVD symbol	IVD symbol added
LOT :150226/5	Version: 200224/1
PI number	No PI number
No manufacturer symbol	Manufacturer symbol added

Angiotensin I RIA

en

Radioimmunoassay of Angiotensin I for the *in vitro* determination of Plasma Renin Activity (PRA) in human plasma

KIPB3518

IN VITRO DIAGNOSTIC

DIAsource ImmunoAssays SA - Rue du Bosquet 2, B-1348 Louvain-la-Neuve, Belgium - Tel: +32 10 84 99 11 - Fax : +32 10 84 99 90

1. PRINCIPLE OF THE ASSAY

The Angiotensin I RIA serves for the quantitative determination of plasma renin activity (PRA) by the radioimmunoassay of the product of the reaction, angiotensin I.

The generation of angiotensin I is the result of the enzymatic cleavage of the renin substrate, angiotensinogen, in plasma samples in the presence of ACE inhibitor (ACE - Angiotensin-Converting Enzyme), an enzymatic inhibitor that blocks the conversion of angiotensin I to angiotensin II.

The immunoassay of angiotensin I is a radioimmunological competition assay. Unknown samples, control, and calibrators are incubated in polyclonal antibody-coated tubes with ¹²⁵I-labeled angiotensin I as tracer. After incubation, the contents of the tubes are aspirated. The bound radioactivity is then determined in a gamma counter. A standard curve is established and unknown values are determined by interpolation from the standard curve.

2. REAGENTS PROVIDED

All reagents of the kit are stable until the expiry date indicated on the kit label, if stored at 2-8°C. Storage conditions for reagents after reconstitution or dilution are indicated in paragraph Assay Procedure. Expiry dates printed on component vial labels apply to the long-term storage by manufacturer only, prior to assembling of the kit. Do not take into account.

Anti-angiotensin polyclonal antibody-coated tubes:
2 x 50 tubes (ready-to-use)

¹²⁵I-labeled angiotensin I: one 11 ml vial (ready-to use)

The vial contains 260 kBq, at the date of manufacture, of ¹²⁵I-labeled angiotensin I in buffer with bovine serum albumin, Pro Clin 300 (<0.06%; see Precautions) and a dye.

Calibrators: six x 1 ml vials (ready-to-use)

The calibrator vials contain from 0 to approximately 30 ng/mL of angiotensin I in buffer with bovine serum albumin and sodium azide (<0.1%; see § Precautions). The exact concentration is indicated on each vial label. The calibrators were calibrated against RP 86/536.

Control sample: one 1 ml vial (ready-to-use)

The vial contains angiotensin I in buffer with bovine serum albumin and sodium azide (<0.1%; see § Precautions). The expected value is in the concentration range indicated in a supplement.

Enzymatic inhibitor; one vial (lyophilized)

It also contains sodium azide (<0.1%, see § Precautions)

Wash Solution (20x): one 50 ml vial

Concentrated solution has to be diluted before use

Ag	125I
----	------

CAL	N
-----	---

CONTROL

NEUTR	SOLN
-------	------

WASH	SOLN	CONC
------	------	------

3. MATERIAL REQUIRED BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- precision micropipets (75 µL; 100 µL).
- adjustable dispensers (200 µL; 300 µL; 2 mL).
- water bath.
- ice bath.
- vortex-type mixer.
- horizontal or orbital shaker.
- aspiration system.
- gamma counter set for ¹²⁵I.

4. PRECAUTIONS

4.1 General remarks:

- Enzymatic inhibitor solution, calibrators, control sample and analyzed samples must be cooled to 2-8°C before pipeting.
- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- Do not mix the reagents from kits of different lots.
- A calibration curve must be established with each assay.
- It is recommended to perform the immunoassay in duplicate.
- Each tube must be used only once

4.2 Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material is subject to the regulations of the country of use.

Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipeting of radioactive solutions by mouth.
- Avoid all contact with radioactive materials by using gloves and laboratory overalls.
- All manipulation of radioactive substances should be done in an appropriate place, distant from corridors and other busy places.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.
- This kit contains ¹²⁵I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

4.3 Sodium azide

Some reagents contain sodium azide as a preservative. Sodium azide may react with lead, copper or brass to form explosive metal azides. Dispose of the reagents by flushing with large amounts of water through the plumbing system.

4.4 ProClin 300

R43 may cause sensitisation by skin contact

4.5 Materials of human origin

All plasma samples should be handled as if capable of transmitting hepatitis or AIDS and waste should be discarded according to the country rules.

5. SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

Plasma samples have to be collected into cold EDTA tubes

- Separate plasma from cells by centrifugation at 2-8°C.
- Keep plasma samples frozen (<-20°C, 1 year maximum) if determination is not to be performed immediately, after aliquoting in order to avoid repeated freezing and thawing.

Note: The temperature of plasma samples must be kept at 2-8°C in the course of sampling. Avoid further manipulation to prevent both formation and decomposition of angiotensin I.

6. ASSAY PROCEDURE

6.1 Preparation and storage of reagents

6.1.1 Preparation of enzymatic inhibitor solution

The content of the vials is reconstituted with the volume of cold distilled water (4°C) indicated on the label and mixed. The reconstituted enzymatic inhibitor may be stored at 2-8°C until the expiry date of the kit.

6.1.2. Preparation of wash solution

Pour the content of the vial into 950 mL of distilled water and homogenize. The diluted solution may be stored at 2-8°C until the expiry date of the kit.

6.2 Enzymatic step – generation of Angiotensin I

6.2.1 Remarks and recommendations

- The enzymatic inhibitor has to be cooled to 4°C before addition to the sample.
- Both incubation temperatures (4°C and 37°C) must be adhered to strictly, even slight variations may cause severe errors in determination.
- The enzymatic incubation time at 37°C should be determined as precisely as possible and kept within narrow limits for the whole set of tubes.
- The promptness of the temperature increase from 4°C to 37°C and the following reverse drop are critical. A circulating water bath is convenient for warming and, the use of an iced-cooled water bath is advisable for cooling.
- The promptness of the temperature increase and drop may be improved by using tubes made of material with good thermal conductivity (glass).
- If low plasma renin activity of the sample is expected, the incubation time of the enzymatic step may be prolonged for up to 3 hours.

6.2.2 Enzymatic step – procedure

Attention: Do not treat the calibrators and the control sample.

- Add 200 µL of pre-cooled enzymatic inhibitor to 200 µL of each plasma sample and mix.
- Split each sample into two 200 µL aliquots.
- Place the first aliquot into an ice-cold water bath in a refrigerator (intended for the determination of background angiotensin I at 4°C).
- Place the second one into the water bath set for 37°C (intended for the determination of generated angiotensin I at 37°C).
- Incubate all aliquots for 1 hour.
- After incubation, cool samples from 37°C to 4°C rapidly using ice water bath

6.3 Immunoassay procedure

Step 1 Additions *	Step 2 Incubation	Step 3 Counting
To antibody coated tubes, add successively: - 75 µL of calibrator, control or sample after enzymatic incubation at 37°C and at 4°C respectively and - 100 µL of tracer.** Mix.	Incubate 2 hours at 18-25°C with shaking (> 280 rpm).	Aspirate carefully the contents of tubes (except the 2 tubes "total cpm"). Wash with 2 mL of wash solution. Aspirate twice. Determine activity (cpm) for 1 min.

* Calibrators, control sample and analyzed samples have to be cooled to 4°C before pipeting. Mix samples gently before they are added.

** Add 100 µL of tracer to 2 additional tubes to obtain total cpm.

7. RESULTS

Results are obtained from the calibration curve by interpolation. The curve serves for the determination of angiotensin I concentrations in samples measured at the same time as the calibrators.

7.1 Calibration curve

The results in the package insert were calculated using a logit-log curve fit (weighted cubic regression) with B/T (%) or B/B0 (%) on vertical axis and the angiotensin I concentrations of calibrators on the horizontal axis (ng/mL). Other data reduction methods may give slightly different results.

Total activity : 68 511 cpm				
Calibrators	Angiotensin I (ng/mL)	cpm (n=3)	B/T (%)	B/B0 (%)
0	0.00	17 444	25.5	100
1	0.30	13 732	20.0	78.7
2	1.00	9 390	13.7	53.8
3	3.00	5 431	7.93	31.1
4	10.0	2 746	4.01	15.7
5	30.0	1 243	1.81	7.13

(Example of calibration curve, do not use for calculation)

7.2 Samples

For the control and samples incubated at 4°C, or at 37°C, locate the B/T (%) or the B/B0 (%) value on the vertical axis and read off the corresponding angiotensin I concentration in ng/mL on the horizontal axis.

7.3 Calculation of plasma renin activity

The determination of plasma renin activity is performed indirectly by the measurement of the in vitro generation of angiotensin I (A-I) per hour. Background A-I, determined on plasma samples incubated at 4°C, is subtracted from the A-I generated at 37°C for the calculation of PRA using the following equation:

$$\text{PRA ng of A-I /mL/hr} = \frac{[A-I(37^\circ\text{C}) - A-I(4^\circ\text{C})] \times 2}{\text{Enzymatic incubation time (hrs)}}$$

Where

A-I (37°C): angiotensin concentration in ng/mL of sample incubated at 37°C
A-I (4°C): angiotensin concentration in ng/mL of sample incubated at 4°C

8. QUALITY CONTROL

Good laboratory practices imply that control samples must be used regularly to ensure the quality of the results obtained. These samples must be processed exactly in the same way as the assay samples, and it is recommended to analyse their results using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address:

tech.support@diasource.be

9. EXPECTED VALUES

It is suggested that each laboratory establishes its own normal values. The PRA values presented below are indicative only.

N	Normal adult	2.5th - 97.5th percentile (ng/mL/hr)	Median (ng/mL/hr)	Min-Max (ng/mL/hr)
38	Early Morning, Supine	0.32-1.84	0.79	0.30-1.90
41	Upright, 2 Hours	0.60-4.18	2.20	0.48-4.88

10. PERFORMANCE CHARACTERISTICS

(For more details, see the data sheet "APPENDIX")

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary

10.1 Sensitivity

10.1.1 Analytical sensitivity: 0.07 ng/mL

10.1.2 Functional sensitivity: 0.20ng/mL

10.2 Specificity

The antibody used in the immunoassay is highly specific for angiotensin I. Extremely low cross reactivity was obtained against angiotensin II.

Moreover, the influence of possible interferences on PRA result is eliminated by subtraction of background.

10.3 Precision

10.3.1 Intra-assay

Samples were assayed in at least 25 replicates in the same series. The coefficients of variation were found below or equal to 11.3%.

10.3.2 Inter-assay

Samples were assayed in duplicate in 10 different series. Coefficients of variation were found below or equal to 20.9%.

10.4 Accuracy

10.4.1 Dependence on time of enzymatic incubation

The samples were incubated with enzymatic inhibitor for 60, 120, and 180 minutes. No significant effect on PRA results was found.

10.4.2 Dilution test

Plasma samples were serially diluted in the zero calibrator. The recovery percentages were obtained between 78% and 99%.

10.4.3 Recovery test

Plasma samples were spiked with known quantities of angiotensin I. The recovery percentages were obtained between 104% and 123%.

10.5 Measurement range (from analytical sensitivity to highest calibrator):

0.07 to approximately 30 ng/mL.

11. LIMITATIONS OF THE METHOD

The non-respect of the instructions in this package insert may affect results significantly.

Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

Do not use hemolyzed, lipemic or icteric samples.

For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays.

Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

Revision date : 2020-02-24

Other translations of this instructions for use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

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