# TRYPSIN RIA

*KIPCE07*

# History

---

## Summary of change :

<b>Previous Version :</b>	<b>Current Version :</b>
190719/1	200224/1
Multilanguage IFU	Addition of the following sentence at the end of the English IFU: "Other translations of this Instructions for Use can be downloaded from our website: <a href="https://www.diasource-diagnostics.com/">https://www.diasource-diagnostics.com/</a> "
No Manufacturer symbol	Manufacturer symbol added
No IVD symbol	IVD symbol added
"LOT"	"Version"
PI number	PI number cleared

# TRYPSIN RIA

en

KIPCE07

IN VITRO DIAGNOSTIC USE

DIAsource ImmunoAssays SA - Rue du Bosquet 2, B-1348 Louvain-la-Neuve, Belgium - Tel: +32 10 84 99 11 - Fax : +32 10 84 99 90

## 1. INTENDED USE

Radioimmunoassay for the in vitro quantitative measurement of human trypsin in serum.

## 2. INTRODUCTION

Differential diagnosis in patients presenting with "acute abdomen" is often difficult, particularly when the symptoms are due to pancreatic disease. Some means of assessing pancreatic function in these circumstances is desirable, and tests should also be suitable for monitoring the progress of the condition. In searching for a diagnostic marker, it is necessary to examine the characteristic products of the exocrine pancreas. In practice, these are the various pancreatic enzymes- principally, amylase, lipase and trypsin which are involved in the intestinal digestion of carbohydrates, fats and proteins. Of these enzymes, only trypsin is produced solely by the pancreas and its measurement may thus be considered to provide an ideal basis for a specific test of pancreatic function.

### 2.1 Chemical and biological properties of trypsin

Trypsin is one of the enzymes involved in the digestion of dietary protein in the small bowel. It functions as an endopeptidase with a marked preferential specificity for peptide bonds arising from the carbonyl groups of L-arginyl and L-lysyl residues. Chemically, trypsin is a protein consisting of 201 amino acids with a molecular weight of 22,900.

In common with other proteolytic enzymes of pancreatic origin, trypsin is synthesized and secreted from the acinar cells of the exocrine pancreas as an inactive precursor or proenzyme, trypsinogen. This occurs in two forms, both of which are converted to enzymatically active trypsin by the action of enterokinase produced by cells in the duodenum. This kinase, in the presence of  $Ca^{++}$  ions, removes the hexapeptide Val-(Asp)<sub>4</sub>-Lys from trypsinogen, thereby unmasking the active centre of the trypsin molecule. In addition trypsin itself can activate the trypsinogen by autocatalysis.

Some trypsin appears to be excreted intact from the gut, as it can be detected and measured in stool. In addition, the trypsinogen, trypsin and probably trypsin bound to a pancreatic inhibitor may find their way from the pancreatic cells into the blood via the interstitial fluid. In blood all tryptic activity is blocked by three specific inhibitors ( $\alpha_1$ -antitrypsin,  $\alpha_2$ -macroglobulin and inter- $\alpha$ -trypsin inhibitor). Nevertheless, normal serum still retains some proteolytic activity as shown by the hydrolysis of synthetic substrates, but this cannot be due solely to trypsin.

## 3. PRINCIPLE OF THE TEST

Trypsin can be measured in duodenal or pancreatic juice by spectrophotometric methods, which depend upon its ability to hydrolyse synthetic substrates. However, as the trypsin is inactive and the enzymatic activity of trypsin in serum is blocked by the inhibitors mentioned above, this approach cannot be used to determine blood levels of the enzyme. If the immunological rather than the enzymatic properties of trypsin could be used in a method for its determination in serum, the foregoing problem would not arise.

The Trypsin RIA kit permits radioimmunological determination of the enzyme trypsin in human serum and other biological fluids. The kit utilizes the principle of competitive protein binding analysis, using a double antibody radioimmunoassay method, free trypsin being separated from antibody-bound trypsin by an anti-rabbit precipitating serum.

During the first incubation, serum trypsin and <sup>125</sup>I-labelled trypsin compete for highly specific antibodies against trypsin, with the result that <sup>125</sup>I-labelled trypsin is bound to the antibodies in inverse proportion to the amount of serum trypsin present. During the second incubation the trypsin or <sup>125</sup>I-trypsin antibody complex reacts with an anti-rabbit-gammaglobulin, causing precipitation. The precipitate can be centrifuged down, the supernatant decanted and the iodine-125 activity in the precipitate measured in a gamma scintillation counter. The activity measured is compared with a calibrator curve prepared under the same conditions and provides an indication of the trypsin concentration in the serum.

The reagents supplied in one Trypsin RIA kit are sufficient for preparation of a calibrator curve and measurement of 42 patients' sera, in duplicate. The trypsin control serum included in the kit provides an analytical means of checking the accuracy of the assay procedure performed in the laboratory.

The assay has a range of 0-1,300 ng/mL serum, with greatest precision in the range 80-500 ng/mL. The reagents in the kit are ready for use, except the calibrators and the wash reagent, and should be stored at 2 - 8°C.

The Trypsin RIA radioimmunoassay kit permits measurement of trypsin as a protein in the serum, and is not affected by the presence of serum inhibitors. The radioimmunoassay is specific for trypsin, trypsinogen and the enzyme's inhibited forms.

There is no internationally accepted calibrator for trypsin at present. The Trypsin RIA values cannot be readily compared with the results of other trypsin assays.

#### 4. REAGENT PROVIDED

Ag	125I	1 vial of <sup>125</sup> I-trypsin (human), < 90 kBq, 22 mL buffer with bovine albumin, rabbit IgG, sodium azide and an orange red dye.
ANTISERUM		1 vial of rabbit anti-trypsin antiserum 10.5 mL buffer with bovine albumin and sodium azide.
CAL	N	7 vials of human trypsin calibrators, per 0.5 mL freeze dried horse serum and sodium azide, concentration in the nominal range of 0-1,300 ng trypsin/mL.*
CONTROL	N	1 vial of human trypsin control serum, 0.5 mL human serum and sodium azide, concentration stated.
PREC	AGENT	1 vial of precipitation reagent, 55 mL of buffer, goat anti-rabbit IgG serum, bovine albumin and sodium azide.
WASH	SOLN	1 tube of wash reagent (1 buffer tablet).
CONC		1 instruction for use.

\* The values shown above are the target values. The real values are indicated on the label.

The dissolved reagents contain sodium azide as preservative. Avoid swallowing and contact with the skin or mucous membranes. Sodium azide may react with lead or copper piping to form highly explosive metal azides. During waste disposal, flush the drains thoroughly to prevent a build-up of these products.

#### 5. MATERIAL REQUIRED BUT NOT PROVIDED

Incubation tubes. Microlitre pipettes with disposable plastic tips (or Dispensettes) 100, 200, and 500 µL. Dispensers 1.0 mL. Measuring cylinders, beaker. Rotary mixer, centrifuge ≥ 1,500 g, gamma scintillation counter.

#### 6. REAGENTS PREPARATION

The kit components, which have been stored at 2 - 8 °C, are brought up to room temperature (RT, 18-25 °C).

##### Calibrators:

Carefully dissolve the calibrators in distilled water, proceeding as follows:

Gently tap the **vials of calibrator** to dislodge any particles which may be adhering to the stopper. Then carefully remove the stoppers, placing them upside down on the work surface. Using a suitable micropipette, dispense precisely 500 µL **distilled water** into each vial, replace the stoppers and allow approximately 10 minutes for the freeze-dried material to **dissolve completely**. Reconstitution and mixing may be accelerated by swirling the contents and/or rolling the vials between the hands.

##### Wash Buffer:

The wash buffer is prepared by dissolving the buffer tablet in 100 mL distilled water.

All the reagents should be stored at 2- 8°C.

On occasions where large numbers of samples are to be assayed, reagents have to be pooled from more than one kit bearing the sample lot number. Only one calibration curve should be employed to interpolate all samples. It is possible to prepare 2 calibration curves with the reagents provided, but a consequently smaller number of serum samples can be measured.

#### 7. SPECIMEN COLLECTION AND PREPARATION

After blood sampling, serum is obtained by the usual methods.

It is pipetted off, assayed directly or stored up to 3 days at 2 - 8°C for later use, or alternatively frozen at -20°C.

#### 8. PROCEDURE

1. Number sufficient incubation tubes (3-5 mL), as given in the Table (calibrators, control serum, serum samples and "total activity" tubes).
2. Pipette 100 µL **calibrator, control serum** and **serum samples** into test tubes that have been prepared.
3. Dispense 200 µL <sup>125</sup>I-**trypsin** into each test tube (including "total activity" tubes).
4. Dispense 100 µL **rabbit anti-trypsin serum** into each test tube (excluding "total activity" tubes), mix the contents of the tubes on the rotary mixer, and incubate for 3h (3-5h) at room temperature (18-25 °C) away from direct light
5. Dispense 500 µL **precipitation reagent** into each test tube (excluding "total activity" tubes), mix the contents of the tubes on the rotary mixer, and incubate for **30 min** (30-60 min) at room temperature away from direct light.

6. Dispense 1mL of **wash reagent** into each test tube (excluding "total activity" tubes), centrifuge the tubes at  $\geq 1,500$  g for 15 min at 20 °C and decant the supernatants.
7. Measure the test tubes for 1 minute in a gamma scintillation counter. 40,000 to 20,000 cpm are to be expected as total activity.

### 9. CALCULATION OF RESULTS


The radioactivity measured in counts per minute (cpm) of the calibrators  $S_0 - S_6$  is plotted against the relevant trypsin concentrations (ng/ml) on prepared graph paper. The "ideal" calibrator curve is constructed through these points. The mean value is calculated from the two measured values of the control serum and patients' sera and the trypsin concentration per millilitre of serum is read off from the calibrator curve.

### 10. TYPICAL DATA

The Figure shows a typical calibrator curve obtained with Trypsin RIA assay.

#### Example of a calibration curve

Tube groups	Mean cpm	B/T x 100	B/Bo x 100	Concentration ng/mL
T	35000			
Calibrator 0	22225	63.5	100	0
Calibrator 1	19713	56.3	88.7	43
Calibrator 2	17002	48.5	76.5	85
Calibrator 3	13535	38.6	60.9	171
Calibrator 4	9779	27.9	44.0	341
Calibrator 5	7445	21.2	33.5	602
Calibrator 6	5929	16.9	26.6	1365
Control	12090	34.5	54.4	220

### 11. CLINICAL RESULTS

#### 11.1. Normal levels

The function of enzymatic inactivated serum trypsin is not known, but trypsin can be detected in the serum as long as the exocrine portion of the pancreas is active.

The serum trypsin concentration of 151 presumed healthy adults has been measured, it ranges from 160 to 600 ng/mL (5<sup>th</sup> and 95<sup>th</sup> percentile respectively) with a median at 315 ng/mL.

#### 11.2. Chronic pancreatitis

The current methods available for the diagnosis of chronic pancreatitis include X-ray visualisation of calcification of the pancreas, determination of changes in amylase, lipase or bicarbonate levels in the duodenal fluid, or changes in the volume of fluid secreted. Some of these diagnostic procedures are very difficult to perform, are unreliable and unpleasant for the patient. The radioimmunoassay is therefore of great importance.

It has been shown that the serum trypsin levels of patients with chronic pancreatitis are lower than 140 ng trypsin/ml in 50 % of cases, which is probably a result of extensive destruction of the acinar cells of the pancreas.

On the other hand, in 20 % of patients with chronic pancreatitis, the serum trypsin level exceeds the normal upper limit of 400 ng/mL. This probably indicates an acute phase of the disease (see 11.3).

The severity of abdominal pain does not appear to correlate with trypsin serum levels, but steatorrhoea is frequently accompanied by a marked fall in serum trypsin.

In patients with chronic pancreatitis but "normal" trypsin basal levels, the serum trypsin increases markedly 20-120 minutes after stimulation with secretin followed by pancreozymin. This appears to be the result of increased permeability of the inflamed pancreatic tissue, since no such rise in serum trypsin is found in healthy patients undergoing the same test.

### 11.3. Acute pancreatitis

Patients with acute pancreatitis invariably have significantly increased serum trypsin concentration. The levels range between 600 and 6,500 ng trypsin/mL serum.

### 11.4. Carcinoma of pancreas

Patients with carcinoma of the pancreas may have moderately elevated serum trypsin levels (450-1,200 ng/mL) but generally there is no difference. If pancreatitis and/or renal insufficiency are excluded, other diagnostic procedures must be undertaken to confirm the diagnosis of carcinoma.

### 11.5. Cystic fibrosis of the pancreas

It is possible to use the radioimmunoassay of trypsin to diagnose cystic fibrosis of the pancreas.

#### Note:

It appears that the exocrine portion of the pancreas does not function as a unit after stimulation, as amylase, lipase and serum trypsin do not always correlate as expected. There may therefore be further possibilities for differential diagnosis. Normal serum trypsin concentrations (range 160-600 ng/mL) do not exclude the possibility of pancreatic impairment. High serum trypsin levels have also been found in renal insufficiency and the renal status of patients should therefore also be checked when interpreting trypsin levels. Bearing in mind the above limitation, trypsin levels below 160 ng/mL and above 600 ng/mL can be regarded as the lower and upper cut-off points of normality respectively when screening patients with abdominal symptoms.

## 12. PERFORMANCES AND LIMITATIONS

### 12.1. Imprecision

This was evaluated with 2 samples assayed 10 times in the same series and in 15 different series.

Within-run			Between-run		
Samples	Mean (ng/mL)	CV (%)	Samples	Mean (ng/mL)	CV (%)
1	550	2.8	3	162	4.9
2	167	3.3	4	483	6.7

### 12.2. Detection limit

The detection limit is defined as being the smallest concentration different from 0 with a confidence interval of 95 %. It has been determined as being 8.0 ng/mL.

## 13. INTERFERENCE

The presence of bilirubin at concentrations of up to 250 mg/L, hemoglobin up to 10 g/L and triglycerides up to 20 g/L have no effect on the assay results. The immuno-assay is protected against heterophilic antibodies. However, we cannot guarantee that this protection is exhaustive.

## 14. PRECAUTIONS AND WARNINGS

### Safety

For *in vitro* diagnostic use only.

This kit contains <sup>125</sup>I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

## 15. BIBLIOGRAPHY

1. Adrian TE, Bestermann HS, Mallinson CN, Galarotis C, Bloom SR. Clin Sci Mol Med. 1978;54:24.
2. Bray PT, Clark GCF, Moody GJ, Thomas JDR. Clin Chim. 1977;80:333.
3. Crossley JR, Smith PA, Elliott RB. Lancet. March 3. 1979;472.
4. Elias E, Redshaw M, Wood T. Lancet II. 1977;66.
5. Fedail SS, Salmon PR, Harvey RF, Brown PB, Read AE. Abstract. 1978;19:A445.
6. Goldberg JM. Clin Chem. 1976;22:638.
7. Lesser PB, Marshaw AL. Gastroenterology. 1975;68:630.
8. Nagasawa S, Han BH, Sugihari H, Suzuki T. J Biochem. 1970;67:821.
9. Temler RS, Felber JP. Biochim Biophys. 1976;445:720.

## 16. SUMMARY OF THE PROTOCOL

### Trypsin RIA assay:

Labelling of the test tubes	Calibrators (µL)							Control serum (µL)	Sera (µL)				Total activity (µL)
	S <sub>0</sub>	S <sub>1</sub>	S <sub>2</sub>	S <sub>3</sub>	S <sub>4</sub>	S <sub>5</sub>	S <sub>6</sub>	C	1	2	3	etc.	T
Calibrators S <sub>0</sub>	100/100	100/100	100/100	100/100	100/100	100/100	100/100						
CAL <sub>1</sub>													
CAL <sub>2</sub>													
CAL <sub>3</sub>													
CAL <sub>4</sub>													
CAL <sub>5</sub>													
CAL <sub>6</sub>													
Control serum								100/100					
Sera									100	100	100	etc.	
<sup>125</sup> I-trypsin	←-----200 µL -----→												
Anti-trypsin serum	←-----100 µL-----→												
	Mix ; incubate 3 h at 18-25 °C												
Precipitation reagent	←-----500 µL -----→												
	Mix ; incubate for 30 minutes at 18-25 °C												
Wash reagent	←-----1 mL -----→												
	Centrifuge for 15 min at ≥ 1,500 g, 20 °C, decant the supernatant												
	Measure the precipitate												
	Measure the solution												

Revision date : 2020-02-24

Other translations of this Instructions for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

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