# **Anti-IA<sub>2</sub> RIA**

***KIPM2050***

# History

---

## Summary of change :

<b>Previous Version :</b> 180223/1	<b>Current Version :</b> 200224/1
Multilanguage IFU	Addition of the following sentence at the end of the English IFU: "Other translations of this Instructions for Use can be downloaded from our website: <a href="https://www.diasource-diagnostics.com/">https://www.diasource-diagnostics.com/</a> "
No Manufacturer symbol	Manufacturer symbol added
No IVD symbol	IVD symbol added
PI number	PI number cleared

Read entire protocol before use.

## Anti-IA<sub>2</sub> RIA

### I. INTENDED USE

Radioligand assay for the determination of autoantibodies to Protein Tyrosine Phosphatase IA<sub>2</sub> in human serum.

### II. GENERAL INFORMATION

- A. **Proprietary name :** DIAsource Anti-IA<sub>2</sub> RIA Kit
- B. **Catalog number :** KIPM2050: 50 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.  
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

**For technical assistance or ordering information contact :**  
Tel : +32 (0)10 84.99.11 Fax : +32 (0)10 84.99.91

### III. CLINICAL BACKGROUND

Type 1 diabetes, also known as insulin-dependent diabetes mellitus (IDDM), results from a chronic autoimmune destruction of the insulin-secreting pancreatic beta cells, probably initiated by exposure of genetically susceptible host to an environmental agent. Autoimmune destruction of beta cells is thought to be completely asymptomatic until 80 - 90 % of the cells are lost. This process may take years to complete and may occur at any time.

During the preclinical phase, this autoimmune process is marked by circulating autoantibodies to beta cell antigens. These autoantibodies are present years before the onset of type 1 diabetes and prior to clinical symptoms. Early studies utilized the immunofluorescence test for islet-cell antibodies (ICA), which has been difficult to standardize and is now replaced by a combination of several radioimmunoassays for antibodies against specific beta cell antigens, such as insulin (IAA), glutamic acid decarboxylase (GAD) and tyrosine phosphatase ICA 512 (IA<sub>2</sub>).

IA<sub>2</sub>, a member of the protein tyrosine phosphatases family is localized in the dense granules of pancreatic beta cells and the second defined recombinant islet cell antigen. IA<sub>2</sub> shares sequence identity with the islet cell antigen 512. The higher frequency of antibodies to IA<sub>2</sub> is explained by the presence of autoantibodies directed to the COOH terminus of IA<sub>2</sub> which is lacking in the ICA512 molecule.

IA<sub>2</sub> autoantibodies are present in the majority of individuals with new-onset type 1 diabetes and in individuals in the pre-diabetic phase of the disease. The appearance of autoantibodies to IA<sub>2</sub> seems to be correlated with the rapid progression to overt type 1 diabetes.

The combination of tests for GAD65 and IA<sub>2</sub> autoantibodies is highly relevant for risk assessment of type 1 diabetes in children and adolescence. The screening for GAD65 and IA<sub>2</sub> autoantibodies detect more than 90 % of subjects at risk for type 1 diabetes and may, therefore, possess the potential to replace ICA technique.

#### IV. PRINCIPLES OF THE METHOD

The DIASource anti-IA<sub>2</sub> kit is a direct assay based on the principle of radioligand assays. Highly purified human recombinant IA<sub>2</sub> (intracellular fragment of IA<sub>2</sub>) is labeled with 125-Iodine. This tracer is used in excess and bound by the IA<sub>2</sub> autoantibodies of the sample.

The DIASource anti-IA<sub>2</sub> tracer meets the highest requirements with regard to purity, enzymologic identity, fast reaction kinetics, cross reactivity at zero level and stability. These are the main prerequisites for the specific binding of the tracer and its exclusive recognition by the IA<sub>2</sub> autoantibodies of the sample.

By adding dissolved Protein A which binds to the Fc moiety of the autoantibodies, sandwich-type complexes are formed. By adding a precipitation buffer the immune complex will be precipitated as a solid phase which facilitates the simple separation of the bound fraction (B) by centrifugation. After removing the supernatant which contains the non-bound tracer by aspiration or decantation, the radioactivity of the remaining precipitate is measured.

The concentration of IA<sub>2</sub> autoantibodies (anti-IA<sub>2</sub>) in the sample is reflected by the specifically bound tracer amount. The radioactive signal (cpm) of the bound fraction (B) is proportional to the autoantibody concentration.

No immune complex is formed if autoantibodies against IA<sub>2</sub> are absent in the sample, as the tracer binds solely to IA<sub>2</sub> autoantibodies, but not to Protein A.

A standard curve with a range of 0.01 - 60 U/ml is established by measuring cpm respectively the binding B/T % of the calibrators 1 - 5. The anti-IA<sub>2</sub> concentration value of the patient's sample is directly read off against this curve.

#### V. REAGENTS PROVIDED

Reagents	50 Tests Kit	Reconstitution		
<table border="1"> <tr> <td>Ag</td> <td>125I</td> </tr> </table> <p>TRACER: <sup>125</sup>Iodine labelled IA<sub>2</sub> (human, recombinant)</p>	Ag	125I	1 vial Lyo <0.1 MBq	<b>Reconstitute</b> with 2.6ml Reconstitution buffer
Ag	125I			
<table border="1"> <tr> <td>CAL</td> <td>N</td> </tr> </table> <p>Calibrators - N = 1 to 5 Anti- IA<sub>2</sub> calibrators in human serum (conc: see vial label)</p>	CAL	N	5 vials 0.25ml	<b>Ready</b> for use
CAL	N			
<table border="1"> <tr> <td>PROTEIN A</td> </tr> </table> <p>Protein A component</p>	PROTEIN A	1 vial Lyo	<b>Reconstitute</b> with 2.6ml Reconstitution buffer	
PROTEIN A				
<table border="1"> <tr> <td>REC</td> <td>SOLN</td> </tr> </table> <p>Reconstitution buffer</p>	REC	SOLN	1 vial 6 ml	<b>Ready</b> for use
REC	SOLN			
<table border="1"> <tr> <td>CONTROL</td> <td>N</td> </tr> </table> <p>Controls - N = 1 or 2 Anti- IA<sub>2</sub> controls in human serum (conc: see vial label)</p>	CONTROL	N	2 vials 0.25ml	<b>Ready</b> for use
CONTROL	N			
<table border="1"> <tr> <td>PREC</td> <td>AGENT</td> </tr> </table> <p>Precipitation buffer</p>	PREC	AGENT	1 vial 55ml	<b>Ready</b> for use
PREC	AGENT			

#### VI. SUPPLIES NOT PROVIDED

- Precision pipettes 20 µl, 100 µl, 1000 µl
- Disposable pipette tips
- Distilled or de-ionized water
- Absorbent paper or paper towel
- Vortexer
- Centrifuge
- γ-counter

#### VII. REAGENT PREPARATION

Allow all of the components to reach room temperature prior to use in the assay, except Precipitation buffer.

##### Tracer:

Reconstitute with 2.6 ml Reconstitution Buffer per vial. Reconstituted tracer remains stable for 2 weeks, stored at 2 - 8 °C.

##### Reconstitution Buffer:

Reconstitution Buffer is ready for use and serves for the reconstitution of the tracer and the Protein A component.

##### Precipitation Buffer:

Precipitation buffer is ready for use and serves for the precipitation and washing of the radioactive immune complex. **Use cold.**

##### Protein A component:

Reconstitute with 2.6 ml Reconstitution Buffer per vial. The reconstituted solution remains stable for 2 weeks stored at 2 - 8 °C.

##### Calibrators:

Ready for use.

##### Controls:

Ready for use.

#### VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

The DIASource anti-IA<sub>2</sub> kit has been designed for 50 determinations. This is sufficient for the analysis of 18 unknown samples as well as for calibrators and control sera, assayed in duplicates.

The expiry date of each component is reported on its respective label and that of the complete kit on the box label.

Upon receipt, all components of the DIASource anti-IA<sub>2</sub> kit have to be kept at 2 - 8 °C, preferably in the original kit box.

#### VIX. SPECIMEN COLLECTION AND PREPARATION

- Blood is taken by venipuncture. After clotting, the serum is separated by centrifugation. Plasma is also suitable for use in the DIASource anti-IA<sub>2</sub> kit. Lipaemic and hemolytic samples should not be employed.

- The samples may be kept at 2 - 8 °C up to three days. Long-term storage requires - 20 °C.

- Repeated freezing and thawing should be avoided. For multiple use, initially aliquot samples and keep at - 20 °C.

#### X. PROCEDURE

**Use conical tubes for the DIASource anti-IA<sub>2</sub> kit only.**

1. Label test tubes appropriately.
2. Pipette into the corresponding tubes according to assay scheme
  - a. 20 µl calibrators,
  - b. 20 µl control sera,
  - c. 20 µl patient's samples, each.
3. Add 50 µl tracer (prepared from Tracer and Reconstitution buffer), each, to all tubes, including those for total radioactivity .  
*Tubes are now separated until radioactivity is measured.*
4. Incubate overnight (at least 18 hours at 4 - 8 °C).
5. Add 50 µl Protein A solution (prepared from Protein A component and Reconstitution buffer), each.  
(Agitate the mixture gently prior to use - please section VII. Reagent Preparation before use).
6. Incubate 60 minutes (at 4 - 8 °C).
7. Add 1 ml cold precipitation buffer, each.
8. Centrifuge the tubes for 30 minutes at 4-8°C at 2000 x g.
9. Aspirate supernatant completely or decant. For removal of any remaining liquid, turn tubes upside down (minimum 5 minutes) and absorb any droplets by carefully tapping on blotting paper.
10. Measure radioactivity of all tubes including.  
Recommended counting time: 1 minute

## XI. CALCULATION OF RESULTS

The standard curve is established by plotting the mean cpm-values of the calibrators 1 - 5 on the ordinate, y-axis, versus their respective anti-IA<sub>2</sub>-concentrations on the abscissa, x-axis. The respective binding rates B, related to the total radioactivity T may be used as well for setting up the standard curve (B/T %).

**The use of log/log processing is recommended for best results !**

The anti-IA<sub>2</sub> concentrations of the controls and the unknown samples are directly read off in U/ml against the respective cpm or B/T values.

The DIASource anti-IA<sub>2</sub> kit may be used also with Computer Assisted Analysis using software able to plot log/log curves with spline smoothing, such as for sandwich-type assays (IRMA).

## XII. TYPICAL DATA


The following data are for illustration only and should never be used instead of the real time calibration curve.

Test tubes	cpm (a)	cpm (b)	cpm (mean)	B/T %	U/ml
Total activity T	60029	58617	59323	100 %	---
Calibrator 1	2526	2606	2566	4.3	0.01
Calibrator 2	4338	4379	4359	7.3	0.5
Calibrator 3	6621	6459	6540	11.0	2.5
Calibrator 4	14858	14822	14840	25.0	12
Calibrator 5	34776	35136	34956	58.9	60
Patient 1	8253	8245	8249	13.9	4.0

Calculation of patient sample 1:

$$\frac{B}{T} (\%) = \frac{8249}{59323} \times 100 = 13,9 \%$$

Typical standard curve:

## XIII. PERFORMANCE AND LIMITATIONS

### A. Detection limit

The most appropriate and statistically reasonable definition of the lower detection limit of any assay is at present the so-called functional assay sensitivity.

This functional assay sensitivity generally represents that concentration which corresponds to the 10 % (within-assay) and to the 20 % (between assay) coefficient of variation in the respective precision profiles of the assay in the lower concentration range. Upon correct and thorough performance of DIASource anti-IA<sub>2</sub> kit, this value is found at approx. 0.8 U/ml.

Anti-IA<sub>2</sub> values below this defined level of functional assay sensitivity do not meet the statistical criteria for reliability according to GLP (Good Laboratory Practice) and therefore cannot be distinguished from zero due to the statistically necessary certainty.

Anti-IA<sub>2</sub> concentrations above approx. 0.8 U/ml, however, fulfill these criteria and are consequently assessed as valid.

Sample	Intra-assay		Inter-assay	
	Mean (U/ml)	CV (%)	Mean (U/ml)	CV (%)
1	0.76	10.4	0.92	17.0
2	5.80	4.2	6.13	2.9
3	24.68	2.6	24.76	1.9

### B. Specificity

The high quality of the tracer (<sup>125</sup>I-IA<sub>2</sub>) ensures the exclusive reaction of anti-IA<sub>2</sub> autoantibodies in direct assay principle. Any detectable cross reaction with autoantibodies to GAD65, Insulin, thyroglobulin, thyroidal peroxidase, to the TSH receptor and Acetylcholine receptor does not exist.

### C. Calibration

The units in the DIASource anti-IA<sub>2</sub> kit are arbitrary units. DIASource anti-IA<sub>2</sub> assay based on <sup>125</sup>I-IA<sub>2</sub> shows good agreement with an established reference assay. Sensitivity and specificity shows 100 % consensus.

### D. Parallelism of standards and serum samples

Dilutions of specimen in anti-IA<sub>2</sub> free human serum are determined according to their expected theoretical values with DIASource anti-IA<sub>2</sub> M. Based on the heterogeneous nature of the autoantibody population in human serum and in view of epitope specificity and affinity of the autoantibodies, in some cases the sample dilution is not leading to the theoretically expected value.

### E. Limitations

Healthy individuals should be tested negative by using the DIASource anti-IA<sub>2</sub> kit. However, IA<sub>2</sub> autoantibodies may be also present in a rare neurological disorder, Stiff-man Syndrome (SMS). In sera from patients with SMS IA<sub>2</sub> autoantibodies appear seldom. Any clinical diagnosis should not be based on the results of in vitro diagnostic method alone. Physicians are supposed to consider all clinical and laboratory findings possible to state a diagnosis.

## XIV. REFERENCE INTERVALS

Anti-IA <sub>2</sub>	
IA <sub>2</sub> antibodies negative	< 1.0 U/ml
IA <sub>2</sub> antibodies borderline	1.0 – 2.0 U/ml
IA <sub>2</sub> antibodies positive	≥ 2.0 U/ml

It is recommended that each laboratory establishes its own normal and pathological reference ranges as usually done for other diagnostic parameters, too. Therefore, the above mentioned reference values provide a guide only to values which might be expected.

Normal range was established by evaluation of data from patients with type 1 diabetes and healthy control subjects.

## XV. PRECAUTIONS AND WARNINGS

### Safety

For *in vitro* diagnostic use only. This kit contains <sup>125</sup>I (half-life: 59 days), emitting ionizing X (27 keV) and γ (35 keV) radiations.

This kit is for *in vitro* use only. Follow the working instructions carefully. This instruction manual is valid only for the present kit with the given composition. An exchange of single components is not in agreement with CE regulations.

The expiration dates stated on the respective labels are to be observed. The same relates to the stability stated for reconstituted reagents.

All reagents should be kept at 2 - 8 °C before use in the original shipping container.

Some of the reagents contain small amounts (< 0.1 %) of sodium azide as a preservative. They must not be swallowed or allowed to come into contact with skin or mucosa. The possible formation of heavy metal azides in the drainage has to be prevented by sufficient rinsing with water.

Source materials derived from human body fluids or organs used in the preparation of this kit were tested and found negative for both Hepatitis an HIV antibody. However, no known test guarantees the absence of such viral agents. Therefore, handle all components and all patient samples as if potentially hazardous.

Since the kit contains radioactive material the following precautions should be observed:

- Do not smoke, eat or drink while handling radioactive material in any room designated for working with radioactive material,
- Always use protective gloves,
- Never pipette radioactive material by mouth,
- Wipe up spills promptly, washing the affected surface thoroughly with a decontaminant,
- Place contaminated tissues, tubes, bench covers, gloves etc. in a specially marked container, discard liquid and solid radioactive waste only as permitted by federal, state or local authorities and regulations.

It is the responsibility of the user of this product to handle radioactive material in accordance to the national rules given by law or other statements of the local authorities.

In any case GLP with all general and individual regulation has to be applied to the use of this kit.

**XVI. SUMMARY OF THE PROTOCOL**

1	Label test tubes*	Cals	Ctrl	Pat.-Sera 1, 2 etc.	T
2	Pipette Calibrators 1 - 5 Control I - II Patient sera	20 µl	20 µl	20 µl	
3	Pipette Tracer (prepared from Tracer and Reconstitution buffer)	50 µl	50 µl	50 µl	50 µl
4	Incubate**	overnight (at least 18 hours) (at 4 - 8 °C)			
5	Pipette Protein A component (prepared from Protein A component and Reconstitution buffer)	50 µl	50 µl	50 µl	
6	Incubate**	60 minutes (at 4-8 °C)			
7	Pipette Cold precipitation buffer	1 ml	1 ml	1 ml	
8	Centrifuge	30 minutes at 4-8°C at 2000 x g			
9	Decant supernatant or Aspirate supernatant	leave tubes upside down on absorbent paper for minimum 5 minutes quantitatively.			
10	Count radioactivity	Counting time: 1 minute			

\* use conical tubes

\*\* Prior to incubation, agitate the tubes briefly in order to ensure homogeneous reaction conditions.

DIAsource Catalogue Nr : KIPM2050	Revision nr : 200224/1
--------------------------------------	---------------------------

Revision date : 2020-02-24

**Other translations of this Instructions for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>**

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