# Reverse T3-RIA

*R-EW-125*

# History

---

## Summary of change:

<b>Previous Version:</b> <b>200724</b>	<b>Current Version:</b> <b>201008</b>
	<p>The following correlation was added:</p> <p><b>6. SPECIMEN COLLECTION AND HANDLING</b></p> <p>Serum or plasma (EDTA and heparin) provides similar results.</p> <p>Y (hep. Li plasma) = 0.95 x (serum) + 0.01 r = 0.95                      n = 20</p> <p>Y (hep. Na plasma) = 0.93 x (serum) + 0.004 r = 0.89                      n = 20</p> <p>Y (EDTA plasma) = 0.88 x (serum) + 0.01 r = 0.97                      n = 20</p>
<p>In the EN, DE and ES versions:</p> <p><b>7. ASSAY PROCEDURE</b></p> <p><b>3. Test procedure</b></p> <p>Centrifuge for:</p> <ul style="list-style-type: none"><li>• 20 minutes at 2000 x g;</li></ul>	<p>Correction in the EN, DE and ES versions:</p> <p><b>7. ASSAY PROCEDURE</b></p> <p><b>3. Test procedure</b></p> <p>Centrifuge for:</p> <ul style="list-style-type: none"><li>• 25 minutes at 2000 x g;</li></ul>

## FOR IN VITRO DIAGNOSTIC USE ONLY

### R-EW-125 125 Determinations

#### 1. INTENDED USE

The kit allows the quantitative determination of reverse  $T_3$  by RIA in human serum or plasma.

#### 2. EXPLANATION OF THE TEST

The presence of 3,3',5' triiodothyronine ( $rT_3$ ) in human blood has been demonstrated some years ago (1). Only 3%  $rT_3$ , structurally corresponding to a  $T_4$  molecule without iodine atom in position 5, is synthesized by thyroid gland. The rest is originated from monodeiodination of  $T_4$ , when the last one links peripheral receptors which consist of the enzymes for the production either of  $T_3$  or of  $rT_3$ . The monodeiodination occurs in many organs, but the liver and the kidney are the major sites of conversion. That is why it is possible you can detect anomalous  $rT_3$  levels owing to pathological situations not involving thyroid function (2). Because  $T_3$  is active, while  $rT_3$  is inactive, the production of  $rT_3$  instead of  $T_3$  starting from  $T_4$  can determine several metabolic modifications. In amniotic fluid  $rT_3$  levels are very high; at the birth high TSH levels determines a rapid increase of  $T_3$  and of  $T_4$  even if in small amount, but  $rT_3$  levels are unchanged for some days. This situation can be due to the small intrathyroidal storage of  $rT_3$ . For the same reason the administration of exogenous TSH does not increase  $rT_3$  levels in normal adults (3). It is very interesting to underline that in several conditions associated with low  $T_3$  values, high  $rT_3$  concentration are detectable (4). Some systemic illnesses such as hepatic cirrhosis and protein malnutrition are named «low- $T_3$  syndrome» (5). These pathologies are associated to low  $T_3$  levels, with concomitant high  $rT_3$  levels. Moreover Chopra demonstrated the same situation in newborn sera (2). This physiological  $rT_3$  increase brings back to normal only within 8-10th days of life, while all thyroidal parameters brings back to normal within the first week of life (6).  $RT_3$  profile during pregnancy has been used to reveal gestational age and fetal well-being, even if clear correlations are still unknown.

#### 3. PRINCIPLES OF THE METHOD

The kit has been designed for quantitative determination, by RIA, of reverse  $T_3$  in serum or plasma, without any preliminary treatment of the samples.

Antigen, whether from standards or from samples competes with radioactive tracer ( $^{125}I$ - $RT_3$ ) for the binding sites of antibody. After incubation, the amount of tracer bound to antibody will be inversely proportional to the amount of antigen present in standards or in samples. B/F separation is performed by addition of polyethylene glycol and the immunocomplex is collected by centrifugation.

The main steps of assay are:

- Incubation of reaction mixture for 3 hours at room temperature (18-25°C).
- Precipitation of immunocomplex by polyethyleneglycol.
- Centrifugation and decantation or aspiration.
- Counting of radioactivity.

#### 4. REAGENTS

**Do not interchange reagent from different kit lots. Lot numbers of reagent should be as stated on the "Certificate of Analysis". Do not use kit component beyond their expiration date.**

The Reverse kit contains sufficient reagents for 125 tubes. On receipt, store the kit at 2-8°C, until the expiration date printed on the label.

Reagents	125 Tests Kit	Colour Code	Reconstitution
<b>ANTISERUM</b> Antiserum: Raised in rabbit, in Tris buffer with bovine serum albumin and sodium azide, <0.1% w/w	1 Vial Lyophilized		<b>Add 12.5 ml distilled water</b>

<b>Ag</b> <sup>125I</sup> TRACER : contains $RT_3$ labelled with $^{125}I$ in Tris buffer with Sodium Azide, <0.1%	1 Vial 13ml 75 kBq	Red	<b>Ready to Use</b>
<b>CAL 0</b> Calibrator zero in phosphate buffer with bovine serum albumin and sodium azide, 0.05% w/w	1 Vial Lyophilized	Yellow	<b>Add 2 mL distilled water</b>
<b>CAL N</b> Calibrator - N = 7 in phosphate buffer with bovine serum albumin and sodium azide, 0.05% w/w	7 Vials Lyophilized	Yellow	<b>Add 1 mL distilled water</b>
<b>CONTROL N</b> Control - N = 1 Contains human serum, sodium azide	1 Vial Lyophilized	Silver	<b>Add 1 mL distilled water</b>
<b>PEG SOLN</b> PEG Solution : contains polyethyleneglycol 20% in phosphate buffer with Tween 20, 0.1% and sodium azide, <0.1% w/w.	2 Vials 62.5ml		<b>Ready to Use</b>

1.  **$RT_3$  Antiserum:** Store at 2-8°C. Reconstitute the contents of the vial with 12.5 ml distilled water and mix gently until the complete solubilization of freeze-dried residue. After reconstitution store at 4°C for 3 days or at -20°C for longer length of time until expiration of the vial label.
2.  **$RT_3$  Calibrators:** Expiry date is stated on the vial labels. Store at 2-8°C. Reconstitute the contents of calibrator zero with 2 ml distilled water and the remaining calibrators with 1 ml distilled water. Mix gently until the complete solubilization of freeze-dried residue. After reconstitution store at 4°C for 3 days or at -20°C for longer length of time until expiration of the vial label. The value of each calibrator vary from batch to batch. **For the exact value, refer to the vial label.**
3.  **$^{125}I$ -Reverse  $T_3$  :** Store at 2-8°C until expiration date printed on the vial label. The reagent is ready to use, allow the reagent to reach temperature (18-25°C) and mix thoroughly by gentle inversion avoiding foam before use.
4. **20% TW PEG Solution:** Store at 2-8°C until the expiration date on the vial label. Do not freeze.
5. **Control:** Store at 2-8°C. The expiry date is reported on vial label. Reconstitute the contents of the vial with 1 ml distilled water and mix gently until the complete solubilization of freeze-dried residue. Store at 2-8°C for 1-2 days or at -20°C for longer length of time.

#### Indications of Reagents Deterioration

1. Moisture in any of the lyophilized reagents.
2. Precipitation in any of the reconstituted reagents except the PEG Solution.

#### 5. WARNING AND PRECAUTIONS FOR USERS

**For *in vitro* diagnostic use.**

**Only experienced laboratory personnel should use this test and handling should be in agreement with GLP.**

**Radioactive Material - Not for Internal or External Use in Humans or Animals.**

This radioactive material may be received, acquired, possessed, and used only by physicians, clinical laboratories or hospitals and only for *in vitro* clinical or laboratory tests not involving internal or external administration of the material, or the radiation there from to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations of each country.

Physical characteristics of  $^{125}I$ : see end of instructions

- Safety Precautions:** The following precautions should be observed in handling radioactive material:
  - Store radioactive materials in a designated area.
  - Do not eat, drink, smoke or apply cosmetics where radioactive materials are being handled.
  - Do not pipette by mouth.
  - Wear gloves when handling radioactive materials and wash hands thoroughly afterwards.
  - Cover working area with disposable absorbent paper.
  - Wipe up all spills immediately and thoroughly and dispose of the contaminated materials as radioactive waste.
  - Dispose of the liquid radioactive waste into the sanitary sewage system if permitted by the local regulations.

- CHEMICAL HAZARD Sodium Azide (NaN<sub>3</sub>) Warning:** Some of the reagents in this kit contain sodium azide as a preservative. For all such reagents, the concentration of sodium azide is <0.1% p/p. Sodium azide may react with lead and copper plumbing to form explosive metal azides. Dispose of all non radioactive reagents by flushing with large amounts of water through the plumbing system.

**H370 very toxic if swallowed.**  
**EUH031 contact with acids liberates very toxic gas.**

- POTENTIAL BIOHAZARDOUS MATERIAL Warning:** This kit may contain some reagents made with human serum or plasma. The serum or plasma used has been tested by an FDA-approved method and found to be non-reactive for HIV-1/2 Antibodies, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, these reagents should be handled at the Biosafety Level 2 as recommended for any potentially infectious human serum or blood specimen in the Centers for Disease Control/National Institutes of Health manual "Biosafety in Microbiological and Biomedical Laboratories," 3<sup>rd</sup> Edition 1993.

## 6. SPECIMEN COLLECTION AND HANDLING

Serum or plasma samples can be used.

**Plasma:** add an anticoagulant, centrifuge for 10 minutes and collect the plasma.

**Serum:** allow the blood to clot at room temperature (18-25°C). Centrifuge for 10 minutes and collect the serum.

**Storage:** if the assay is performed within 24 hours from collection, store at 4°C; for longer storage, store at -20°C, avoiding to thaw the sample more than once; if necessary aliquot. If after thawing the sample is turbid, we suggest to centrifuge it before the assay. Grossly lipemic or badly hemolyzed samples must not be used.

**Dilutions:** the samples can be assayed undiluted. Eventual dilutions must be performed using calibrator zero supplied with the kit.

Serum or plasma (EDTA and heparin) provides similar results.

$$Y (\text{hep. Li plasma}) = 0.95 \times (\text{serum}) + 0.01 \quad r = 0.95 \quad n = 20$$

$$Y (\text{hep. Na plasma}) = 0.93 \times (\text{serum}) + 0.004 \quad r = 0.89 \quad n = 20$$

$$Y (\text{EDTA plasma}) = 0.88 \times (\text{serum}) + 0.01 \quad r = 0.97 \quad n = 20$$

## 7. ASSAY PROCEDURE

### 1. Materials Provided

The Reverse T<sub>3</sub> kit is sufficient for 125 tubes and contains the following reagents:

Reagent	Number
Reverse T <sub>3</sub> Antiserum	1
Reverse T <sub>3</sub> Calibrators	8
<sup>125</sup> I Reverse T <sub>3</sub>	1
20% Tween Peg Solution	2
Control	1

### 2. Equipment required but not supplied with the kit

- Distilled water.
- Plastic test tubes (about 1 x 7 cm).
- Graduates pipettes.
- Automatic micropipettes with disposable tips (0.1 ml).

- Vortex mixer.
- Gamma-Counter.
- Multi-sample centrifuge
- Automatic dispenser for addition of PEG.

### 3. Test procedure

Allow the reagents to reach room temperature (18-25°C) and mix gently before the use.

- For each assay prepare 5 groups of tubes:
    - 2 tubes for total radioactivity.
    - 2 tubes for NSB (non specific binding).
    - 2 tubes for B0 ("0" concentration of cold antigen).
    - 2 tubes for each calibrator.
    - 2 tubes for each sample and Control.
  - Pipette into the respective tube 0.1 ml of sample, calibrator and Control .
  - Pipette into B0 and NSB respectively 0.1 and 0.2 ml of calibrator zero.
  - Add into all tubes 0.1 ml of <sup>125</sup>I-RT<sub>3</sub>.
  - Add into all tubes except NSB and Total Radioactivity tubes 0.1 ml of RT<sub>3</sub> Antiserum.
  - Mix and incubate for 3 hours at room temperature (18-25°C).
  - Add into all tubes except into total radioactivity tubes 1.0 ml of 20% TW PEG Solution.
  - Mix on vortex and centrifuge all tubes except total Radioactivity tubes.
- Centrifuge for:
- 25 minutes at 2000 x g;
  - 20 minutes at 2500 x g;
  - 15 minutes at 3500 x g;
- $$g = 1.1178 \times 10^{-5} \times r \times n^2; \quad (g=\text{gravity}; n = \text{number of turns/minute}; r = \text{radius in cm}).$$
- Eliminate carefully the supernatant by decantation or aspiration (except total radioactivity).
  - Count in Gamma-Counter for at least 1 minute.

### 4. Assay Procedure Scheme

The volumes are expressed in ml

Reagents	Tubes	Sample	B0	Calibrator	NSB	Total Radioactivity
Samples or Serotest	0.1	-	-	-	-	-
Calibrator zero	-	0.1	-	-	0.2	-
Calibrators	-	-	-	0.1	-	-
<sup>125</sup> I-RT <sub>3</sub>	0.1	0.1	0.1	0.1	0.1	0.1
RT <sub>3</sub> Antiserum	0.1	0.1	0.1	0.1	-	-
Mix and incubate for 3 hours at room temperature (18-25°C).						
20% TW PEG Solution		1.0	1.0	1.0	1.0	-
Mix on vortex and centrifuge. Aspirate or decant (except from total radioactivity). Count all tubes in Gamma-Counter for at least 1 minute.						

## 8. RESULTS

NSB mean counts may be subtracted from all mean counts except from total radioactivity. B0 mean counts are used to calculate the % binding in absence of cold antigen (maximum binding).

$$\frac{\text{B0 mean counts}}{\text{Total radioactivity mean counts}} \times 100 = \% \text{ binding}$$

% relative binding of calibrator and samples is calculated with the following formula:

$$\frac{\text{Mean counts (calibrator, sample or Serotest)}}{\text{B0 mean counts}} \times 100 = \% \text{ relative binding}$$

Draw the dose-response curve by plotting the % relative binding of each calibrator (y axis) against the relative concentration (x axis) using semilogarithmic paper.

By interpolating on the standard curve the % relative binding of each sample, one obtains the sample antigen concentration in ng/ml. For diluted samples it is necessary to multiply the value read on standard curve by dilution factor.

For data computerization we suggest to refer to bibliography (8,9). The user can choose, according to assay characteristic, the calculation method and the graphic expression which allows the best data treatment.

### 1. Example of calculation - Representative Standard Curve

Tubes	Mean cpm	B/T (%)	B/B0 (%)	Concentrations (ng/ml)
Total Radioactivity	19649,1	-	-	-
NSB	652,5	3,32	-	-
B0	7869,2	36,73	-	-
C.1 0.02 ng/ml	7322,0	-	92,42	-
C.2 0.06 ng/ml	6499,8	-	81,02	-
C.3 0.12 ng/ml	5354,3	-	65,15	-
C.4 0.33 ng/ml	4059,3	-	47,21	-
C.5 0.65 ng/ml	2941,4	-	31,72	-
C.6 1.15 ng/ml	2302,8	-	22,87	-
C.7 2.14 ng/ml	1703,4	-	14,56	-
Control	3141,6	-	33,40	0.57

Do not use this example instead of data obtained during the assay. For the exact value, refer to the vial label.

### INTERNAL QUALITY CONTROL

- If the results obtained for Control are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises

### 9. LIMITATIONS

1. To obtain reliable results it is necessary to follow carefully the instruction reported in this package insert.
2. Do not use the reagents beyond the expiry date stated on vial labels.
3. Do not mix reagents from different lots.
4. After reconstitution, store the reagents at 4°C for 3 days in the stoppered vials. After reconstitution and storage at 4°C, the reagents can be used only once. For longer than 3 days storage, freeze the reagents at -20°C. After freezing, allow the reagents to reach room temperature (18-25°C) gradually, avoiding to use water-bath at 37°C. Avoid to thaw the reagents more than once.
5. Samples contaminated with radioactivity may give inaccurate results.

### 10 EXPECTED VALUES

	Concentration (ng/ml)
Normal men and women : 77 samples	0.17-0.44

We suggest that each laboratory establishes its own normal range. Conversion into international units or nmol/l: nmol/l = ng/ml • 1.54.

### 11. PERFORMANCE CHARACTERISTICS

#### 1. Accuracy

Recovery studies were performed by adding rT<sub>3</sub> serum or plasma samples. RT<sub>3</sub> values were determined before and after addition, and the % recovery of added rT<sub>3</sub> calculated.

Sample	rT <sub>3</sub> added (ng/ml)	rT <sub>3</sub> measured (ng/ml)	% Recovery
1	0.83	0.70	84.4
	0.36	0.34	95.1
	0.13	0.11	81.1
	0.06	0.06	94.0
	0.02	0.02	107.9
2	0.96	0.81	85.1
	0.52	0.47	90.0
	0.15	0.14	91.5
	0.07	0.06	80.1

#### 2. Dilution

In a dilution study plasma samples were diluted with calibrator zero and the recovery determined.

Sample	Dilution factor	rT <sub>3</sub> measured (ng/ml)	% Recovery
1	1	1.79	100.0
	2	0.94	95.3
	4	0.44	102.1
	8	0.21	106.1
2	1	1.65	100.0
	2	0.48	97.4
	4	0.25	93.7
	8	0.12	100.4

#### 3. Precision

Intra-Assay coefficient of variation was evaluated using three samples measured 12 times in the same assay. The within assay variability is shown below:

Sample	Mean ng/ml	Standard Deviation ng/ml	C. V. %
A	0.67	0.03	4.1
B	0.26	0.01	4.0
C	0.10	0.01	9.7

Inter-Assay coefficient of variation (C.V.) was evaluated using three samples run in 10 different assays. The between assay variability is shown below:

Sample	Mean ng/ml	Standard Deviation ng/ml	C. V. %
A	0.62	0.04	5.7
B	0.28	0.02	7.4
C	0.14	0.01	10.5

#### 4. Sensitivity

The LoB (Limit of Blank) was calculated by measuring the blank 20 times and was calculated as the mean + 2 Standard Deviations of the distribution of the test values ; the calculated LoB is 0.014 ng/ml.

The LoD (Limit of Detection) was calculated as the LoB + 1.645 standard deviation of a low concentration sample tested in 10 different run .The calculated LoD is 0.025 ng/ml.

The LoQ (Limit of Quantification) was calculated by testing 5 low values samples, 9 times. The calculated LoQ is 0.035 ng/ml.

#### 5. Specificity

The specificity has been evaluated by the interference of the following compounds, similar to rT<sub>3</sub> % of interference has been calculated according to Abraham's method (x/y • 100) where x and y are respectively the weight of rT<sub>3</sub> and of interfering compound which causes a decrease of 50% in the binding capacity:

Compound	% Cross-Reactivity
Reverse T <sub>3</sub>	100
L-Thyroxine (T <sub>4</sub> )	0.22
3,5,3'- Triiodothyronine (T <sub>3</sub> )	0.024
3,3'- Diiodo-L-Thyronine	15
3,5- Diiodo-L-Thyronine	0.0002
3,5,3'- Triiodothyroacetic Acid (Triac)	0.023
3- Monoiodothyrosine (MIT)	0.0001
3,5- Diiodothyrosine (DIT)	0.0048
Acetylsalicylic Acid	<0.0001
Phenylbutazone	<0.0001

The potentially interfering effects of hemoglobin at 500 mg/dl, of bilirubin at 100 mg/dl and of triglycerides at 250 mg/dl have been evaluated. The results of this test do not demonstrate any significant interference as shown in the table below.

<b>Interfering substance</b>	<b>Serum sample (ng/ml)</b>	<b>Plasma sample (ng/ml)</b>
-	0.53	0.18
Haemoglobin	0.55	0.18
Bilirubin	0.53	0.15
-	0.58	0.15
Triglycerides	0.51	0.17

Other translations of this Instruction for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

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