Fecal Rotavirus and Adenovirus Duo Antigen Test Kit

RAPEPKT926

History

Summary of change :

Previous Version :	Current Version :
110225/1	200224/1
No History	History added
Old Diasource logo	New Diasource logo
Multilanguage IFU	Addition of the following sentence at the end of the English IFU: "Other translations of this Instructions for Use can be downloaded from our website: https://www.diasource-diagnostics.com/ "
No IVD symbol	IVD symbol added
LOT : 110225/1	Version: 200224/1
PI number : 1701148	No PI number
No manufacturer symbol	Manufacturer symbol added

Fecal Rotavirus and Adenovirus Duo Antigen Kit **en**

Rapid Immunochromatographic Test Device for the detection of
Rotavirus and Adenovirus Antigens in Feces

RAPEPKT926

IN VITRO DIAGNOSTIC

DIAsource ImmunoAssays SA - Rue du Bosquet 2, B-1348 Louvain-la-Neuve, Belgium - Tel: +32 10 84 99 11 - Fax : +32 10 84 99 91

INTENDED USE

This rotavirus and adenovirus Duo antigen test kit is intended for the direct qualitative detection of the presence of rotavirus and/or adenovirus antigens in patient fecal samples. The test might be used as an aid for the diagnosis of an adenovirus and/or rotavirus infection resulting in acute gastroenteritis.

SUMMARY OF PHYSIOLOGY

Rotaviruses are the main cause of acute gastroenteritis and diarrhea, especially in children under the age of two years. Rotaviruses have been identified in almost 40% of the feces of children with gastroenteritis. Rotavirus is the cause of up to 50% of the hospitalized cases of diarrhea in infant and young children. If not treated, the infection may result in severe dehydration and disorders of body electrolyte balance. Therefore, it can be mortal in risk populations such as children, the elderly or immunosuppressed individuals. Rotavirus is transmitted by oral-fecal contact with an incubation period of 1-3 days. Characteristic symptoms include vomiting, hydrodiarrhoea for between 3 and 8 days, high temperature and stomach pains. A large amount of rotavirus particles is shed during infection.

Adenoviruses are one of the main causes of acute gastroenteritis and diarrhea, especially in children under the age of two years. Adenoviruses have been identified in almost 12% of the feces of children with gastroenteritis. It was reported that adenovirus is the second leading cause of the hospitalized cases of diarrhea in infant and young children. If not treated, the infection may result in severe dehydration and disorders of body electrolyte balance. Therefore, it can be mortal in risk populations such as children, the elderly or immunosuppressed individuals. Adenovirus is transmitted by oral-fecal contact, but can result from the inhalation of aerosols as well. Its incubation period lasts 5 to 8 days. Characteristic symptoms include vomiting, hydrodiarrhoea, high temperature and stomach pains. There are 41 known human adenoviruses primarily differentiated by serology and DNA analysis. Morphologically, the viruses are non-enveloped icosahedral structured with a diameter of about 80 nm.

Diagnosis of gastroenteritis with rotavirus and/or adenovirus infection can be established based on the detection of the virus particles by electron microscopy or the virus antigen by specific immunoassay methods.

ASSAY PRINCIPLE

This is a two-in-one test including a rotavirus antigen test strip and an adenovirus antigen test strip that are back-to-back positioned in one test tube.

The rotavirus and adenovirus Duo antigen rapid test strips employ two group paired monoclonal antibodies for either rotavirus antigen or adenovirus antigen. Dye-conjugated monoclonal antibodies against antigen VP6 of group A rotavirus, and solid-phase specific rotavirus antibodies together form one of the antibody pairs. The antibodies specifically against human adenovirus antigens and a solid-phase coated monoclonal antibody specifically against adenovirus hexon antigens form the other pair of antibodies. In this test the specimen is first treated with an extraction solution to extract rotavirus and/or adenovirus antigens from the feces. Following extraction, the only step required is to screw the rotavirus and adenovirus Duo test strip tube into the sample collection device. As the sample extraction flows through the chamber and reach the test strips, the colored particles migrate. In the case of a positive result the specific antibodies present on the membrane will capture the colored particles. Different colored lines will be visible, depending upon the virus content of the sample. These lines, after 5 minutes of incubation at room temperature, are used to interpret the result.

REAGENTS: Preparation and Storage

This test kit must be stored at 2 – 8°C upon receipt. For the expiration date of the kit refer to the label on the kit box. All components are stable until this expiration date.

Prior to use allow all reagents to come to room temperature if the kit is stored at refrigerated condition.

- **TUBE** 30 Fecal specimen collection device: containing sampling tube, sampling lid and pre-added extraction solution (1.1 ml Tris Buffer containing 0.1% bovine serum albumin) in

the sampling tube. This device should be stored at 2 to 8°C. Do not freeze.

- **STRIP** 30 Test strip: one dipstick for the Fecal Rotavirus and Adenovirus Duo antigen test is assembled in a transparent housing and sealed in a foil pouch with desiccant. It should remain in its original sealed pouch until ready for use. The test strip should be stored at 2 to 30°C. Do not freeze.
- 30 Disposable spoons for sample collection
- Instruction for use.

SAFETY PRECAUTIONS

The reagents are for professional use only. Source material from which reagents of bovine serum was derived in the contiguous 48 United States. It was obtained only from donor health animals maintained under veterinary supervision and found free of contagious diseases. Wear gloves while performing this test and handle these reagents and patient samples as if they are infectious. Do not get in eyes, on skin, or on clothing. Do not ingest. On contact, flush with copious amounts of water for at least 15 minutes. When the assay procedure is completed, dispose of specimens (biohazard materials) carefully after autoclaving for at least one hour. Alternatively, treat with a 0.5 or 1% solution of sodium hypochlorite for one hour before disposal. Use Good Laboratory Practices.

MATERIALS REQUIRED BUT NOT PROVIDED

1. Disposable pipette for watery sample collection
2. Positive Control

SPECIMEN COLLECTION


1. Stool specimens can be collected at any time of the day.
2. Collect a random sample of feces in a clean, dry cup or toilet paper.
3. Unscrew the sampling lid and keep the sampling tube in a vertical position to prevent the loss of any extraction solution (figure 1).
4. Insert and twist the tip of the sampling lid into the stool specimen at two or more different sites (figure 2).
5. Collect fecal sample that is stuck to the surface of the sampling lid. Do not intentionally collect any separate and large pieces of fecal sample into the tube.
6. Replace the sampling lid into the tube and secure tightly (figure 3).
7. The specimen is ready for testing, transportation or storage. It can be stored at 2-8°C for up to 14 days and at room temperature for up to 5 days.

Note: Two specimens from three consecutive bowel movements are recommended from American Cancer Society. Specimen should not be collected during digital rectal examination.

TEST PROCEDURE

1. Bring the sealed foil pouch test strips and collected specimens to room temperature.
2. Shake the sampling tube vigorously to ensure a good liquid suspension.

3. Position the sampling tube upside down vertically and let it settle for about 1 minute.
4. Remove the test strip from the sealed foil pouch (figure A).
5. Screw the test strip in a vertical position into the sampling tube by **breaking** the bottom seal of the sampling tube. Secure tightly! (Figure B)
6. Allow the solution to flow into the bottom space of the test strip and keeping the device in a vertical position.
7. Read test result at 5 minutes. Do not interpret test result after 10 minutes.

(3) Rotavirus and Adenovirus Antigen Positive

Test lines of both rotavirus and adenovirus test strips are observed.


- **Negative:**

If both rotavirus and adenovirus test lines have no colored band and the control line displays a red/pink colored band, the test result is negative.

INTERPRETATION OF RESULTS

- **Valid Test:**


The appearance of control line on the test reading membrane area indicates that the test is valid. The control line is located at the upper half of the test membrane area.

- **Positive:**

If two red/pink colored bands are visible on either strip within 10 minutes, the test result is positive for the respective Adenovirus and/or Rotavirus.

(1) Adenovirus Antigen Positive (Green tape covered strip)

(2) Rotavirus Antigen Positive (Blue tape covered strip)

QUALITY CONTROL

Good laboratory practices recommend the use of appropriate controls. There are two types of controls for the Rotavirus and Adenovirus Duo antigen test, the internal procedural control and external controls.

1. **Internal procedural control:** Each fecal rotavirus and adenovirus Duo antigen test contains of a built-in procedural control. A coloured control line will appear if the test has been performed correctly and the reagents are reactive. This control does not ensure that the test line antibody is accurately detecting the presence or absence of viral antigen in the test fecal sample.
2. **External controls:** It is recommended to use external positive controls. The external positive controls are not provided with this kit. External controls are used to assure that the test line antibodies are reactive. However, external controls will not detect an error in the performance of a particular test with a patient sample. It is recommended that the external control be tested once per kit.

You should always follow local, state, federal guidelines for running quality control.

LIMITATION OF THE PROCEDURE

1. The test should be used only for the detection of rotavirus and/or adenovirus antigens in fecal samples.
2. The test is qualitative and no quantitative interpretation should be made with respect to the intensity of the positive line, when reporting the result
3. Over two hundred samples were evaluated to assure the correct performance of the test. The correlation of the results with other techniques (ELISA) was satisfactory. However, interferences in the performance of the tests should not be excluded.
4. No cross-reactions with other viruses or substances were observed during the evaluation of the test. A negative result does not totally exclude a possible rotavirus or adenovirus infection. The significance of the results must be evaluated in relation to the patient's clinical symptoms.
5. As with all diagnostic tests, the definitive clinical diagnosis must not be based on the result of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated. The fecal rotavirus and adenovirus Duo antigen test is designed for the aid of clinical diagnosis and should not replace other diagnostic procedures.

PERFORMANCE

The sensitivity and specificity of this adenovirus antigen test device are studied with 212 clinical samples and compared with an adenovirus antigen ELISA test. The results are shown below.

CARE Strip \ ELISA	Positive	Negative
	Positive	61
Negative	1	149
Total	62	150

Sensitivity: 98% (61/62 = 98.4%)

Specificity: 99% (149/150 = 99.3%)

Accuracy: 99% (210/212 = 99.1%)

Inter-series and intra-series accuracy: 100 %

The sensitivity and specificity of this rotavirus antigen test device are studied with 206 clinical samples and compared with a rotavirus antigen ELISA test. The results are shown below.

ELISA	Positive	Negative	Total
CARE Strip	68	2	70
Positive	2	134	136
Negative	70	136	206
Total			

Specificity: 98.5 % (134/136)

Sensitivity: 97.1 % (68/70)

Accuracy: 98.1% (202/206)

Inter-series and intra-series accuracy: 100 %

Interference: Cross reactivity has been evaluated and found to be negative compared to positive specimens of *Cryptosporidium parvum* and *Giardia lumbli*a.

LITERATURE

1. Nishio O, Ooseto M, Takagi K, Yamasita Y, Ishihara Y, Isomura S. Enzyme-linked immunosorbent assay employing monoclonal antibodies for direct identification of enteric adenoviruses (Ad40,41) in feces. Microbiol Immunol. 1990;34(10):871-7.
2. Vizzi E, Ferraro D, Cascio A, Di Stefano R, Arista S. Detection of enteric adenoviruses 40 and 41 in stool specimens by monoclonal antibody-based enzyme immunoassays. Res Virol. 1996 Nov-Dec;147(6):333-9.
3. Singh-Naz N, Rodriguez WJ, Kidd AH, Brandt CD. Monoclonal antibody enzyme-linked immunosorbent assay for specific identification and typing of subgroup F adenoviruses. J Clin Microbiol. 1988 Feb;26(2):297-300.
4. Uhnnoo I, Wadell G, Svensson L, Johansson ME. Importance of enteric adenoviruses 40 and 41 in acute gastroenteritis in infants and young children. J Clin Microbiol. 1984 Sep;20(3):365-72.
5. Uhnnoo I, Svensson L, Wadell G. Enteric adenoviruses. Baillieres Clin Gastroenterol. 1990 Sep;4(3):627-42.
6. Shinozaki T, Araki K, Fujita Y, Kobayashi M, Tajima T, Abe T. Epidemiology of enteric adenoviruses 40 and 41 in acute gastroenteritis in infants and young children in the Tokyo area. Scand J Infect Dis. 1991;23(5):543-7.
7. D. Van Beers , M. DE Foor , R. Viehoff , D. Col , M. Venuti and T. Set-up of a new rapid immunochromatographic diagnostic test for a Rotavirus detection. Leclipteux. Progress in Clinical Virology III , Bologne , Septembre 1997.
8. Sneyers et al. Detection of rotavirus in faecal specimens with a monoclonal antibody enzyme-linked immunosorbent assay : comparison with polyclonal antibody enzyme-immunoassays and a latex agglutination test. Comp. Immun. Microbiol. Infect. Dis. , vol 12 , n°4 , pp 95-104 , 1989
9. I. Van der Donck et al. Comparison of Three Rapid Immunoassays for the Detection of Rotavirus Antigen in Stool Samples ESCV Winter Meeting 1999, Rotterdam, the Netherlands
10. I. Wilhelmi et al. Evaluacion de tres Metodos de Deteccion de Rotavirus en Heces 6th Congreso Nacional de Virologia, Madrid, 26th Oct. 99

Consult instructions for use

Manufacturer

Storage temperature

Contains sufficient for n tests

Use by

In vitro diagnostic medical device

Batch code

Test Strip

Catalogue number

Collection device

Other translations of this Instructions for Use can be downloaded from our website:

<https://www.diasource-diagnostics.com/>

Revision date : 2020-02-24

Rapid Immunochromatographic Test Device for the detection of Rotavirus and Adenovirus Antigens in Feces


RAPEPKT926 - Instructions for Fecal Sample Collection

IN VITRO DIAGNOSTIC

DIAsource ImmunoAssays SA - Rue du Bosquet 2, B-1348 Louvain-la-Neuve, Belgium - Tel: +32 10 84 99 11 - Fax : +32 10 84 99 91


1

1. Take the sampling tube and unscrew the sampling lid, keeping the sampling tube in a **vertical position** to prevent loss of solution.


2

1. Hold the sampling lid by the **Thumb Grip**.
2. Use the **tip** of the sampling lid to collect a small amount of fecal sample at two or more sites. Only take the fecal sample that sticks to the sampling lid tip (never intentionally place any separate piece of fecal sample into the tube).

3

1. Insert and screw the sampling lid back into the sampling tube in a **vertical position**. Do not spill any solution from the tube.
2. Tightly seal the lid with the tube.

Rotavirus and Adenovirus Duo Antigen Kit

en

Rapid Immunochromatographic Test Device for the detection of Rotavirus and Adenovirus Antigens in Feces


RAPEPKT926 - Instructions for Test Procedures

IN VITRO DIAGNOSTIC

DIAsource ImmunoAssays SA - Rue du Bosquet 2, B-1348 Louvain-la-Neuve, Belgium - Tel: +32 10 84 99 11 - Fax : +32 10 84 99 91


A

1. Shake to dissolve the stool into solution.
2. Turn the sampling tube **upside down vertically**.
3. Remove the test strip from foil pouch.

B

1. Insert and screw the test strip **in a vertical position** into the sampling tube by breaking the bottom seal of the sampling tube.
2. Allow the solution to flow into the bottom space of test strip, keeping the device **in a vertical position**.
3. You may soon see a red fluid moving across the white area of the test strip. Read test result after 5 minutes.

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