

ELISA SHAKER

The ELISA Shaker, allows shaking and incubation simultaneously.

CAT#: DIA4000

Special Features

- Simultaneous shaking & incubation operation.
- Buzzer indication on completion of incubation
- Indication of Remaining time
- Current temperature of incubation on display, on pressing TEMP Key.
- Separate Timer ON indication on Keypad



Technical Specification

Operating Modes	Shaker, Incubator, Shaking & Incubator
-----------------	--

Temperature Control

Temperature Range	37° C to 42° C
-------------------	----------------

Resolution	1° C
------------	------

Incubation Time	1 to 999 min
-----------------	--------------

Shaker

Frequency	400 to 700 RPM
-----------	----------------

Amplitude	2 mm
-----------	------

Operating Position	On horizontal flat, rigid and vibration free surface
--------------------	--

Operating Conditions

Temperature	From + 18° C to 32° C
-------------	-----------------------

Relative Humidity	Up to 80%
-------------------	-----------

Storage Conditions

Temperature	From -10° C to 35° C
-------------	----------------------

Relative Humidity	Up to 80%
-------------------	-----------

Enclosure	ABS Fire retardant
-----------	--------------------

Size (cm)	28 X 25 X 16 (l X b X h)
-----------	--------------------------

Weight	3 Kgs (Approx)
--------	----------------