

Rabbit Crude anti Testosterone



PRODUCT SPECIFICATION

Cat # 1302535-V

Batch*:	<i>[Batch dependent]</i>	
Size	<i>[≥ 100 µL]</i>	
Format:	Undiluted neat serum	
Purification:	None	
Preservative:	NaN3 <0,1%	
Immunogen:	Testosterone 3-CMO-BSA/KLH	
Host:	Rabbit	
Cross-reactivity:	Testosterone:	100%
	17OH Progesterone:	<0.01%
	Androstenedione:	0.3%
	Cortisol:	< 0.01%
	DHEA:	<0.01%
	DHEA-SO4:	<0.01%
	DHT:	0.3%
	Estradiol:	<0.01%
	Progesterone:	0.01%
	Applications:	Immunoassays
Storage:	-20°C	
Recommended dilution:	1/800.000	
Expiry date:	<i>[Batch dependent]</i>	

* Only the figures before the last slash should be taken into account for the batch identification.

Date

Signature