


mAb Hu Leptin Detecting, F(ab)'2, HRP Conjugated



PRODUCT SPECIFICATION

Cat # **5122818**

Batch*:	<i>[Batch dependent]</i>
Size	<i>[≥ 0,1 mL]</i>
Format:	Liquid
Buffer:	Tris-HCL 0,05 M BSA 5 g/L Thymol 0,5 g/L
Applications:	Immunoassays
Storage:	-20°C
Expiry date:	<i>[Batch dependent]</i>

* Only the figures before the last slash should be taken into account for the batch identification.

Date

Signature