

mAb Hu PG Capture, F(ab)'2, Purified Unconjugated



PRODUCT SPECIFICATION

Cat # **5314627**

Clone:	1R1 14A6 3B2
Batch*:	<i>[Batch dependent]</i>
Size	<i>[≥ 1 mg]</i>
Format:	F(ab)'2 Unconjugated, Liquid
Concentration:	<i>[≥ 1 mg/mL, Batch dependent]</i>
Buffer:	BBS 0,1M pH 8,0
Preservative:	None
Purification:	Protein A affinity chromatography
Isotype:	IgG1
Immunogen:	Human recombinant PG
Mouse / fusion partners:	Balb/c (Charles River)/ NSO
Specifications :	Recognizes natural and recombinant Human PG
Applications:	Immunoassays
Storage:	2-8°C
Expiry date:	<i>[Batch dependent]</i>

* Only the figures before the last slash should be taken into account for the batch identification.

Date

Signature