

## mAb PTH Capture Purified Unconjugated



### PRODUCT SPECIFICATION

Cat # **5314926**

Clone:	OBP 1
Batch*:	<i>[Batch dependent]</i>
Size	<i>[≥ 1 mg]</i>
Format:	Purified Unconjugated, Liquid
Concentration:	<i>[≥ 1 mg/mL, Batch dependent]</i>
Buffer:	PBS 0.01M 0,85% NaCl pH 7,2
Preservative:	NaN3 0,05%
Purification:	Protein A affinity chromatography
Isotype:	IgG1
Immunogen:	Synthetic PTH [1-34]
Host:	Mouse
Specifications :	No cross reactivity with PTH [39-84]
Storage:	2-8°C
Expiry date:	<i>[Batch dependent]</i>

\* Only the figures before the last slash should be taken into account for the batch identification.

Date

Signature