


mAb Hu IGF-1, Purified Unconjugated



PRODUCT SPECIFICATION

Cat # **5315836**

Clone:	OBP 6
Batch*:	<i>[Batch dependent]</i>
Size	<i>[≥ 1 mg]</i>
Format:	Purified Unconjugated, Liquid
Concentration:	<i>[≥ 1 mg/mL, Batch dependent]</i>
Buffer:	PBS 15mM NaCl 0,15M pH 7,2
Preservative:	NaN3 0,09%
Purification:	Protein A affinity chromatography
Isotype:	IgG1 K
Immunogen:	IGF-1
Host:	Mouse
Specifications :	No cross-reactivity with IGF-2
Storage:	2-8°C
Expiry date:	<i>[Batch dependent]</i>

* Only the figures before the last slash should be taken into account for the batch identification.

Date

Signature