


mAb Hu TSH Capture, Purified Unconjugated



PRODUCT SPECIFICATION

Cat # 5318806

Clone:	103B 2/2B4 1A8
Batch*:	<i>[Batch dependent]</i>
Size	<i>[≥ 1 mg]</i>
Format:	Purified Unconjugated, Liquid
Concentration:	<i>[≥ 1 mg/mL, Batch dependent]</i>
Buffer:	BBS 0,1M pH 8,0
Preservative:	None
Purification:	Protein A affinity chromatography
Isotype:	IgG1 K
Immunogen:	Human recombinant TSH
Mouse / fusion partners:	Balb/c (Charles River)/ NSO
Specifications :	Recognizes natural and recombinant Human TSH, epitope III (1)
Applications:	Immunoassays
Storage:	2-8°C
Expiry date:	<i>[Batch dependent]</i>

* Only the figures before the last slash should be taken into account for the batch identification.

Date

Signature

(1) Donadio S et al.

Both core and terminal glycosylation alter epitope expression in thyrotropin and introduce discordances in hormone measurements.

Clin Chem Lab Med 2005;43(5):519–530