


## mAb Hu 1,25(OH)<sub>2</sub>-Vit D Capture, Purified Unconjugated



### PRODUCT SPECIFICATION

Cat # **5319256**

|                          |   |
|--------------------------|---|
| Clone:                   | EG4   |
| Batch*:                  | <i>[Batch dependent]</i>                          |
| Size                     | <i>[≥ 1 mg]</i>                                   |
| Format:                  | Purified Unconjugated, Liquid                     |
| Concentration:           | <i>[≥ 1 mg/mL, Batch dependent]</i>               |
| Buffer:                  | PBS 0,1M pH 7,4                                   |
| Preservative:            | NaN <sub>3</sub> 0,5 g/L                          |
| Purification:            | Protein A affinity chromatography                 |
| Isotype:                 | NA  |
| Immunogen:               | 1,25(OH) <sub>2</sub> -Vit D – BSA / KLH          |
| Mouse / fusion partners: | Balb/c (Charles River)/ NSO                       |
| Specifications :         | Recognizes 1,25(OH) <sub>2</sub> – D <sub>3</sub> |
| Applications:            | Immunoassays                                      |
| Storage:                 | 2-8°C   |
| Expiry date:             | <i>[Batch dependent]</i>                          |

\* Only the figures before the last slash should be taken into account for the batch identification.

Date

Signature