

Crude Serum anti Aldosterone



PRODUCT SPECIFICATION

Cat # **5331225**

Batch*:	<i>[Batch dependent]</i>
Size	<i>[≥ 100 µL]</i>
Format:	Undiluted neat serum
Aspect:	Yellow clear liquid
Preservative:	NaN ₃ 0,2 g/L
Purification:	None
Immunogen:	Aldosterone 3-CMO-BSA/KLH
Host:	Rabbit L307A

Cross-reactivity :

11-Deoxycortisol	<0.05%
11-Deoxycorticosterone	0.30%
Androstenedione	<0.05%
Androsterone	<0.05%
Cortisol	<0.05%
Dihydrotestosterone (DHT)	<0.05%
Testosterone	<0.05%
17-Hydroxyprogesterone	<0.05%
Progesterone	<0.05%
DHEA	<0.05%
Estriol (E3)	<0.05%
Estrone	<0.05%

Applications:	Immunoassays
Recommended dilution:	1/100.000
Storage:	-20°C
Expiry date:	<i>[Batch dependent]</i>

* Only the figures before the last slash should be taken into account for the batch identification.

Date

Signature